Chemical Analysis Results for Potable Water from ISS Expeditions 21 through 25

John E. Straub II¹, Debrah K. Plumlee², and John R. Schultz³ Wyle Integrated Science & Engineering Group, Houston, Texas, 77058

J. Torin McCoy⁴
NASA Johnson Space Center, Houston, Texas, 77058

The Johnson Space Center Water and Food Analytical Laboratory (WAFAL) performed detailed ground-based analyses of archival water samples for verification of the chemical quality of the International Space Station (ISS) potable water supplies for Expeditions 21 through 25. Over a 14-month period the Space Shuttle visited the ISS on four occasions to complete construction and deliver supplies. The onboard supplies of potable water available for consumption by the Expeditions 21 to 25 crews consisted of Russian ground-supplied potable water, Russian potable water regenerated from humidity condensate, and US potable water recovered from urine distillate and condensate. Chemical archival water samples that were collected with U.S. hardware during Expeditions 21 to 25 were returned on Shuttle flights STS-129 (ULF3), STS-130 (20A), STS-131 (19A), and STS-132 (ULF4), as well as on Soyuz flights 19-23. This paper reports the analytical results for these returned potable water archival samples and their compliance with ISS water quality standards.

Nomenclature

CE	Capillary Electrophoresis
CWC	Contingency Water Container
DAI	Direct Aqueous Injection
DMSD	Dimethysilanediol

DWEL Drinking Water Exposure Limit
EPA Environmental Protection Agency
GC/MS Gas Chromatography/Mass Spectrometry

HA Health Advisory
IC Ion Chromatography

ICP/MS Inductively Coupled Plasma/Mass Spectrometry

ISE Ion Selective Electrode
ISS International Space Station
JSC Johnson Space Center
LC Liquid Chromatography

LC/MS Liquid Chromatography/Mass Spectrometry

LCV Leuco Crystal Violet

MCL Maximum Contaminant Level

MORD Medical Operations Requirements Document

N/A Not Applicable NA Not Analyzed

NASA National Aeronautics & Space Administration

NTU Nephelometric Turbidity Unit PWD Potable Water Dispenser RIP Rack Interface Panel

¹ Senior Engineer, Water and Food Analytical Laboratory, NASA Johnson Space Center, Mail Stop: Wyle/HEF/37A.

² Senior Scientist, Water and Food Analytical Laboratory, NASA Johnson Space Center, Mail Stop: Wyle/HEF/37A.

³ Senior Supervisor, Water and Food Analytical Laboratory, NASA Johnson Space Center, Mail Stop: Wyle/HEF/37A.

⁴ Environmental Operations Lead, Environmental Factors Branch, NASA Johnson Space Center, Mail Stop: SF2

SM Service Module

SRV-K System for Regeneration of Condensate Water SVO-ZV System for Water Storage and Dispensing SWEG Spacecraft Water Exposure Guideline

TDS Total Dissolved Solids
THM Trihalomethanes

TOCA Total Organic Carbon Analyzer

U.S. United States

UPA Urine Processor Assembly

UV Ultraviolet

WAFAL Water and Food Analytical Laboratory

WPA Water Processor Assembly WRS Water Recovery System

Introduction

While the International Space Station (ISS) recently celebrated the 10-year anniversary of human occupation by some 25 different expedition crews, the Johnson Space Center's (JSC) Water and Food Analytical Laboratory (WAFAL) quietly continued its key role of performing ground-based chemical analyses of returned archival water samples to verify that the ISS potable water remained safe for crew consumption. Chemical analysis results for samples collected during Expeditions 1-20 have been previously published. This paper presents and discusses the analytical results from chemical analyses of the archival potable water samples collected and returned during Expeditions 21 - 25, as detailed in Table 1. Only chemical analysis results are discussed herein, as the JSC Microbiology Laboratory has responsibility for microbiological monitoring of the ISS potable water supplies.

The WAFAL sent a representative to the Kennedy Space Center in Florida to retrieve and return the ISS potable water samples that returned on Shuttle flights STS-129 (ULF3), STS-130 (20A), STS-131 (19A), and STS-32(ULF4). These samples were packed into coolers along with ice packs and transported to Houston via commercial air transport. The ISS water samples that returned on Soyuz 19 through Soyuz 23 were retrieved from the vehicle by Russian specialists and subsequently transferred to a NASA representative for cold packing and transportation to Houston via NASA jet. The WAFAL representative took possession of the return after their arrival in Houston for delivery to the JSC laboratory.

Once the samples arrived in the WAFAL, allocation was performed based upon return sample volume. The samples collected into U.S. 1-L Teflon® sample bags typically contained sufficient sample volume (> 500 mL) to support full chemical characterization using the standard and custom analytical methods identified in Table 2.

Water samples collected for return on Soyuz were typically less than 500 mL volume, requiring elimination of some analyses and/or reductions in sensitivity of other analyses performed. During data analysis and reporting, the analytical results for each sample are compared as appropriate against either the Russian Segment potable water quality requirements found in the *ISS Medical Operations Requirement Document* (MORD)⁹ or the U.S. Segment potable water quality requirements found in the *System Specification for the ISS* document.¹⁰

Background

Over the 14-month period beginning October 2009, the Expeditions 21 to 25 crews onboard the International Space Station had access to 3 different sources of potable water: U.S. potable water recovered from urine distillate and humidity condensate, Russian ground-supplied potable water, and Russian potable water regenerated from humidity condensate. These different water supplies were available via water systems located in the U.S. Segment and the Russian Segment of the ISS.

U.S. Segment Water Recovery System

The U.S. Segment Water Recovery System (WRS) was delivered on STS-126 (ULF2) in November of 2008 and began processing a combined wastewater feed consisting of both urine distillate and humidity condensate. During an initial checkout period of about 90-days the product water was tested to verify quality and was subsequently approved for crew consumption beginning in May 2009 prior to the start of ISS 6-person crew operations.

Expedition	Table 1. Summa Flight No.	Samples	Sample Type	Sample Collection Date	Sample Receipt
pea_e	1.1.5	Received	Sumpre 1, pe	Sumple Content Suit	Date
	STS-129 (ULF3)	1	PWD Ambient	11/10/2009	11/30/2009
		1	PWD Hot		
		1	SVO-ZV		
		1	SRV-K Hot		
	Subtotal:	4			•
21	Soyuz 19	1	PWD Ambient	10/20/2009	12/14/2009
		1	PWD Hot		
		1	SVO-ZV		
	Subtotal:	3	•		<u>-</u>
	Total:	7			
	STS-130 (20A)	1	PWD Ambient	1/6/2010	2/22/2010
		2	SVO-ZV	1/6/2010, 2/3/2010	
		1	SRV-K Hot	2/3/2010	
		1	SRV-K Warm	1/6/2010	
22	Subtotal:	5			
22	Soyuz 20	1	SVO-ZV	3/3/2010	3/26/2010
		1	SRV-K Hot		
		1	SRV-K Warm		
	Subtotal:	3	•		
	Total:	8			
	STS-131 (19A)	1	PWD Ambient	3/3/2010	4/22/2010
		1	PWD Hot		
		1	PWD Aux Port	2/25/2010	
		1	WPA RIP		
	Subtotal:	4	•		•
	STS-132 (ULF4)	2	PWD Ambient	4/26/2010, 5/18/2010	5/27/2010
		2	SVO-ZV	4/26/2010, 5/18/2010	
23		1	SRV-K Hot	5/18/2010	
		2	SRV-K Warm	4/26/2010, 5/18/2010	
	Subtotal:	7			
	Soyuz 21	1	PWD Hot	3/31/2010	6/3/2010
		1	SVO-ZV		
		1	SRV-K Hot		
	Subtotal:	3			-
	Total:	14			
	Soyuz 22	2	PWD Ambient	7/14/2010, 9/15/2010	9/26/2010
24		1	PWD Hot	8/25/2010	
44		1	WPA RIP	7/29/2010	
	Total:	4			
	Soyuz 23	2	PWD Hot	10/19/2010, 11/23/2010	11/29/10
25		1	SRV-K Warm	11/23/10	
	Total:	3			

The WRS includes a urine processor assembly (UPA) that uses a distillation process to purify pretreated urine. The resulting urine distillate is stored in a wastewater tank where it is combined with humidity condensate recovered from the ISS atmosphere. The WRS also includes a water processor assembly (WPA) that uses adsorption/ion exchange and thermal catalytic oxidation to remove contaminants from the combined wastewater feed stream. After iodine biocide is added, the WPA product water is then stored for delivery to the potable water bus. The potable water dispenser (PWD) receives water directly from the potable bus and provides either hot or ambient water for crew use, after removing the iodine biocide at the point of use. Iodinated potable water transferred from the Space Shuttle can also be added directly to the WPA storage tank in the event that make-up water is needed.

A total organic carbon analyzer (TOCA) was also delivered with the WRS in November of 2008. In Since that time the TOCA has served as a real-time tool for monitoring the total organic carbon (TOC) content of the WPA product water on a weekly basis. The TOCA includes a dedicated hose for direct sampling of water from the WPA product tank. It has also been used monthly to analyze samples collected from the PWD use-point in Teflon® sample bags. In-flight TOCA results provide insight into the total amount of organic contamination in the WPA water without identification of specific organic constituents.

Table 2. Water Analytical Methods

Parameter	Method
pH & conductivity	Potentiometric
Total Dissolved Solids	Gravimetric
Turbidity	Nephelometric
Iodine & iodide	Leuco crystal violet (LCV)
Fluoride	Ion chromatography (IC)
Metals/Minerals	Inductively coupled plasma/mass spectrometry (ICP/MS)
Inorganic anions & cations	Ion chromatography (IC)
Total organic carbon (TOC)	Ultraviolet or heated persulfate oxidation
Alcohols & glycols	Direct Injection gas chromatography/mass spectrometry (GC/MS)
Volatile organics	GC/MS with a purge & trap concentrator
Semi-volatile organics	GC/MS after liquid/liquid extraction
Organic acids & amines	Capillary electrophoresis (CE)
Urea/Caprolactam	Liquid chromatography (LC) with UV diode array detector
Formaldehyde	GC/MS after derivatization & extraction
Glycerol	LC/MS/MS
Dimethylsilanediol	Direct injection GC/MS or LC with refractive index detector

Russian Segment Water Systems

The Russian condensate water recovery system (SRV-K), which is located in the Service Module (SM), treats humidity condensate recovered from SM cabin air into potable water as previously described. ¹⁻³ Condensate from the U.S. Segment that has been stored in a Contingency Water Container (CWC) can also be transferred using a Condensate Feed Unit and processed by the SRV-K. After the condensate is filtered to remove particulate matter, it flows through a catalytic filter reactor, phase separator, and multifiltration beds to remove organic and inorganic contaminants. Prior to storage the processed water flows through a conditioning bed where silver biocide is added, as well as minerals (calcium, magnesium, and fluoride) to improve palatability. The conditioned product water is pumped from storage to the SRV-K galley where it is pasteurized then made available to the crews via two dispenser ports (hot and warm). Should the demand for potable water exceed the availability of condensate for processing, the crews can install a container of stored potable water to a connection located downstream of the conditioning bed and upstream of the galley for use as make-up water.

The Russian Segment stored potable water system or SVO-ZV is also located in the Service Module as previously described. ¹⁻³ The SVO-ZV provides the crew access to Russian ground-supplied potable water (Rodnik water) that has been launched on Russian Progress vehicles and delivered to the ISS. During pre-launch preparation, groundwater from Korolev, Russia is filtered and silver is added electrolytically as a biocide before storage on Progress in 210-liter Rodnik tanks.

Discussion of Analytical Results

Results from chemical analyses of the SRV-K (regenerated), SVO-ZV (stored), and U.S. WPA archival potable water samples collected during Expeditions 21 through 25 are summarized in Appendices 1, 2 and 3, respectively. Specification limits, if any, are included in these data tabulations for comparison with the analytical results. A discussion of findings obtained from the chemical analysis of archival samples collected during each expedition and their compliance with ISS requirements follows.

EXPEDITION 21

A total of 7 chemical archival potable water samples, including 1 SRV-K hot, 2 SVO-ZV, 2 PWD hot, and 2 PWD ambient, were collected during Expedition 21 as detailed in Table 1. The 4 samples collected on November 10, 2009 were returned on STS-129 (ULF3) and received in the WAFAL on November 30, 2009. The 3 samples collected on October 20, 2009 were returned on Soyuz 19 and received in the WAFAL on December 14, 2009. All of the samples were collected using U.S. 1-liter Teflon® water sample bags. The 4 samples that were returned on Shuttle all had sufficient sample volume to support full chemical characterization. Due to limited sample volume, turbidity was not analyzed on the PWD hot sample and solids were not analyzed on any of the 3 samples collected for Soyuz return.

ISS US SEGMENT:

WPA Processed Water Samples

All chemical parameters measured for the 4 PWD samples were within limits specified in Table LXX of SSP $41000.^{10}$ Nickel levels ranged from 17 to 33 ug/L (see Figure 1). The total iodine (I) was <0.05 mg/L and meets the point of consumption limit of 0.2 mg/L total I (see Figure 2). The iron levels have continued to be at low levels after the PWD disinfection flush on 3/23/09 (see Figure 3). The TOC values ranged from 0.14 to 0.30 mg/L and are well within the 3 mg/L limit (see Figure 4). Trace levels of methyl sulfone (34 to 111 μ g/L) and formaldehyde (<2 to 2 μ g/L) were identified; however, these levels were toxicologically insignificant.

ISS RUSSIAN SEGMENT:

SRV-K Potable Water Sample

All chemical parameters measured for the SRV-K potable water sample were within the requirements found in Table D-1 of the ISS MORD, Revision C. The nickel level of 34 μ g/L was well below the limit of 100 μ g/L (see Figure 5). The silver biocide level of 14 μ g/L was low, indicating that heating of the water by the pasteurization unit continues to be the main source of microbial control in the SRV-K galley. As shown in Figure 6, an updated plot of the TOC trend in the SRV-K water samples, the TOC level was 0.17 mg/L and well below the specification limit. Trace levels of acetone (4 μ g/L) chloroform (0.7 μ g/L), diethyl ether (4 μ g/L), acetaldehyde (3 μ g/L), dibutylphthalate (8 μ g/L), and N-phenyl-2-naphthylamine (4 μ g/L) were identified. Their levels, however, were toxicologically insignificant.

SVO-ZV Potable Water Samples

All chemical parameters measured for the 2 SVO-ZV water samples were within the ISS MORD requirements except for a turbidity level of 1.6 NTU in the October 20, 2009 sample, which slightly exceeded the 1.5 NTU limit. The concern with elevated turbidity in the SVO-ZV samples is that particulates causing the turbidity can shield bacteria from the silver biocide. The dissolved silver level of 81 μ g/L in the October 20 sample indicates that a small amount of the particulates are due to colloidal silver which may help to mitigate this concern. The total silver levels of 28 to 111 μ g/L continued to be lower than the expected range of 400 to 500 μ g/L, thereby increasing the risk of microbial growth in the water. Updated plots of manganese, and of turbidity, total and colloidal silver (Ag_{colloidal} = Ag_{total} - Ag_{dissolved}), and formate in SVO-ZV samples are presented in Figures 7 and 8. The TOC levels in the 2 SVO-ZV samples were 0.36 mg/L and 2.70 mg/L, well below the 20 mg/L limit. Low levels of chloroform (1.6 to 51.9 μ g/L) and a trace level of bromodichloromethane (1.9 μ g/L) were identified; however, these levels were not considered to be toxicologically significant.

EXPEDITION 22

A total of 8 chemical archival potable water samples, including 2 SRV-K hot, 2 SRV-K warm, 3 SVO-ZV, and 1 PWD ambient, were collected during Expedition 22 as detailed in Table 1. The 5 samples collected on January 6, 2010 and February 3, 2010 were returned on STS-130 (20A) and received in the WAFAL on February 22, 2010. The 3 samples collected on March 3, 2010 were returned on Soyuz 20 and received in the WAFAL on March 26, 2010. All of the samples were collected in U.S. 1-liter Teflon® water sample bags. The 5 samples that were returned on Shuttle all had sufficient sample volume to support full chemical characterization. Due to limited sample volume, solids were not analyzed on the 3 samples that were returned on Soyuz.

ISS US SEGMENT:


WPA Processed Water Sample


All chemical parameters measured for the PWD ambient sample were within limits specified in SSP 41000. The nickel level was 10 μ g/L (refer to Figure 1). The total iodine (I) level was <0.05 mg/L and met the point of consumption limit of 0.2 mg/L (Figure 2). The iron level continued to be at a low level after the PWD disinfection flush on 3/23/09 (Figure 3). The TOC value was 0.16 mg/L and is well within the 3 mg/L limit (see Figure 4). Trace levels of methyl sulfone (34 μ g/L) and formaldehyde (4 μ g/L) were identified; however, these levels were toxicologically insignificant.


ISS RUSSIAN SEGMENT:

SRV-K Potable Water Samples

All chemical parameters measured for the 4 SRV-K warm & hot potable water samples were within the ISS MORD requirements. The nickel levels (26-48 µg/L) were within specifications (Figure 5). The silver level (added










Figure 5 - Nickel Levels in SRV-K Water Samples ISS Flight 4A to Soyuz 23




Figure 6 - Total, Formate, & Non-formate Organic Carbon in SRV-K Potable Water ISS Flight 4A to Soyuz 23




Figure 7 - Manganese Levels in SVO-ZV Water Samples ISS Flight 5A to Soyuz 21





Figure 8 - Turbidity, Formate, Total & Colloidal Silver in SVO-ZV Water Samples



biocide) continued to be low in the SRV-K samples (2-90 μ g/L), indicating that heating of the water by the galley continues to be the main source of microbial control. The TOC levels ranged from 0.19-1.32 mg/L (Figure 6). Trace levels of bromodichloromethane (1.2 μ g/L), 1,2-dichloroethane (0.6 μ g/L), chloroform (28 μ g/L), diethyl ether (4 μ g/L), o-xylene (1.5 μ g/L), and dibutylphthalate (5 μ g/L) were identified; however, these levels were not considered to be toxicologically significant.

SVO-ZV Potable Water Samples

All chemical parameters measured for the 3 SVO-ZV water samples were within the ISS MORD requirements except for turbidity (2.0 - 3.8 NTU) which exceeded the 1.5 NTU limit. The concern with elevated turbidity in the SVO-ZV samples is that particulates causing the turbidity can shield bacteria from the silver biocide. The dissolved silver levels ranged from 2-90 μ g/L (Figure 8) indicating that a small amount of the particulates are due to colloidal silver which may mitigate this concern. The total silver levels in the January 6 and March 3 samples of 173 and 276 μ g/L, respectively, (Figure 8) were both below the expected range of 400-500 μ g/L, increasing the risk of microbial growth in the water. On the other hand, the total silver level in the February 3 sample was 405 μ g/L and within the expected range. Manganese levels in the 3 samples ranged from 35-36 μ g/L and were well below the ISS MORD limit (Figure 7). The TOC levels ranged from 3.26-3.48 mg/L, all well within specifications. Although chloroform (89.2 μ g/L) exceeded the 80 μ g/L EPA limit in the March 3rd sample, it was within the 6,500 μ g/L Spacecraft Water Exposure Guideline (SWEG) for total trihalomethanes and thus did not pose a crew health risk. Only trace levels of bromodichloromethane (3.0-3.1 μ g/L) and chloroform (68-77 μ g/L) were identified in the other 2 samples.

EXPEDITION 23

A total of 14 chemical archival potable water samples, including 2 SRV-K hot, 2 SRV-K warm, 3 SVO-ZV, 3 PWD ambient, 2 PWD hot, 1 PWD Auxiliary Port, and 1 WPA Rack Interface Panel (RIP) were collected during Expedition 23 as detailed in Table 1. The 4 samples collected on February 25, 2010 and March 3, 2010 were returned on STS-131 (19A) and received in the WAFAL on April 22, 2010. The 7 samples collected on April 26, 2010 and May 18, 2010 were returned on STS-132 (ULF4) and received in the WAFAL on May 27, 2010. The 3 samples collected on March 31, 2010 were returned on Soyuz 21 and received in the WAFAL on June 3, 2010. All of the samples were collected using U.S. 1-liter Teflon® water sample bags. The 5 samples that were returned on Shuttle all had sufficient sample volume to support full chemical characterization. Due to limited sample volume, solids were not analyzed on any of the 3 samples collected for Soyuz return.

ISS US SEGMENT:

WPA Processed Water Samples

All chemical parameters measured for the 7 potable water samples taken from the PWD, WPA RIP, and PWD Auxiliary Port met requirements in SSP 41000. The nickel levels (see Figure 1) ranged from 26-270 μ g/L, with the highest level coming from the Aux port. The total iodine (I) was <0.05 mg/L for the PWD ambient and hot samples and meets the requirement limit of 0.2 mg/L at the points of consumption (Figure 2). Total iodine levels were 2.24 and 2.72 mg/L in the PWD Aux and WPA RIP samples, both within the range of 2-4 mg/L that is the desired level for WPA product water. Iron was not detected in any of the 7 samples. The TOC levels ranged from 0.16-0.26 mg/L, well within the 3 mg/L limit (see Figure 4). Trace levels of acetone (11-39 μ g/L), iodomethane (5 μ g/L), methyl sulfone (34-98 μ g/L), o-xylene (1.5 μ g/L) and formaldehyde (2-6 μ g/L) were identified in the samples. These trace levels were toxicologically insignificant.

ISS RUSSIAN SEGMENT:

SRV-K Potable Water Samples

All chemical parameters measured for the 4 SRV-K water samples met requirements listed in the ISS MORD with the exception of the nickel level of 116 μ g/L in the March 31, 2010 SRV-K hot sample, which slightly exceeded the MORD specification of 100 μ g/L, but was well within the established SWEG of 300 μ g/L (Figure 5). The total silver level in the March 31 sample was 131 μ g/L, which is within the acceptable biocidal range (>100 μ g/L). The silver biocide levels returned to typical low levels in the SRV-K samples collected on April 26 and May18 (5-9 μ g/L), indicating that heating of the water in the galley was the main source of microbial control at that time. The TOC levels in the SRV-K samples ranged from 0.18-0.51 μ g/L. Traces of 1,4-dichlorobenzene (1.4 μ g/L), dibutylphthalate (6 μ g/L), 2-methylthiobenzothiazole (4 μ g/L), formaldehyde (2-6 μ g/L), and chloroform (1.7 μ g/L) were identified in the SRV-K samples, which were not toxicologically significant levels. The TOC levels in the 4 SRV-K water samples ranged from 0.18-0.51 μ g/L and were well below the limit (Figure 6).

SVO-ZV Potable Water Samples

All chemical parameters measured for the 3 SVO-ZV water samples met requirements listed in the ISS MORD. Manganese ranged from 27-35 μ g/L, well below the specification limit (Figure 7). The total silver levels of 132-187 μ g/L were within the acceptable biocidal range (>100 μ g/L), but on the lower end of that range which increases the risk of microbial growth (see Figure 8). The TOC levels in the SVO-ZV water samples ranged from 0.48-3.53 mg/L. The only organics identified in the samples were traces of bromodichloromethane (2.6 μ g/L), formaldehyde (3-4 μ g/L), and chloroform (1.4-47.8 μ g/L). Although the chloroform level of 47.8 μ g/L in the March 31, 2010 sample is higher than the historical average of 13.5 μ g/L in SVO-ZV samples, it falls well below the EPA maximum contaminant limit for trihalomethanes of 80 μ g/L. Accordingly, none of the organics detected were at levels sufficient to pose a crew health risk.


EXPEDITION 24

A total of 4 chemical archival potable water samples, including 2 PWD ambient, 1 PWD hot, and 1 WPA RIP were collected during Expedition 24 as detailed in Table 1. The 4 samples collected on July 14, 2010, July 29, 2010, August 25, 2010, and September 15, 2010 were returned on Soyuz 22 and received in the WAFAL on September 26, 2010. All of the samples were collected using U.S. 1-liter Teflon® water sample bags. Due to limited sample volume, turbidity was not analyzed on the WPA RIP sample, and solids were not analyzed on any of the samples.

ISS US SEGMENT:

PWD Potable Water Samples

All chemical parameters measured for the 3 potable water samples from the PWD met requirements in SSP 41000. The nickel levels (see Figure 1) ranged from 4-29 μ g/L. The total iodine (I) was <0.05 mg/L and met the requirement of 0.2 mg/L at points of consumption (Figure 2). The iron levels remained below the detection limit (Figure 3). The TOC levels ranged from 1.11-2.19 mg/L, within the 3 mg/L limit, but trending up beginning July 14, 2010 (see Figure 4). A comparison of in-flight versus archival TOC results is presented in Figure 9. The archival TOC results confirm the trend seen with the TOCA in-flight analyses, although the TOCA results, on average, are 45% lower than the archive results (Figure 9).



Routine analyses for target organics showed only trace levels of 2-butanone ($<2-4~\mu g/L$), methyl sulfone ($<4-36~\mu g/L$) and formaldehyde ($5-8~\mu g/L$) in the samples that were toxicologically insignificant. As the target list did not reveal the contaminant responsible for the increase in TOC, efforts to look for unknown compounds were initiated. A method for glycerol was developed but this compound was not detected in any of the potable samples. Further

work, however, showed an unknown peak in the gas chromatography/mass spectrometry (GC/MS) analysis for glycols. This unknown was identified as dimethylsilanediol (DMSD) by GC/MS analysis. ¹¹ A method for analysis of DMSD was then developed, which provided estimated DMSD levels from 4.7-8.2 mg/L. The results from the DMSD analysis should be considered estimates only as the synthesized standard used for calibration was not traceable. Using the estimated DMSD levels, the organic carbon accountability for the potable water samples ranged from 98-111% indicating that DMSD was the primary contributor to the elevated TOC in the WPA product water. An interim SWEG of 25 mg/L was subsequently established for DMSD. Accordingly, in-flight levels did not pose a crew health risk. ¹²

Non-routine silicon analysis by inductively coupled plasma/mass spectrometry (ICP/MS) was also conducted to confirm that the unknown contained silicon and whether there were additional silicon containing contaminants in addition to the DMSD. The results show that the silicon levels ranged from 1.21-2.13 mg/L. These levels are slightly less than values predicted using the estimated DMSD levels, confirming that most of the silicon was accounted for in the DMSD analysis.

More detailed descriptions of the extensive work performed to identify and quantify the DMSD and of the multidisciplinary teamwork used to evaluate associated risk and investigate potential sources and root cause can be found in other recent publications. ^{12,13}

WPA Processed Water Sample

The TOC level in the product water sample from the WPA RIP was 2.18 mg/L, which also reflects an upward trend (Figure 4). Methyl sulfone and formaldehyde were identified at toxicologically insignificant levels of 150 μ g/L and 8 μ g/L, respectively. Glycerol was not detected and the silicon level was 1.77 mg/L. The estimated DMSD concentration was 7.3 mg/L, and well below the interim SWEG of 25 mg/L. Organic carbon accountability for the sample improved to 90%, with DMSD as the primary contributor.

EXPEDITION 25

A total of 3 chemical archival potable water samples, including 2 PWD hot and 1 SRV-K warm, were collected during Expedition 25 as detailed in Table 1. The 3 samples were collected on October 19, 2010 and November 23, 2010, returned on Soyuz 23 and received in the WAFAL on November 29, 2010. All of the samples were collected using U.S. 1-liter Teflon® water sample bags. Due to limited sample volume, total dissolved solids were not analyzed on any of the samples and turbidity was not analyzed on the October 19 PWD sample or the SRV-K sample.

ISS US SEGMENT:

PWD Potable Water Samples

All chemical parameters measured for the 2 potable water samples from the PWD met ISS quality requirements listed in SSP 41000. The nickel levels of 17 and 23 µg/L were well within specifications (see Figure 1). The total iodine (I) levels were below the detection limit (<0.05 mg/L) and met the ISS limit of 0.2 mg/L at points of consumption (Figure 2). The iron levels remained below the detection limit (Figure 3). The TOC levels of 2.51 mg/L on October 19 and 0.15 mg/L on November 23 indicate that the TOC level in the WPA product water continued to rise until mid-October, but never reached the 3 mg/L limit before falling precipitously (see Figure 4). These results confirm the TOCA in-flight data trend and an offset between TOCA results and archive results (Figure 9). The primary contributor to the elevated TOC again proved to be DMSD. The DMSD level in the October 19 sample of 8.49 mg/L is the highest seen in product water to date, but still well below the recently established SWEG of 25 mg/L. As expected, DMSD was not detected (<2 mg/L) in the November 23 sample collected after the TOC of the WPA product water returned to a nominal low level. Again, these DMSD results should be considered estimates as the synthesized standard used for calibration of the DMSD method was not traceable. Routine analyses for other target organics detected only trace levels of acetone (< 2 and 40 µg/L), 2-butanone (< 2 and 24 µg/L), and methyl sulfone (48 and 50 µg/L) in the samples. Non-routine silicon analysis by ICP/MS was also performed to confirm the presence of DMSD and determine if there were additional silicon containing contaminants present in the samples. The results indicated silicon levels of 0.94 and 0.22 mg/L in the October 19 and November 23 samples, respectively. These levels are slightly less than values predicted using the estimated DMSD levels, thereby confirming that most of the silicon in the samples was being accounted for in the DMSD analysis.

ISS RUSSIAN SEGMENT:

SRV-K Potable Water Sample

All chemical parameters measured for the SRV-K warm water sample were within the ISS MORD requirements. The nickel level of 32 μ g/L was well below the specification limit (Figure 5). The silver biocide level of 75 μ g/L was below the acceptable biocidal range (>100 μ g/L) indicating that heating of the water by the pasteurization unit was the main source of microbial control in the SRV-K galley at the time. As shown in Figure 6, the TOC level was 0.38 mg/L and well below the specification limit. No specific organics were identified in the sample.

Conclusions

The chemical analysis results for the archival potable water samples collected and returned during Expeditions 21 through 25 indicate that the ISS potable water supplies complied with ISS quality specifications and were acceptable for crew consumption.

Turbidity levels exceeded the ISS MORD limit of 1.5 NTU in 4 of 8 SVO-ZV samples collected during Expeditions 21 through 25. Although elevated turbidity in itself does not pose a direct crew health risk, the concern is that particulate matter causing the turbidity could shield bacteria from the silver biocide that is added for microbial control. Analytical results for dissolved silver indicate that some of the particulates contributing to elevated turbidity are colloidal silver, which has biocidal properties and may help to mitigate the turbidity concern. The silver biocide levels in 2 of 8 SVO-ZV samples were below the acceptable biocidal range (>100 μ g/L), which indicates increased risk of microbial growth. Continued close monitoring of the SVO-ZV silver level is therefore recommended to determine if remedial action is required.

Silver biocide levels in all but one SRV-K sample remained typically low ($<100~\mu g/L$), indicating that heating of the water by the pasteurization unit continued to be the primary means of microbial control in the SRV-K galley. The nickel level of 116 $\mu g/L$ in the March 31, 2010 SRV-K hot sample slightly exceeded the ISS MORD specification of 100 $\mu g/L$, but was well within the established SWEG of 300 $\mu g/L$.

From early July to mid-October the TOC level in the WPA product water showed a continuous rise, approaching the 3 mg/L limit, before falling back to nominal low levels in mid-November. The analytical results for archival PWD and WPA product water samples confirm the TOCA in-flight data trend and an offset between TOCA results and archive results. DMSD was eventually identified as the primary contributor to the elevated TOC. The November 23 PWD sample, collected after the TOC of the WPA product water returned to a nominal low level, did not contain DMSD. It is recommended that WPA troubleshooting support continue as well as support for the multidisciplinary effort to establish root cause and the environmental source(s) of the DMSD.

Appendix

The chemical analysis results for the archival potable water samples collected from the Russian Segment SRV-K (regenerated water) system during Expeditions 21 through 25 are presented in Appendix 1. Analytical results for the chemical archive samples collected from the Russian Segment SVO-ZV (stored water) system during these 5 expeditions are presented in Appendix 2. Appendix 3 contains the chemical analysis results for U.S. Segment Water Processor Assembly product water samples collected during Expeditions 21-25.

Acknowledgments

All of the work described in this paper was performed at the JSC Water and Food Analytical Laboratory (WAFAL) under NASA contract NAS9-02078. The authors wish to acknowledge the ISS Expeditions 21-25 crewmembers for collecting and packing the chemical archival samples for return to the ground. Mickie Benoit of the WAFAL is acknowledged for coordinating the retrieval and delivery of the returned ISS samples to WAFAL. The following WAFAL chemists are gratefully acknowledged for performing chemical analyses of the chemical archival water samples: Jim Alverson, Lydia Ding, Mike Kuo, Esther Liu, Jeff Rutz, and Dawn Zapp. Finally, our Russian colleagues, Peter Andreichuk and Elena Zapryagaylo of RSC-Energia, Yuri Sinyak and Sergei Harin of the Institute of Biomedical Problems, and Leonid Bobe of NIICHIMMASH are also acknowledged.

References

- ¹Straub, J. E., Plumlee, D. K., and Schultz, J. R., "ISS Expeditions 16 thru 20: Chemical Analysis Results for Potable Water", AIAA-2010-6042, AIAA Proceedings of the 2010 International Conference on Environmental Systems, 2010.
- ²Straub, J. E., Plumlee, D. K., and Schultz, J. R., "Chemical Analysis Results for Potable Water Returned from ISS Expeditions 14 and 15", *SAE International Journal of Aerospace*. 1(1): 556-577, 2008.
- ³Straub, J. E., Plumlee, D. K., and Schultz, J. R., "Sampling and Chemical Analysis of Potable Water for ISS Expeditions 12 and 13". *SAE International Journal of Aerospace*, 2007.
- ⁴Straub, J. E., Plumlee, D. K., and Schultz, J. R., "ISS Expeditions 10 & 11 Potable Water Sampling and Chemical Analysis Results", *SAE International Journal of Aerospace*, 2006.
- ⁵Straub, J. E., Plumlee, D. K., and Schultz, J. R., "Chemical Analysis of ISS Potable Water from Expeditions 8 and 9", *SAE International Journal of Aerospace*, 2005.
- ⁶Straub, J. E., Plumlee, D. K., and Schultz, J. R., "ISS Potable Water Sampling and Chemical Analysis: Expeditions 6 & 7", SAE International Journal of Aerospace, 2004.
- ⁷Plumlee, D. K. and Schultz, J. R., "ISS Potable Water Sampling and Chemical Analysis: Expeditions 4 & 5", *SAE International Journal of Aerospace*, 2003.
- ⁸Plumlee, D. K. and Schultz, J. R., "Chemical Sampling and Analysis of ISS Potable Water: Expeditions 1-3", *SAE International Journal of Aerospace*, 2002.
- ⁹ISS Medical Operations Requirements Document, SSP 50260, Revision C, NASA Johnson Space Center, February 2006, Section 7.2 and Appendix D, Table D-1.
- ¹⁰System Specification for the International Space Station, SSP 41000 BP, National Aeronautics and Space Administration, September 28, 2009, Table LXX.
- ¹¹Bentley, N, et al., "Second-Generation International Space Station Total Organic Carbon Analyzer Verification Testing and On-Orbit Performance Results", AIAA-2010-6043, *AIAA Proceedings of the 2011 International Conference on Environmental Systems*, 2010.
- Systems, 2010.

 12 Rutz, J. A. et al., "Discovery and Identification of Dimethylsilanediol as a Contaminant in ISS Potable Water", AIAA Technical Paper #1021758. 41 threnational Conference on Environmental Systems, 2011.

 13 McCoy, J. T. et al., "The Story Behind the Numbers: Lessons Learned from the Integration of Monitoring Resources in
- ¹³McCoy, J. T. et al., "The Story Behind the Numbers: Lessons Learned from the Integration of Monitoring Resources in Addressing an ISS Water Quality Anomaly", AIAA Technical Paper #1020578 . 41st International Conference on Environmental Systems, 2011.

Appendix 1. ISS SRV-K Potable Water (Regenerated) Summary for Expeditions 21 through 25

Mission					ISS ULF3/Exp. 21	ISS 20A	/Exp. 22	Sovuz 20	0/Exp. 22	Soyuz 21/Exp. 23	IS	SS ULF4/Exp.	23	Soyuz 23/Exp. 25
Sample Location			Potable Water		SRV-K Hot	SRV-K Warm	SRV-K Hot	SRV-K Warm	SRV-K Hot	SRV-K Hot		SRV-K Warm	SRV-K Hot	SRV-K Warm
Sample Description Sample Date Analysis/Sample ID	Units	Test Conducted by	Maximum Contaminant Level (MCL)	Maximum Contaminant Level Source	Potable Water 11/10/2009 20091130011	Potable Water 1/6/2010 20100222006	Potable Water 2/3/2010 20100222008	Potable Water 3/3/2010 20100326002	Potable Water 3/3/2010 20100326003	Potable Water 3/31/2010 20100603001	Potable Water 4/26/2010 20100527010	Potable Water 5/18/2010 20100527011	5/18/2010	Potable Water 11/23/2010 20101129003
Physical Characteristics						1								
Physical Characteristics	77 **	TI C	7.7.0.0	MORD	7.01	7.00	7.00	7.45	7.04	0.00	0.71	0.40	0.05	r 00
DH Conductivity	pH unit		5.5-9.0	MORD	7.31	7.09	7.00	7.45	7.34	6.33	6.51	6.43	6.35	5.82
Turbidity	μS/cm NTU	U.S. U.S.	1.5*	MORD	15 0.4	14 0.2	189 0.4	49 0.5	50 0.2	187 0.5	0.2	0.3	21 0.4	124 NA
Total Dissolved Solids	mg/L	U.S.	100 (1,000*)	MORD	<5	15	112	NA	NA	NA	<5	NA	<5	NA NA
Total Dissolved Solids	mg/L	0.5.	100 (1,000)	Moltb	ν,	10	112	1471	1411	1471	ν,	14/1	ν,	101
Iodine (LCV)						1								
Total I	mg/L	U.S.	0.05	MORD	< 0.05	< 0.05	< 0.05	< 0.05	< 0.05	< 0.05	< 0.05	< 0.05	< 0.05	< 0.05
Anions (IC/ISE)														
Bromide	mg/L	U.S.			< 0.5	< 0.5	< 0.5	< 0.5	< 0.5	< 0.5	< 0.5	< 0.5	< 0.5	< 0.5
Chloride	mg/L	U.S.	250	MORD	< 0.15	< 0.15	4.47	1.09	1.14	4.89	0.21	0.18	0.19	4.37
Fluoride	mg/L	U.S.	1.5/4	MORD/EPA	<0.1	< 0.1	<0.1	<0.1	<0.1	<0.1	<0.1	< 0.1	<0.1	<0.1
Nitrate as Nitrogen (NO3-N)	mg/L	U.S.	10	MORD/EPA	< 0.11	< 0.11	0.16	< 0.11	< 0.11	< 0.11	<0.11	< 0.11	<0.11	0.38
Nitrite as Nitrogen (NO2-N)	mg/L	U.S.	1	EPA	NA	NA	NA	NA	NA	NA	NA	NA	NA	NA
Phosphate as P (PO4-P)	mg/L	U.S.			< 0.24	< 0.24	< 0.24	< 0.24	< 0.24	< 0.24	< 0.24	< 0.24	< 0.24	< 0.02
Sulfate	mg/L	U.S.	250	MORD	<0.75	< 0.75	13.5	0.94	1.02	19.2	< 0.75	< 0.75	< 0.75	11.6
G-H (GO)						1								
Cations (IC)	/T	U.S.	2/1	MORD/SWEG	<0.002	<0.002	< 0.002	<0.002	<0.002	0.018	<0.002	<0.002	< 0.002	<0.002
Ammonia as Nitrogen (NH3-N)	mg/L	U.S.	2/1	MORD/SWEG	<0.002 <0.002	<0.002	<0.002	<0.002 <0.002	<0.002	<0.002	<0.002	<0.002	<0.002	<0.002 <0.002
Lithium	mg/L	0.3.			<0.002	<0.002	<0.002	<0.002	<0.002	<0.002	<0.002	<0.002	<0.002	<0.002
Metals (ICP/MS)						+								
Calcium	mg/L	U.S.	100	MORD	2.42	2.34	28.0	9.36	9.61	21.1	3.23	3.17	3.17	18.1
Magnesium	mg/L	U.S.	50	MORD	0.03	0.04	4.21	0.16	0.16	5.34	0.17	0.21	0.21	4.14
Potassium	mg/L	U.S.		on	0.04	<0.01	0.72	0.01	<0.01	1.71	0.02	0.04	0.03	1.03
Sodium	mg/L	U.S.			0.01	0.02	2.12	0.01	0.01	3.88	0.13	0.15	0.15	5.12
Aluminum	μg/L	U.S.			<2	<2	19	3	4	53	6	<4	4	101
Antimony	μg/L	U.S.	6	EPA	<2	<2	<2	<2	<2	<2	<2	<4	<2	<4
Arsenic	μg/L	U.S.	10	MORD/EPA	<1	<1	<1	<1	<1	<1	<1	<2	<1	<2
Barium	μg/L	U.S.	1,000/10,000	MORD/SWEG	<1	<1	11	<1	<1	8	<1	<2	<1	12
Beryllium	μg/L	U.S.	4	EPA	<1	<1	<1	<1	<1	<1	<1	<2	<1	<2
Cadmium	μg/L	U.S.	5/22	MORD/SWEG	<1	<1	<1	<1	<1	1	<1	<2	<1	<2
Chromium	μg/L	U.S.	100	MORD/EPA	<5	<5	<5	<5	<5	<5	<5	<10	<5	<10
Cobalt	μg/L	U.S.			<1	<1	<1	<1	<1	<1	<1	<2	<1	<2
Copper	μg/L	U.S.	1,000/1,300	MORD/EPA	11	4	7	2	9	9	7	15	13	5
Iron	μg/L	U.S.	300	MORD	6	6	39	17	15	36	< 5	<10	<5	19
Lead	μg/L	U.S.	50/15	MORD/EPA	<1	<1	<1	<1	<1	1	<1	<2	<1	<2
Manganese	μg/L	U.S.	50/300	MORD/SWEG	<1	<1	13	<1	<1	28	1	<2	1	6
Mercury	μg/L	U.S.	2	MORD/EPA	< 0.5	< 0.5	< 0.5	< 0.5	<0.5	<0.5	< 0.5	<0.5	< 0.5	<1
Molybdenum	μg/L	U.S.	40	EPA HA MORD/SWEG	<1 34	<1 48	<1 26	<1 43	<1 39	<1	<1 36	<2 68	<1 46	<2 32
Nickel	μg/L	U.S. U.S.	100/300 10/50	MORD/SWEG MORD/EPA		48 <1		43 <1	39 <1	116 <1	36 <1	68 <2		32 <2
Selenium Silver	μg/L μg/L	U.S.	500/400	MORD/SWEG	<1 14	<1 7	<1 108	<1 7	3	<1 131	5	<2 7	<1 9	<2 75
Silver, Dissolved	μg/L μg/L	U.S.	300/400	WIORD/SWEG	5	2	90	<4	<4	58	<2	<4	3	57
Zinc Zinc	μg/L μg/L	U.S.	5,000/2,000	MORD/SWEG	59	69	27	27	21	54	34	129	42	18
Lanc	μ8/ L	0.0.	0,000/2,000	MORD/SWEG	33	33	~'	ω, r	~1	J1	34	160	- 16	10
Silicon (ICP/MS)	1					1	1		1			1	1	
Silicon (ICP/MS)	μg/L	U.S.			NA	NA	NA	NA	NA	NA	NA	NA	73	793
	1.2					1			1			<u> </u>		1
Total Organic Carbon (Sievers)						<u> </u>								
Total Inorganic Carbon	mg/L	U.S.			2.79	2.64	18.1	6.33	6.44	17.2	3.62	3.54	3.58	14.2

NA=Not analyzed; MI=Matrix interference

^{*}MORD limit 1.5 mg/L (Russian method)

^{**}limit does not include contribution from formate #TDS allowable limit after mineralization SWEG - 1000 days (5-2006)

Appendix 1. ISS SRV-K Potable Water (Regenerated) Summary for Expeditions 21 through 25

Mission					ISS ULF3/Exp. 21	ISS 20A	/Exp. 22	Sovuz 20)/Exp. 22	Soyuz 21/Exp. 23	IS	SS ULF4/Exp.	23	Soyuz 23/Exp. 25
Sample Location			Potable Water		SRV-K Hot	SRV-K Warm	SRV-K Hot	SRV-K Warm	SRV-K Hot	SRV-K Hot			SRV-K Hot	SRV-K Warm
			Maximum	Maximum										
Sample Description		Test	Contaminant	Contaminant	Potable Water	Potable Water	Potable Water	Potable Water	Potable Water	Potable Water	Potable Water	Potable Water	Potable Water	Potable Water
Sample Date		Conducted	Level	Level	11/10/2009	1/6/2010	2/3/2010	3/3/2010	3/3/2010	3/31/2010	4/26/2010	5/18/2010	5/18/2010	11/23/2010
Analysis/Sample ID	Units	by	(MCL)	Source	20091130011	20100222006	20100222008	20100326002	20100326003	20100603001	20100527010		20100527012	20101129003
Total Organic Carbon	mg/L	U.S.	20**	MORD	0.17	0.21	1.32	0.21	0.19	0.51	0.18	0.27	0.27	0.38
Volatile Organics														
Acetone	μg/L	U.S.	15,000	SWEG	4	<2	<2	<2	<2	<2	<2	<2	<2	<6
Acryloniltrile	μg/L	U.S.			<2	<2	<2	<2	<2	<2	<2	<2	<2	<6
Allyl chloride (3-Chloropropene)	μg/L	U.S.			<2	<2	<2	<2	<2	<2	<2	<2	<2	<6
Benzene	μg/L	U.S.	5	EPA	< 0.4	< 0.4	< 0.4	< 0.4	< 0.4	< 0.4	< 0.4	< 0.4	< 0.4	<1.2
Bromobenzene	μg/L	U.S.			< 0.4	< 0.4	< 0.4	< 0.4	< 0.4	< 0.4	< 0.4	< 0.4	< 0.4	<1.2
Bromochloromethane	μg/L	U.S.	90	EPA HA	<4	<4	<4	<4	<4	<4	<4	<4	<4	<12
Bromodichloromethane	μg/L	U.S.	THM 80	EPA	< 0.4	< 0.4	1.2	< 0.4	< 0.4	< 0.4	< 0.4	< 0.4	< 0.4	<1.2
Bromoform	μg/L	U.S.	THM 80	EPA	<2	<2	<2	<2	<2	<2	<2	<2	<2	<6
Bromomethane	μg/L	U.S.	10	EPA HA	<2	<2	<2	<2	<2	<2	<2	<2	<2	<6
2-Butanone (Methyl ethyl ketone)	μg/L	U.S.	4,000	EPA HA	<2	<2	<2	<2	<2	<2	<2	<2	<2	<6
n-Butylbenzene	μg/L	U.S.			< 0.4	< 0.4	< 0.4	< 0.4	< 0.4	< 0.4	< 0.4	< 0.4	< 0.4	<1.2
sec-Butylbenzene	μg/L	U.S.			< 0.4	< 0.4	< 0.4	< 0.4	< 0.4	< 0.4	< 0.4	< 0.4	< 0.4	<1.2
tert-Butylbenzene	μg/L	U.S.			<0.4	<0.4	< 0.4	<0.4	<0.4	<0.4	<0.4	<0.4	<0.4	<1.2
Carbon disulfide	μg/L	U.S.			<2	<2	<2	<2	<2	<2	<2	<2	<2	<6
Carbon tetrachloride	μg/L	U.S.	5	EPA	<0.4	<0.4	< 0.4	<0.4	<0.4	< 0.4	<0.4	< 0.4	<0.4	<1.2
Chloroacetonitrile	μg/L	U.S.			<10	<10	<10	<10	<10	<10	<10	<10	<10	<30
Chlorobenzene	μg/L	U.S.	100	EPA	<0.4	<0.4	<0.4	<0.4	<0.4	<0.4	<0.4	<0.4	<0.4	<1.2
1-Chlorobutane (Butyl chloride)	μg/L	U.S.			<0.4	< 0.4	<0.4	<0.4	<0.4	<0.4	< 0.4	< 0.4	<0.4	<1.2
Chloroethane	μg/L	U.S.		OTTER (ED.)	<2	<2	<2	<2	<2	<2	<2	<2	<2	<6
Chloroform	μg/L	U.S.	6,500/THM 80	SWEG/EPA	0.7	< 0.4	27.6	<0.4	<0.4	1.7	<0.4	<0.4	<0.4	<1.2
Chloromethane 2-Chlorotoluene	μg/L	U.S. U.S.	30	EPA HA	<2	<2	<2	<2 <0.4	<2	<2	<2 <0.4	<2	<2 <0.4	<6 <1.2
4-Chlorotoluene	μg/L	U.S.	100	EPA HA EPA HA	<0.4 <0.4	<0.4 <0.4	<0.4 <0.4	<0.4	<0.4 <0.4	<0.4 <0.4	<0.4	<0.4 <0.4	<0.4	<1.2
Dibromochloromethane	μg/L μg/L	U.S.	100 THM 80	EPA HA EPA	<0.4	<0.4	<0.4	<0.4	<0.4	<0.4	<0.4	<0.4	<0.4	<1.2
1.2-Dibromo-3-chloropropane (DBCP)	μg/L μg/L	U.S.	0.2	EPA	<0.4	<2	<0.4	<0.4	<0.4 <2	<0.4	<0.4	<0.4	<0.4	<1.2 <6
1,2-Dibromoethane (EDB)	μg/L μg/L	U.S.	0.05	EPA	<0.4	<0.4	<0.4	<0.4	<0.4	<0.4	<0.4	<0.4	<0.4	<1.2
Dibromomethane	μg/L μg/L	U.S.	0.03	ыл	<0.4	<0.4	<0.4	<0.4	<0.4	<0.4	<0.4	<0.4	<0.4	<1.2
1,2-Dichlorobenzene	μg/L	U.S.	600	EPA	<0.4	<0.4	<0.4	<0.4	<0.4	<0.4	<0.4	<0.4	<0.4	<1.2
1,3-Dichlorobenzene	μg/L	U.S.	600	EPA HA	<0.4	<0.4	<0.4	<0.4	<0.4	<0.4	<0.4	<0.4	<0.4	<1.2
1,4-Dichlorobenzene	μg/L	U.S.	75	EPA	<0.4	<0.4	1.3	<0.4	<0.4	<0.4	<0.4	1.4	<0.4	<1.2
trans-1,4-Dichloro-2-butene	μg/L	U.S.	-		<0.4	<0.4	< 0.4	<0.4	<0.4	<0.4	<0.4	<0.4	<0.4	<1.2
Dichlorodifluoromethane	μg/L	U.S.	1,000	EPA HA	<2	<2	<2	<2	<2	<2	<2	<2	<2	<6
1,1-Dichloroethane	μg/L	U.S.			< 0.4	< 0.4	< 0.4	< 0.4	< 0.4	< 0.4	< 0.4	< 0.4	< 0.4	<1.2
1,2-Dichloroethane	μg/L	U.S.	5	EPA	< 0.4	< 0.4	< 0.4	< 0.4	0.6	< 0.4	< 0.4	< 0.4	< 0.4	<1.2
1,1-Dichloroethene	μg/L	U.S.	7	EPA	< 0.4	< 0.4	< 0.4	< 0.4	< 0.4	< 0.4	< 0.4	< 0.4	< 0.4	<1.2
cis1,2-Dichloroethene	μg/L	U.S.	70	EPA	< 0.4	< 0.4	< 0.4	< 0.4	< 0.4	< 0.4	< 0.4	< 0.4	< 0.4	<1.2
trans-1,2-Dichloroethene	μg/L	U.S.	100	EPA	< 0.4	< 0.4	< 0.4	< 0.4	< 0.4	< 0.4	< 0.4	< 0.4	< 0.4	<1.2
1,2-Dichloropropane	μg/L	U.S.	5	EPA	< 0.4	< 0.4	< 0.4	< 0.4	< 0.4	< 0.4	< 0.4	< 0.4	< 0.4	<1.2
1,3-Dichloropropane	μg/L	U.S.			< 0.4	< 0.4	< 0.4	< 0.4	< 0.4	< 0.4	< 0.4	< 0.4	< 0.4	<1.2
2,2-Dichloropropane	μg/L	U.S.			< 0.4	< 0.4	< 0.4	< 0.4	< 0.4	< 0.4	< 0.4	< 0.4	< 0.4	<1.2
1,1-Dichloropropanone	μg/L	U.S.			<2	<2	<2	<2	<2	<2	<2	<2	<2	<6
1,1-Dichloropropene	μg/L	U.S.			< 0.4	< 0.4	< 0.4	< 0.4	< 0.4	< 0.4	< 0.4	< 0.4	< 0.4	<1.2
cis-1,3-Dichloropropene	μg/L	U.S.			<0.4	<0.4	< 0.4	<0.4	<0.4	<0.4	<0.4	<0.4	< 0.4	<1.2
trans-1,3-Dichloropropene	μg/L	U.S.			<2	<2	<2	<2	<2	<2	<2	<2	<2	<6
Diethyl ether	μg/L	U.S.			4	4	<2	<2	<2	<2	<2	<2	<2	<6
Ethylbenzene	μg/L	U.S.	700	EPA	<0.4	< 0.4	<0.4	<0.4	<0.4	<0.4	< 0.4	< 0.4	<0.4	<1.2
Ethyl methacrylate	μg/L	U.S.		TTD 4	<2	<2	<2	<2	<2	<2	<2	<2	<2	<6
Hexachlorobutadiene	μg/L	U.S.	1	EPA HA	<2	<2	<2	<2	<2	<2	<2	<2	<2	<6
Hexachloroethane	μg/L	U.S.	1	EPA HA	<2	<2	<2	<2	<2	<2	<2	<2	<2	<6
2-Hexanone	μg/L	U.S.			<2	<2	<2	<2	<2	<2	<2	<2	<2	<6

NA=Not analyzed; MI=Matrix interference

^{*}MORD limit 1.5 mg/L (Russian method)

^{**}limit does not include contribution from formate #TDS allowable limit after mineralization SWEG - 1000 days (5-2006)

Appendix 1. ISS SRV-K Potable Water (Regenerated) Summary for Expeditions 21 through 25

Mission					ISS ULF3/Exp. 21	ISS 20A	/Exp. 22	Sovuz 20	0/Exp. 22	Soyuz 21/Exp. 23	I IS	SS ULF4/Exp.	23	Soyuz 23/Exp. 25
Sample Location			Potable Water		SRV-K Hot	SRV-K Warm	SRV-K Hot	SRV-K Warm	SRV-K Hot	SRV-K Hot		SRV-K Warm	SRV-K Hot	SRV-K Warm
Sample Location					Sitt it flot	DICV IC VVIIII	Sit V It Hot	SKV K VVIIII	Sit V It Hot	SKV K HOT	Siev it wain	Sit V It Wallin	DRV R Hot	Sitv it vvaim
G 1 D 1 1			Maximum	Maximum	D . 11	D . 11 ***	n . 11 ***.	B - 11 *** -	D . 11	D . 11 *** .		n . 11	n . 11	D . 11
Sample Description		Test	Contaminant	Contaminant	Potable Water	Potable Water	Potable Water	Potable Water	Potable Water	Potable Water	Potable Water			Potable Water
Sample Date		Conducted	Level	Level	11/10/2009	1/6/2010	2/3/2010	3/3/2010	3/3/2010	3/31/2010	4/26/2010	5/18/2010	5/18/2010	11/23/2010
Analysis/Sample ID	Units	by	(MCL)	Source	20091130011	20100222006	20100222008	20100326002	20100326003	20100603001	20100527010	20100527011	20100527012	20101129003
Iodomethane	μg/L	U.S.			<2	<2	<2	<2	<2	<2	<2	<2	<2	<6
Isopropylbenzene (Cumene)	μg/L	U.S.	4,000	EPA DWEL	<0.4	<0.4	<0.4	<0.4	< 0.4	< 0.4	<0.4	< 0.4	<0.4	<1.2
4-Isopropyltoluene (Cymene)	μg/L	U.S.			<0.4	< 0.4	< 0.4	<0.4	< 0.4	< 0.4	< 0.4	< 0.4	< 0.4	<1.2
Methacrylonitrile	μg/L	U.S.			<2	<2	<2	<2	<2	<2	<2	<2	<2	<6
Methyl acrylate	μg/L	U.S.			<2	<2	<2	<2	<2	<2	<2	<2	<2	<6
Methyl-t-butylether (MTBE)	μg/L	U.S.			<2	<2	<2	<2	<2	<2	<2	<2	<2	<6
Methylene chloride (Dichloromethane)	μg/L	U.S.	15,000/5	SWEG/EPA	<0.4	<0.4	<0.4	<0.4	<0.4	< 0.4	<0.4	<0.4	<0.4	<1.2
Methyl methacrylate	μg/L	U.S.			<2	<2	<2	<2	<2	<2	<2	<2	<2	<6
4-Methyl-2-pentanone	μg/L	U.S.			<0.4	<0.4	< 0.4	<0.4	<0.4	<0.4	< 0.4	< 0.4	< 0.4	<1.2
Naphthalene	μg/L	U.S.	100	EPA HA	<0.4	<0.4	<0.4	<0.4	<0.4	<0.4	<0.4	<0.4	<0.4	<1.2
Nitrobenzene	μg/L	U.S.			<2	<2	<2	<2	<2	<2	<2	<2	<2	<6
2-Nitropropane	μg/L	U.S.			<2	<2	<2	<2	<2	<2	<2	<2	<2	<6
Pentachloroethane	μg/L	U.S.			<2	<2	<2	<2	<2	<2	<2	<2	<2	<6
Propionitrile (Ethyl cyanide)	μg/L	U.S.			<10	<10	<10	<10	<10	<10	<10	<10	<10	<30
n-Propylbenzene	μg/L	U.S.	100	770.4	<0.4	<0.4	< 0.4	<0.4	<0.4	<0.4	<0.4	<0.4	<0.4	<1.2
Styrene	μg/L	U.S.	100	EPA	<0.4	<0.4	<0.4	<0.4	<0.4	<0.4	<0.4	<0.4	< 0.4	<1.2
1,1,1,2-Tetrachloroethane	μg/L	U.S.	70	EPA HA	<0.4	<0.4	< 0.4	<0.4	<0.4	<0.4	<0.4	< 0.4	<0.4	<1.2
1,1,2,2-Tetrachloroethane	μg/L	U.S.	0.3	EPA HA	<0.4	<0.4	<0.4	<0.4	<0.4	<0.4	<0.4	< 0.4	< 0.4	<1.2
Tetrachloroethene	μg/L	U.S.	5	EPA	<0.4	<0.4	<0.4	<0.4	<0.4	<0.4	<0.4	<0.4	<0.4	<1.2
Tetrahydrofuran	μg/L	U.S.	1.000	TDA	<2	<2	<2	<2	<2	<2	<2	<2	<2	<6
Toluene	μg/L	U.S.	1,000	EPA	<0.4	<0.4	< 0.4	<0.4	<0.4	<0.4	<0.4	<0.4	<0.4	<1.2
1,2,3-Trichlorobenzene	μg/L	U.S.	mo.	TD 4	<0.4	<0.4	< 0.4	<0.4	<0.4	<0.4	<0.4	<0.4	<0.4	<1.2
1,2,4-Trichlorobenzene	μg/L	U.S.	70	EPA	<0.4	<0.4	< 0.4	<0.4	<0.4	<0.4	<0.4	<0.4	<0.4	<1.2
1,1,1-Trichloroethane	μg/L	U.S.	200	EPA	<0.4	<0.4	< 0.4	<0.4	<0.4	<0.4	<0.4	<0.4	<0.4	<1.2
1,1,2-Trichloroethane Trichloroethene	μg/L	U.S. U.S.	5 5	EPA EPA	<0.4 <0.4	<0.4 <0.4	<0.4 <0.4	<0.4 <0.4	<0.4 <0.4	<0.4 <0.4	<0.4 <0.4	<0.4 <0.4	<0.4 <0.4	<1.2 <1.2
Trichlorofluoromethane	μg/L	U.S.	2,000	EPA HA	<0.4	<0.4	<0.4	<0.4	<0.4	<0.4	<0.4	<0.4	<0.4	<1.2
	μg/L	U.S.	40	EPA HA	<0.4	<0.4	<0.4	<0.4	<0.4	<0.4	<0.4	<0.4	<0.4	<1.2
1,2,3-Trichloropropane 1,2,4-Trimethylbenzene	μg/L		40	ЕГА ПА	<0.4									<1.2
1,3,5-Trimethylbenzene	μg/L	U.S. U.S.			<0.4	<0.4 <0.4	<0.4 <0.4	<0.4 <0.4	<0.4 <0.4	<0.4 <0.4	<0.4 <0.4	<0.4 <0.4	<0.4 <0.4	<1.2
Vinyl Acetate	μg/L μg/L	U.S.			<0.4	<0.4	<0.4	<0.4	<2	<2	<2	<2	<0.4	<6
Vinyl Chloride	μg/L μg/L	U.S.	2	EPA	<2	<2	<2	<2	<2	<2	<2	<2	<2	<6
m&p-Xylene	μg/L μg/L	U.S.	~	EPA	<0.4	<0.4	<0.4	<0.4	<0.4	<0.4	<0.4	<0.4	<0.4	<1.2
o-Xylene	μg/L μg/L	U.S.	Total Xylenes 10,000 Total Xylenes 10,000	EPA	<0.4	1.5	<0.4	<0.4	<0.4	<0.4	<0.4	<0.4	<0.4	<1.2
0-Aylene	μg/L	0.3.	Total Xylenes 10,000	ErA	<0.4	1.3	<0.4	<0.4	<0.4	<0.4	<0.4	<0.4	<0.4	<1.2
Volatile Organics - Non-Targets (Tentatively Ident	fled Cor	npounds (>/=	80% match quality)										
Acetaldehyde	μg/L	U.S.	<u> </u>	,	3	not found	not found	not found	not found	not found	not found	not found	not found	not found
Butvraldehyde (Butanal)	μg/L	U.S.			not found	not found	not found	not found	not found	not found	not found	not found	not found	not found
1.8-Cineole	μg/L	U.S.			not found	not found	not found	not found	not found	not found	not found	not found	not found	not found
Cyclohexanone	μg/L	U.S.			not found	not found	not found	not found	not found	not found	not found	not found	not found	not found
Difluorodimethylsilane	μg/L				not found	not found	not found	not found	not found	not found	not found	not found	not found	not found
Dimethoxymethane (Formal)	μg/L	U.S.			not found	not found	not found	not found	not found	not found	not found	not found	not found	not found
2,5-Dimethylfuran	μg/L	U.S.			not found	not found	not found	not found	not found	not found	not found	not found	not found	not found
2,6-Dimethyl-1,7-octadiene	μg/L	U.S.			not found	not found	not found	not found	not found	not found	not found	not found	not found	not found
Dimethyl sulfide (Thiobismethane)	μg/L				not found	not found	not found	not found	not found	not found	not found	not found	not found	not found
p-Dioxane	μg/L	U.S.			not found	not found	not found	not found	not found	not found	not found	not found	not found	not found
1,3-Dioxolane (Ethylene glycol formal)	μg/L	U.S.			not found	not found	not found	not found	not found	not found	not found	not found	not found	not found
Ethyl acetate	μg/L	U.S.			not found	not found	not found	not found	not found	not found	not found	not found	not found	not found
Fluorotrimethylsilane	μg/L	U.S.			not found	not found	not found	not found	not found	not found	not found	not found	not found	not found
2-Heptanone	μg/L μg/L	U.S.			not found	not found	not found	not found	not found	not found	not found	not found	not found	not found
_														
3-Hexanone	μg/L	U.S.			not found	not found	not found	not found	not found	not found	not found	not found	not found	not found
3-Hydroxy-3-methylbutyric acid	μg/L	U.S.			not found	not found	not found	not found	not found	not found	not found	not found	not found	not found

NA=Not analyzed; MI=Matrix interference

^{*}MORD limit 1.5 mg/L (Russian method)

^{**}limit does not include contribution from formate #TDS allowable limit after mineralization SWEG - 1000 days (5-2006)

Appendix 1. ISS SRV-K Potable Water (Regenerated) Summary for Expeditions 21 through 25

Mission					ISS ULF3/Exp. 21	ISS 20A	/Exp. 22	Sovuz 20	0/Exp. 22	Soyuz 21/Exp. 23	IS	SS ULF4/Exp.	23	Soyuz 23/Exp. 25
Sample Location			Potable Water	1	SRV-K Hot	SRV-K Warm	SRV-K Hot	SRV-K Warm	SRV-K Hot	SRV-K Hot		SRV-K Warm	SRV-K Hot	SRV-K Warm
•			Maximum	Maximum										
Sample Description		Test	Contaminant	Contaminant	Potable Water									
Sample Description Sample Date		Conducted	Level	Level	11/10/2009	1/6/2010	2/3/2010	3/3/2010	3/3/2010	3/31/2010	4/26/2010	5/18/2010	5/18/2010	11/23/2010
Analysis/Sample ID	Units	by	(MCL)	Source	20091130011	20100222006	20100222008	20100326002	20100326003	20100603001	20100527010	20100527011	20100527012	20101129003
, ,			(MCL)	Source										
Isobutyronitrile	μg/L	U.S.			not found	not found not found	not found	not found	not found	not found not found				
Isooctanol	μg/L	U.S.			not found		not found	not found	not found					
L-Menthol	μg/L	U.S.			not found									
Menthone	μg/L	U.S.			not found									
Methyl acetate cis-1-Methyl-4-(1-methylethenyl)-cyclohexane	μg/L	U.S. U.S.			not found	not found	not found	not found not found	not found not found	not found not found	not found not found	not found	not found	not found
3-Methyl-2-pentanone	μg/L μg/L	U.S.			not found not found	not found not found	not found not found	not found	not found	not found	not found	not found not found	not found not found	not found not found
		U.S.			not found									
2-Methyl-1-propene alpha-Methyl styrene	μg/L	U.S.			not found	not found not found	not found	not found	not found	not found				
2-Nonanone	μg/L													
2-Octanone	μg/L	U.S. U.S.			not found not found									
2-Octanone 2-Pentanone	μg/L	U.S.			not found not found									
	μg/L	U.S.												
1,1,2-Trichloro-1,2,2-trifluoroethane (Freon 113) Trimethylsilanol	μg/L	U.S.			not found									
Trinieuryisiianoi	μg/L	U.S.			not found									
Extractable Organics														
Acetophenone	μg/L	U.S.			<8	<8	<8	<16	<16	<16	<8	<16	<8	<16
Benzaldehyde	μg/L μg/L	U.S.			<4	<4	<4	<8	<8	<8	<4	<8	<4	<8
Benzoic acid	μg/L μg/L	U.S.			<12	<12	<12	<24	<24	<24	<12	<24	<12	<48
Benzothiazole	μg/L μg/L	U.S.			<4	<4	<4	<8	<8	<8	<4	<8	<4	<8
Benzyl alcohol	μg/L μg/L	U.S.			<4	<4	<4	<8	<8	<8	<4	<8	<4	<8
Benzyl actorol Benzyl butyl phthlate	μg/L μg/L	U.S.	7.000	EPA DWEL	<4	<4	<4	<8	<8	<8	<4	<8	<4	<8
2-Butoxyethanol	μg/L μg/L	U.S.	7,000	EFA DVVEL	<8	<8	<8	<16	<16	<16	<8	<16	<8	<16
2-(2-Butoxyethanol	μg/L μg/L	U.S.			<8	<8	<8	<16	<16	<16	<8	<16	<8	<16
2-(2-Butoxyethoxy)ethalol 2-(2-Butoxyethoxy)ethyl acetate	μg/L μg/L	U.S.			<4	<4	<4	<8	<8	<8	<4	<8	<4	<8
n-Butylpalmitate	μg/L μg/L	U.S.			<8	<8	<8	<16	<16	<16	<8	<16	<8	<16
Butylated hydroxyanisole (BHA)	μg/L μg/L	U.S.			<4	<4	<4	<8	<8	<8	<4	<8	<4	<8
N-Butylbenzenesulfonamide	μg/L μg/L	U.S.			<4	<4	<4	<8	<8	<8	<4	<8	<4	<8
3-tert-Butylphenol	μg/L μg/L	U.S.			<12	<12	<12	<24	<24	<24	<12	<24	<12	<24
Caffeine	μg/L	U.S.			<4	<4	<4	<8	<8	<8	<4	<8	<4	<8
tris-2-Chloroethyl phosphate	μg/L	U.S.			<4	<4	<4	<8	<8	<8	<4	<8	<4	<8
Cholesterol	μg/L	U.S.			<32	<32	<32	<64	<64	<64	<32	<64	<32	<64
o-Cresol (2-Methylphenol)	μg/L	U.S.			<4	<4	<4	<8	<8	<8	<4	<8	<4	<8
Cyclododecane	μg/L μg/L	U.S.			<4	<4	<4	<8	<8	<8	<4	<8	<4	<8
Decamethylcyclopentasiloxane	μg/L	U.S.			<4	<4	<4	<8	<8	<8	<4	<8	<4	<8
Decanoic acid	μg/L	U.S.			<8	<8	<8	<16	<16	<16	<8	<16	<8	<24
2,6-Di-t-butyl-1,4-benzoquinone	μg/L	U.S.			<4	<4	<4	<8	<8	<8	<4	<8	<4	<8
2,4-Di-t-butylphenol	μg/L	U.S.			<4	<4	<4	<8	<8	<8	<4	<8	<4	<8
1,4 Diacetylbenzene	μg/L	U.S.			<4	<4	<4	<8	<8	<8	<4	<8	<4	<8
N,N-Dibutylformamide	μg/L	U.S.			<4	<4	<4	<8	<8	<8	<4	<8	<4	<8
Dibutyl phthalate	μg/L	U.S.	40,000/4,000	SWEG/EPA DWEL	8	<4	5	<8	<8	<8	<4	<8	6	<8
Dibutylamine	μg/L	U.S.	Dialkylamines 300	SWEG	<4	<4	<4	<8	<8	<8	<4	<8	<4	<8
N,N-Diethyl-m-toluamide	μg/L	U.S.	, j	— —	<4	<4	<4	<8	<8	<8	<4	<8	<4	<8
Diethylphthalate	μg/L	U.S.	30,000	EPA DWEL	<4	<4	<4	<8	<8	<8	<4	<8	<4	<8
Diethylene glycol monoethyl ether	μg/L	U.S.	,		<4	<4	<4	<8	<8	<8	<4	<8	<4	<8
N,N-Diethylformamide	μg/L	U.S.			<12	<12	<12	<24	<24	<24	<12	<24	<12	<24
Diiodomethane (Methyl iodide)	μg/L	U.S.			<4	<4	<4	<8	<8	<8	<4	<8	<4	<8
Diisopropyl adipate	μg/L	U.S.			<4	<4	<4	<8	<8	<8	<4	<8	<4	<8
Dimethyl phthalate	μg/L	U.S.			<4	<4	<4	<8	<8	<8	<4	<8	<4	<8
N,N-Dimethyl acetamide	μg/L	U.S.			<4	<4	<4	<8	<8	<8	<4	<8	<4	<8
N,N-Dimethylbenzylamine	μg/L	U.S.	Dialkylamines 300	SWEG	<4	<4	<4	<8	<8	<8	<4	<8	<4	<8
N,N-Dimethylformamide	μg/L	U.S.	J		<8	<8	<8	<16	<16	<16	<8	<16	<8	<16
Dipropylene glycol methyl ether	μg/L	U.S.			<4	<4	<4	<8	<8	<8	<4	<8	<4	<8
1 10	L F-0'~			Į.	•									

NA=Not analyzed; MI=Matrix interference

^{*}MORD limit 1.5 mg/L (Russian method)

^{**}limit does not include contribution from formate #TDS allowable limit after mineralization SWEG - 1000 days (5-2006)

Appendix 1. ISS SRV-K Potable Water (Regenerated) Summary for Expeditions 21 through 25

Mission					ISS ULF3/Exp. 21	ISS 20A	/Exp. 22	Sovuz 20	D/Exp. 22	Soyuz 21/Exp. 23	15	SS ULF4/Exp.	23	Sovuz 23/Exp. 25
Sample Location			Potable Water		SRV-K Hot	SRV-K Warm	SRV-K Hot	SRV-K Warm	SRV-K Hot	SRV-K Hot		SRV-K Warm		SRV-K Warm
			Maximum	Maximum										
Sample Description		Test	Contaminant	Contaminant	Potable Water	Potable Water	Potable Water	Potable Water	Potable Water	Potable Water	Potable Water	Potable Water	Potable Water	Potable Water
Sample Date		Conducted	Level	Level	11/10/2009	1/6/2010	2/3/2010	3/3/2010	3/3/2010	3/31/2010	4/26/2010	5/18/2010	5/18/2010	11/23/2010
Analysis/Sample ID	Units	by	(MCL)	Source	20091130011	20100222006	20100222008	20100326002	20100326003	20100603001	20100527010	20100527011	20100527012	20101129003
Dodecamethylcyclohexasiloxane	μg/L	U.S.	(MOL)	Dource	<4	<4	<4	<8	<8	<8	<4	<8	<4	<8
2-Ethoxyethanol	μg/L μg/L	U.S.			<4	<4	<4	<8	<8	<8	<4	<8	<4	<16
2-Ethyl-1-hexanol	μg/L μg/L	U.S.			<4	<4	<4	<8	<8	<8	<4	<8	<4	<8
2-Ethylhexanoic acid	μg/L	U.S.			<4	<4	<4	<8	<8	<8	<4	<8	<4	<16
bis-2-Ethylhexyl adipate	μg/L	U.S.	400	EPA	<4	<4	<4	<8	<8	<8	<4	<8	<4	<8
bis-2-Ethylhexyl phthalate (Dioctyl phthlate)	μg/L	U.S.	20.000/6	SWEG/EPA	<4	<4	<4	<8	<8	<8	<4	<8	<4	<8
4-Ethylmorpholine	μg/L	U.S.	20,000/0	SVIEG/EIT	<4	<4	<4	<8	<8	<8	<4	<8	<4	<8
1-Formylpiperidine	μg/L	U.S.			<4	<4	<4	<8	<8	<8	<4	<8	<4	<8
Heptanoic acid	μg/L	U.S.			<4	<4	<4	<8	<8	<8	<4	<8	<4	<24
2-Heptanone	μg/L	U.S.			<4	<4	<4	<8	<8	<8	<4	<8	<4	<8
gamma-Hexalactone	μg/L	U.S.			<4	<4	<4	<8	<8	<8	<4	<8	<4	<8
Hexanoic acid	μg/L	U.S.			<8	<8	<8	<16	<16	<16	<8	<16	<8	<24
2-Hexanol	μg/L	U.S.			<4	<4	<4	<8	<8	<8	<4	<8	<4	<8
2-Hydroxybenzothiazole	μg/L	U.S.			<4	<4	<4	<8	<8	<8	<4	<8	<4	<8
Ibuprofen	μg/L	U.S.			<4	<4	<4	<8	<8	<8	<4	<8	<4	<24
Iodoform	μg/L	U.S.			<4	<4	<4	<8	<8	<8	<4	<8	<4	<8
Isophorone	μg/L	U.S.	100	EPA HA	<4	<4	<4	<8	<8	<8	<4	<8	<4	<8
4-Isopropylphenol	μg/L	U.S.			<4	<4	<4	<8	<8	<8	<4	<8	<4	<8
Lauramide	μg/L	U.S.			<4	<4	<4	<8	<8	<8	<4	<8	<4	<8
Lauric acid (Dodecanoic acid)	μg/L	U.S.			<120	<120	<120	<240	<240	<240	<120	<240	<120	<240
p-Menth-1-en-8-ol (alpha-Terpineol)	μg/L	U.S.			<4	<4	<4	<8	<8	<8	<4	<8	<4	<8
2-Mercaptobenzothiazole	μg/L	U.S.	30,000	SWEG	<40	<40	<40	<80	<80	<80	<40	<80	<40	<80
2-Methyl-2,4-pentanediol	μg/L	U.S.	·		<4	<4	<4	<8	<8	<8	<4	<8	<4	<8
1-Methyl-2-pyrrolidinone	μg/L	U.S.			<4	<4	<4	<8	<8	<8	<4	<8	<4	<8
Methyl-4-hydroxybenzoate	μg/L	U.S.			<4	<4	<4	<8	<8	<8	<4	<8	<4	<8
Methyl sulfone	μg/L	U.S.			<4	<4	<4	<8	<8	<8	<4	<8	<4	<8
2-Methyl butyric acid	μg/L	U.S.			<12	<12	<12	<24	<24	<24	<12	<24	<12	<24
2-Methylthiobenzothiazole	μg/L	U.S.			<4	<4	<4	<8	<8	<8	<4	<8	4	<8
Monomethyl phthalate	μg/L	U.S.			<4	<4	<4	<8	<8	<8	<4	<8	<4	<8
Myristic acid	μg/L	U.S.			<24	<24	<24	<48	<48	<48	<24	<48	<24	<64
(+)-Neomenthol	μg/L	U.S.			<4	<4	<4	<8	<8	<8	<4	<8	<4	<8
Nicotine	μg/L	U.S.			<4	<4	<4	<8	<8	<8	<4	<8	<4	<8
Nonadecane	μg/L	U.S.			<4	<4	<4	<8	<8	<8	<4	<8	<4	<8
Nonanoic acid	μg/L	U.S.			<12	<12	<12	<24	<24	<24	<12	<24	<12	<24
1-Octadecanol	μg/L	U.S.			<12	<12	<12	<24	<24	<24	<12	<24	<12	<24
Octamethylcyclotetrasiloxane	μg/L	U.S.			<4	<4	<4	<8	<8	<8	<4	<8	<4	<8
Octanoic acid	μg/L	U.S.			<8	<8	<8	<16	<16	<16	<8	<16	<8	<48
4-tert-Octylphenol	μg/L	U.S.			<4	<4	<4	<8	<8	<8	<4	<8	<4	<8
Oleic acid	μg/L	U.S.			<40	<40	<40	<80	<80	<80	<40	<80	<40	<80
Oxindole	μg/L	U.S.			<4	<4	<4	<8	<8	<8	<4	<8	<4	<8
Palmitic acid	μg/L	U.S.			<120	<120	<120	<240	<240	<240	<120	<240	<120	<240
Palmitoleic acid	μg/L	U.S.			NA	NA	NA	NA	NA	NA	NA	NA	NA	NA
Pentacosane	μg/L	U.S.			<4	<4	<4	<8	<8	<8	<4	<8	<4	<8
sec-Phenethyl alcohol	μg/L	U.S.			<4	<4	<4	<8	<8	<8	<4	<8	<4	<8
Phenol	μg/L	U.S.	1,000/4,000	MORD/SWEG	<4	<4	<4	<8	<8	<8	<4	<8	<4	<8
2-Phenoxyethanol	μg/L	U.S.			<4	<4	<4	<8	<8	<8	<4	<8	<4	<8
N-Phenyl-2-naphthylamine	μg/L	U.S.	260,000	SWEG	4	<4	<4	<8	<8	<8	<4	<8	<4	<8
2-Phenyl-2-propanol	μg/L	U.S.			<4	<4	<4	<8	<8	<8	<4	<8	<4	<8
2-Phenylacetic acid	μg/L	U.S.			<16	<16	<16	<32	<32	<32	<16	<32	<16	<32
Phenethyl alcohol	μg/L	U.S.			<4	<4	<4	<8	<8	<8	<4	<8	<4	<8
2-Phenylphenol	μg/L	U.S.			<4	<4	<4	<8	<8	<8	<4	<8	<4	<8
Salicyclic Acid	μg/L	U.S.			<32	<32	<32	<64	<64	<64	<32	<64	<32	<64
trans-Squalene	μg/L	U.S.			<8	<8	<8	<16	<16	<8	<8	<16	<8	<16
****					+	·				· · · · · · · · · · · · · · · · · · ·	· · · · · · · · · · · · · · · · · · ·			

NA=Not analyzed; MI=Matrix interference

^{*}MORD limit 1.5 mg/L (Russian method)

^{**}limit does not include contribution from formate #TDS allowable limit after mineralization SWEG - 1000 days (5-2006)

Appendix 1. ISS SRV-K Potable Water (Regenerated) Summary for Expeditions 21 through 25

Mission	1				ISS ULF3/Exp. 21	ISS 20A	/Exp. 22	Sovuz 20	0/Exp. 22	Sovuz 21/Exp. 23	IS	S ULF4/Exp.	23	Soyuz 23/Exp. 25
Sample Location			Potable Water		SRV-K Hot	SRV-K Warm	SRV-K Hot	SRV-K Warm	SRV-K Hot	SRV-K Hot		SRV-K Warm	SRV-K Hot	SRV-K Warm
bampic Bookton					510 11 1100	511 11 1141111	510 11 1100	5111 11 1111111	511 1 1101	5111 11101	514 11 VIII.II	510 7 11 7741111	5111 11101	
Samuela Danastastan		TT4	Maximum Contaminant	Maximum	D . 11 11/	D . 11 377 .	D . 11 11/	D . 11 W.	D . 11 327 .	D . 11 W.	D . 11 11/	D . 11 117 .	D . 11 377 .	D . 11 117 .
Sample Description		Test		Contaminant	Potable Water	Potable Water	Potable Water	Potable Water	Potable Water	Potable Water	Potable Water	Potable Water 5/18/2010	Potable Water 5/18/2010	Potable Water 11/23/2010
Sample Date		Conducted	Level	Level	11/10/2009	1/6/2010	2/3/2010	3/3/2010	3/3/2010	3/31/2010	4/26/2010			
Analysis/Sample ID	Units	by	(MCL)	Source	20091130011	20100222006	20100222008	20100326002	20100326003	20100603001	20100527010	20100527011	20100527012	20101129003
Stearic acid	μg/L	U.S.			NA	NA	NA	NA	NA	NA	NA	NA	NA	NA
1-Tetradecanol	μg/L	U.S.			<4	<4	<4	<8	<8	<8	<4	<8	<4	<8
Tetramethylsuccinonitrile	μg/L	U.S.			<4	<4	<4	<8	<8	<8	<4	<8	<4	<8
Tetramethyl thiourea	μg/L	U.S.			<4	<4	<4	<8	<8	<8	<4	<8	<4	<8
Tetramethylurea	μg/L	U.S.			<4	<4	<4	<8	<8	<8	<4	<8	<4	<8
Thymol	μg/L	U.S.			<4	<4	<4	<8	<8	<8	<4	<8	<4	<8
1,3,5-Triallyl-1,3,5-triazine-2,4,6(1H,3H,5H)-trione	μg/L	U.S.			<4	<4	<4	<8	<8	<8	<4	<8	<4	<8
Tributylamine	μg/L	U.S.	Trialkylamines 400	SWEG	<4	<4	<4	<8	<8	<8	<4	<8	<4	<8
Tributyl phosphate	μg/L	U.S.			<4	<4	<4	<8	<8	<8	<4	<8	<4	<8
Triethyl phosphate	μg/L	U.S.			<8	<8	<8	<16	<16	<16	<8	<16	<8	<16
2,2,4-Trimethyl-1,3-pentanediol diisobutyrate	μg/L	U.S.			<8	<8	<8	<16	<16	<16	<8	<16	<8	<16
Tripropylene glycol monomethyl ether	μg/L	U.S.			<4	<4	<4	<8	<8	<8	<4	<8	<4	<8
Undecanoic acid	μg/L	U.S.			<24	<24	<24	<48	<48	<48	<24	<48	<24	<48
2-Undecanone	μg/L	U.S.			<4	<4	<4	<8	<8	<8	<4	<8	<4	<8
Valeric acid (Pentanoic acid)	μg/L	U.S.			<24	<24	<24	<48	<48	<48	<24	<48	<24	<48
Vanillin	μg/L	U.S.			<8	<8	<8	<16	<16	<16	<8	<16	<8	<16
Alcohols (DAI/GC/MS)														1
1-Butanol	μg/L	U.S.			<100	<100	<100	<100	<100	<100	<100	<100	<100	<100
2-Butanol	μg/L	U.S.			<100	<100	<100	<100	<100	<100	<100	<100	<100	<100
Ethanol	μg/L	U.S.			<100	<100	<100	<100	<100	<100	<100	<100	<100	<100
Methanol	μg/L	U.S.			<100	<100	<100	<100	<100	<100	<100	<100	<100	<100
2-Methyl-1-butanol	μg/L	U.S.			<100	<100	<100	<100	<100	<100	<100	<100	<100	<100
2-Methyl-2-butanol	μg/L	U.S.			<100	<100	<100	<100	<100	<100	<100	<100	<100	<100
3-Methyl-1-butanol (Isopentanol)	μg/L	U.S.			<100	<100	<100	<100	<100	<100	<100	<100	<100	<100
2-Methyl-1-propanol	μg/L	U.S.			<100	<100	<100	<100	<100	<100	<100	<100	<100	<100
2-Methyl-2-propanol	μg/L	U.S.			<100	<100	<100	<100	<100	<100	<100	<100	<100	<100
1-Pentanol (Amyl alcohol)	μg/L	U.S.			<100	<100	<100	<100	<100	<100	<100	<100	<100	<100
2-Pentanol (sec-Amyl alcohol)	μg/L	U.S.			<100	<100	<100	<100	<100	<100	<100	<100	<100	<100
3-Pentanol	μg/L	U.S.			<100	<100	<100	<100	<100	<100	<100	<100	<100	<100
1-Propanol	μg/L	U.S.			<100	<100	<100	<100	<100	<100	<100	<100	<100	<100
2-Propanol (Isopropanol)	μg/L	U.S.			<100	<100	<100	<100	<100	<100	<100	<100	<100	<100
Chronic (DAT/CC/MS)						-	-				-			
Glycols (DAI/GC/MS)	./*	TI C	10000/14000	MODD/EDA ***	1000	1000	1000	1000	1000	1000	1000	1000	1000	1000
1,2-Ethanediol (Ethylene glycol)	μg/L	U.S.	12000/14000	MORD/EPA HA	<1000	<1000	<1000	<1000	<1000	<1000	<1000	<1000	<1000	<1000
1,2-Propanediol (Propylene glycol)	μg/L	U.S.			<500	<500	<500	<500	<500	<500	< 500	<500	<500	<500
Glycerol (LC/MS/MS)	<u> </u>					-								
	/T	U.S.			NA	NA	NA	NA	NA	NA	NA	NA	NA	<300
Glycerol (1,2,3-Propanetriol)	μg/L	U.S.			NA.	NA	INA	INA	INA	NA	INA	NA	NA	<300
Silanols (LC/MS/MS) (R & D Method -NIST traces	hla eto=	dard not are:	lahla)											
Dimethylsilanediol (DMSD)	μg/L	U.S.	uudic)		NA	NA	NA	NA	NA	NA	NA	NA	NA	<400
Dimensionaliculor (DWSD)	μg/L	0.3.			IVA	IVA	INA	INA	INA	INA	INA	INA	IVA	<400
Carboxylates (CE)						 	 							
Acetate	μg/L	U.S.			<125	<125	<125	<125	<125	<125	<125	<125	<125	<125
Formate	μg/L μg/L	U.S.	2.500.000	SWEG	<125	<125	<125	<125	<125	<125	<125	<125	<125	<125
Glycolate	μg/L μg/L	U.S.	۵,300,000	SWEG	<125	<125	<125	<125	<125	<125	<125	<125	<125	<125
Glyoxylate	μg/L μg/L	U.S.			<125	<125	<125	<125	<125	<125	<125	<125	<125	<125
Lactate	μg/L μg/L	U.S.			<1000	<1000	<1000	<1000	<1000	<1000	<1000	<1000	<1000	<1000
Oxalate	μg/L μg/L	U.S.			<125	<125	MI	<125	<125	<125	<125	<125	<125	<125
Propionate	μg/L μg/L	U.S.			<125	<125	<125	<125	<125	<125	<125	<125	<125	<125
	rs-	0.0.			120	-120	120	1200	1220	1200	1220	1220	1200	-120
L.	!	1			l	1	1		1		1			

NA=Not analyzed; MI=Matrix interference

^{*}MORD limit 1.5 mg/L (Russian method)

^{**}limit does not include contribution from formate #TDS allowable limit after mineralization SWEG - 1000 days (5-2006)

Appendix 1. ISS SRV-K Potable Water (Regenerated) Summary for Expeditions 21 through 25

No.														
Mission					ISS ULF3/Exp. 21	ISS 20A	/Exp. 22	Soyuz 20)/Exp. 22	Soyuz 21/Exp. 23	IS	S ULF4/Exp.	23	Soyuz 23/Exp. 25
Sample Location			Potable Water		SRV-K Hot	SRV-K Warm	SRV-K Hot	SRV-K Warm	SRV-K Hot	SRV-K Hot	SRV-K Warm	SRV-K Warm	SRV-K Hot	SRV-K Warm
			Maximum	Maximum										
Sample Description		Test	Contaminant	Contaminant	Potable Water	Potable Water	Potable Water	Potable Water	Potable Water	Potable Water	Potable Water	Potable Water	Potable Water	Potable Water
Sample Date		Conducted	Level	Level	11/10/2009	1/6/2010	2/3/2010	3/3/2010	3/3/2010	3/31/2010	4/26/2010	5/18/2010	5/18/2010	11/23/2010
Analysis/Sample ID	Units	by	(MCL)	Source	20091130011	20100222006	20100222008	20100326002	20100326003	20100603001	20100527010	20100527011	20100527012	20101129003
Aldehydes														
Formaldehyde	μg/L	U.S.	12,000/1,000	SWEG/EPA HA	<2	<2	<2	<2	<2	6	3	2	5	<5
Amines (CE)														
Ethylamine	μg/L	U.S.	Monoalkylamines 2000	SWEG	<125	<125	<125	<125	<125	<125	<125	<125	<125	<125
Methylamine	μg/L	U.S.	Monoalkylamines 2000	SWEG	<125	<125	<125	<125	<125	<125	<125	<125	<125	<125
n-Propylamine	μg/L	U.S.	Monoalkylamines 2000	SWEG	<125	<125	<125	<125	<125	<125	<125	<125	<125	<125
Trimethylamine	μg/L	U.S.	Trialkylamines 400	SWEG	<125	<125	<125	<125	<125	<125	<125	<125	<125	<125
Non-volatiles (LC/UV-VIS)														
Urea	μg/L	U.S.			<800	<800	<800	<800	<800	<800	<800	<800	<800	<800
Caprolactam	μg/L	U.S.	100,000	SWEG	<4	<4	<4	<8	<8	<8	<4	<8	<8	<16
Organic Carbon Recovery	percent	U.S.			9.26	1.92	0.53	0.00	0.08	0.51	0.68	0.56	3.08	0.00
Unaccounted Organic Carbon	mg/L	U.S.			0.15	0.20	1.31	0.21	0.19	0.50	0.17	0.27	0.26	0.38

Mission	1				Sovue 10/Eve 21	ISS ULF3/Exp. 21	100 204	/Exp. 22	Sover 20/Evn 22	Soyuz 21/Exp. 23	ICCIII	4/Exp. 23
TATE OF THE PARTY					SVO-ZV	SVO-ZV	SVO-ZV	SVO-ZV	SVO-ZV	SVO-ZV	SVO-ZV	SVO-ZV
Sample Location Sample Description		Test	Potable Water Maximum Contaminant	Maximum Contaminant	SVO-ZV Potable Water	SVO-ZV Potable Water	SVO-ZV Potable Water		SVO-ZV Potable Water	SVO-ZV Potable Water		
Sample Description Sample Date		Conducted	Contaminant Level	Level	10/20/2009	11/10/2009	1/6/2010	2/3/2010	3/3/2010	3/31/2010	Potable Water 4/26/2010	5/18/2010
Analysis/Sample ID	Units	by	(MCL)	Source	20091214003	20091130010	20100222005		20100326001	20100603002	20100527009	
Anatysis/Sample 1D	Umis	Бу	(MCL)	Source	20091214003	20091130010	20100222003	20100222007	20100320001	20100003002	20100327009	20100327013
Physical Characteristics												
pH	pH units	U.S.	5.5-9.0	MORD	7.00	7.00	6.69	6.96	7.14	6.50	6.29	6.46
Conductivity	μS/cm	U.S.			334	304	362	362	363	300	251	269
Turbidity	NTU	U.S.	1.5*	MORD	1.6	0.3	2.0	3.8	3.5	NA	1.1	1.3
Total Dissolved Solids	mg/L	U.S.	100 (1,000*)	MORD	NA	180	215	222	NA	NA	NA	141
Iodine (LCV)												
Total I	mg/L	U.S.	0.05	MORD	< 0.05	< 0.05	< 0.05	< 0.05	< 0.05	< 0.05	< 0.05	< 0.05
Anions (IC/ISE)												
Bromide	mg/L	U.S.			< 0.5	<0.5	< 0.5	< 0.5	< 0.5	< 0.5	< 0.5	< 0.5
Chloride	mg/L	U.S.	250	MORD	9.75	8.91	9.79	9.81	10.1	10.2	8.61	8.68
Fluoride	mg/L	U.S.	1.5/4	MORD/EPA	<0.1	<0.1	<0.1	<0.1	<0.1	0.1	<0.1	<0.1
Nitrate as Nitrogen (NO3-N)	mg/L	U.S.	10	MORD/EPA	0.33	<0.11	0.45	0.49	0.57	0.62	<0.11	<0.11
Nitrite as Nitrogen (NO2-N)	mg/L	U.S.	1	EPA	NA	NA	NA	NA	NA	NA	NA	NA
Phosphate as P (PO4-P)	mg/L	U.S.			< 0.24	< 0.24	< 0.24	< 0.24	< 0.24	< 0.24	< 0.24	< 0.24
Sulfate	mg/L	U.S.	250	MORD	38	28.4	40.7	40.8	41.8	42.1	25.5	26.3
Cations (IC)												
Ammonia as Nitrogen (NH3-N)	mg/L	U.S.	2/1	MORD/SWEG	< 0.002	< 0.002	< 0.002	< 0.002	< 0.002	< 0.002	< 0.002	< 0.002
Lithium	mg/L	U.S.			< 0.002	< 0.002	< 0.002	< 0.002	< 0.002	< 0.002	< 0.002	< 0.002
Metals (ICP/MS)												
Calcium	mg/L	U.S.	100	MORD	44.4	39.5	48.6	49.5	47.4	49.9	33.1	32.9
Magnesium	mg/L	U.S.	50	MORD	11.1	10.1	11	11.1	11.2	12.1	8.65	8.19
Potassium Sodium	mg/L mg/L	U.S. U.S.			2.29 6.37	3.01 7.51	2.05 6.33	2.1 6.34	2.27 6.09	2.2 6.28	2.1 9.62	2.04 8.98
Aluminum	mg/L μg/L	U.S.			54	103	36	41	44	105	170	166
Antimony	μg/L μg/L	U.S.	6	EPA	<2	<2	<2	<2	<2	<2	<2	<2
Arsenic	μg/L μg/L	U.S.	10	MORD/EPA	<1	<1	<1	<1	<1	<1	<1	<1
Barium	μg/L	U.S.	1,000/10,000	MORD/SWEG	25	7	33	33	36	36	6	6
Beryllium	μg/L	U.S.	4	EPA	<1	<1	<1	<1	<1	<1	<1	<1
Cadmium	μg/L	U.S.	5/22	MORD/SWEG	<1	3	<1	<1	<1	<1	1	1
Chromium	μg/L	U.S.	100	MORD/EPA	<5	<5	<5	<5	<5	<5	<5	<5
Cobalt	μg/L	U.S.			<1	<1	<1	<1	<1	<1	<1	<1
Copper	μg/L	U.S.	1,000/1,300	MORD/EPA	2	<1	4	4	4	4	2	1
Iron	μg/L	U.S.	300	MORD	57	76	54	68	102	60	50	56
Lead	μg/L	U.S.	50/15	MORD/EPA	<1	<1	<1	<1	<1	<1	<1	<1
Manganese	μg/L	U.S.	50/300	MORD/SWEG	33	36	35	35	35	35	27	27
Mercury	μg/L	U.S.	2	MORD/EPA	<0.5	<0.5	<0.5	<0.5	<0.5	< 0.5	<0.5	< 0.5
Molybdenum Nickel	μg/L	U.S.	40	EPA HA MORD/SWEG	<1	<1	<1	<1 r	<1	<1	<1	<1 1
Nickel Selenium	μg/L	U.S. U.S.	100/300	MORD/SWEG MORD/EPA	4 <1	3 <1	24	5	5	5	2	
Silver	μg/L μg/L	U.S.	10/50 500/400	MORD/SWEG	<1 111	<1 28	<1 173	<1 405	<1 276	<1 187	<1 132	<1 160
Silver, Dissolved	μg/L μg/L	U.S.	300/400	MORD/SWEG	81	18	56	110	69	90	86	125
Zinc	μg/L μg/L	U.S.	5,000/2,000	MORD/SWEG	30	8	132	27	24	29	7	5
	µg/L	0.0.	0,000/2,000	o	30	J	102	~'		20	' ' ' ' ' ' ' ' ' ' ' ' ' ' ' ' ' ' '	,
L	1			0	1		·	·	1	I.	·	·

NA=Not analyzed; MI=Matrix interference

^{*}MORD limit 1.5 mg/L (Russian method)

^{**}limit does not include contribution from formate #TDS allowable limit after mineralization

Mission					Sovuz 19/Evp. 21	ISS ULF3/Exp. 21	155 204	/Exp. 22	Sovuz 20/Evp. 22	Soyuz 21/Exp. 23	ISS III E	4/Exp. 23
MANAGEME					SVO-ZV	SVO-ZV	SVO-ZV	SVO-ZV	SVO-ZV	SVO-ZV	SVO-ZV	SVO-ZV
Sample Location			Potable Water		SVU-ZV	SVU-ZV	SVU-ZV	3VU-ZV	SVU-ZV	SVU-ZV	SVU-ZV	SVU-ZV
Sample Location			Maximum	Maximum								
Sample Description		Test	Contaminant	Contaminant	Potable Water	Potable Water	Potable Water	Potable Water	Potable Water	Potable Water	Potable Water	Potable Water
Sample Date		Conducted	Level	Level	10/20/2009	11/10/2009	1/6/2010	2/3/2010	3/3/2010	3/31/2010	4/26/2010	5/18/2010
Analysis/Sample ID	Units	by	(MCL)	Source	20091214003	20091130010	20100222005	20100222007	20100326001	20100603002	20100527009	20100527013
Silicon (ICP/MS)												
Silicon (ICP/MS)	μg/L	U.S.			NA	NA	NA	NA	NA	NA	NA	NA
Total Organic Carbon (Sievers)												
Total Inorganic Carbon	mg/L	U.S.			31.5	29.7	31.7	32.9	31.4	31.9	25.8	25.6
Total Organic Carbon	mg/L	U.S.	20**	MORD	2.70	0.36	3.38	3.26	3.48	3.53	0.48	0.48
Volatile Organics												
3	/T	U.S.	15,000	SWEG	<2	<2	<2	<2	<2	<6	<2	<2
Acetone Acryloniltrile	μg/L μg/L	U.S.	13,000	SWEG	<2 <2	<2 <2	<2 <2	<2 <2	<2 <2	<6	<2 <2	<2 <2
Allyl chloride (3-Chloropropene)	μg/L μg/L	U.S.			<2 <2	<2 <2	<2 <2	<2 <2	<2 <2	<6	<2 <2	<2 <2
Benzene	μg/L μg/L	U.S.	5	EPA	<2. <0.4	<2 <0.4	<2 <0.4	<0.4	<2 <0.4	<0 <1.2	<2 <0.4	<2 <0.4
Bromobenzene	μg/L μg/L	U.S.	<u> </u>	EPA	<0.4	<0.4	<0.4	<0.4	<0.4	<1.2	<0.4	<0.4
Bromochloromethane	μg/L μg/L	U.S.	90	EPA HA	<4	<4	<4	<4	<4	<1.2	<4	<4
Bromodichloromethane	μg/L μg/L	U.S.	THM 80	EPA HA EPA	1.9	<0.4	3	3.1	3.5	2.6	<0.4	<0.4
Bromoform	μg/L μg/L	U.S.	THM 80	EPA	<2	<2	<2	<2	<2	<6	<2	<2
Bromomethane	μg/L	U.S.	10	EPA HA	<2	<2	<2	<2	<2	<6	<2	<2
2-Butanone (Methyl ethyl ketone)	μg/L	U.S.	4,000	EPA HA	<2	<2	<2	<2	<2	<6	<2	<2
n-Butylbenzene	μg/L	U.S.	1,000	LITTIE	<0.4	<0.4	<0.4	<0.4	<0.4	<1.2	<0.4	<0.4
sec-Butylbenzene	μg/L	U.S.			<0.4	<0.4	<0.4	<0.4	<0.4	<1.2	<0.4	<0.4
tert-Butylbenzene	μg/L	U.S.			<0.4	<0.4	<0.4	<0.4	<0.4	<1.2	<0.4	<0.4
Carbon disulfide	μg/L	U.S.			<2	<2	<2	<2	<2	<6	<2	<2
Carbon tetrachloride	μg/L	U.S.	5	EPA	< 0.4	< 0.4	< 0.4	< 0.4	< 0.4	<1.2	< 0.4	< 0.4
Chloroacetonitrile	μg/L	U.S.			<10	<10	<10	<10	<10	<30	<10	<10
Chlorobenzene	μg/L	U.S.	100	EPA	< 0.4	< 0.4	< 0.4	< 0.4	< 0.4	<1.2	< 0.4	< 0.4
1-Chlorobutane (Butyl chloride)	μg/L	U.S.			< 0.4	< 0.4	< 0.4	< 0.4	< 0.4	<1.2	< 0.4	< 0.4
Chloroethane	μg/L	U.S.			<2	<2	<2	<2	<2	<6	<2	<2
Chloroform	μg/L	U.S.	6,500/THM 80	SWEG/EPA	51.9	1.6	67.9	76.8	89.2	47.8	1.4	1.6
Chloromethane	μg/L	U.S.	30	EPA HA	<2	<2	<2	<2	<2	<6	<2	<2
2-Chlorotoluene	μg/L	U.S.	100	EPA HA	< 0.4	< 0.4	< 0.4	< 0.4	< 0.4	<1.2	< 0.4	< 0.4
4-Chlorotoluene	μg/L	U.S.	100	EPA HA	< 0.4	< 0.4	< 0.4	< 0.4	< 0.4	<1.2	< 0.4	< 0.4
Dibromochloromethane	μg/L	U.S.	THM 80	EPA	< 0.4	< 0.4	< 0.4	< 0.4	< 0.4	<1.2	< 0.4	< 0.4
1,2-Dibromo-3-chloropropane (DBCP)	μg/L	U.S.	0.2	EPA	<2	<2	<2	<2	<2	<6	<2	<2
1,2-Dibromoethane (EDB)	μg/L	U.S.	0.05	EPA	< 0.4	< 0.4	< 0.4	< 0.4	< 0.4	<1.2	< 0.4	< 0.4
Dibromomethane	μg/L	U.S.			< 0.4	< 0.4	< 0.4	< 0.4	< 0.4	<1.2	< 0.4	< 0.4
1,2-Dichlorobenzene	μg/L	U.S.	600	EPA	< 0.4	< 0.4	< 0.4	< 0.4	< 0.4	<1.2	< 0.4	< 0.4
1,3-Dichlorobenzene	μg/L	U.S.	600	EPA HA	< 0.4	< 0.4	< 0.4	< 0.4	< 0.4	<1.2	< 0.4	< 0.4
1,4-Dichlorobenzene	μg/L	U.S.	75	EPA	< 0.4	< 0.4	< 0.4	< 0.4	< 0.4	<1.2	< 0.4	< 0.4
trans-1,4-Dichloro-2-butene	μg/L	U.S.			< 0.4	< 0.4	< 0.4	< 0.4	< 0.4	<1.2	< 0.4	< 0.4
Dichlorodifluoromethane	μg/L	U.S.	1,000	EPA HA	<2	<2	<2	<2	<2	<6	<2	<2
1,1-Dichloroethane	μg/L	U.S.			<0.4	<0.4	<0.4	<0.4	<0.4	<1.2	<0.4	<0.4
1,2-Dichloroethane	μg/L	U.S.	5	EPA	< 0.4	<0.4	< 0.4	<0.4	<0.4	<1.2	<0.4	< 0.4
1,1-Dichloroethene	μg/L	U.S.	7	EPA	< 0.4	<0.4	< 0.4	<0.4	<0.4	<1.2	<0.4	<0.4
cis1,2-Dichloroethene	μg/L	U.S.	70	EPA	<0.4	<0.4	<0.4	<0.4	<0.4	<1.2	<0.4	<0.4
trans-1,2-Dichloroethene	μg/L	U.S.	100	EPA	< 0.4	<0.4	<0.4	<0.4	<0.4	<1.2	<0.4	<0.4
1,2-Dichloropropane	μg/L	U.S.	5	EPA	<0.4	<0.4	< 0.4	<0.4	<0.4	<1.2	<0.4	<0.4
1,3-Dichloropropane	μg/L	U.S.			<0.4	<0.4	< 0.4	<0.4	<0.4	<1.2	<0.4	< 0.4
2,2-Dichloropropane	μg/L	U.S.			< 0.4	< 0.4	< 0.4	< 0.4	< 0.4	<1.2	< 0.4	< 0.4

NA=Not analyzed; MI=Matrix interference
*MORD limit 1.5 mg/L (Russian method)
**limit does not include contribution from form

^{**}limit does not include contribution from formate #TDS allowable limit after mineralization SWEG - 1000 days (5-2006)

Appendix 2. ISS SVO-ZV Potable Water Summary for Expeditions 21 through 25

Mission					Sovuz 10/Evn 21	ISS ULF3/Exp. 21	ISS 20A	/Exp. 22	Sovuz 20/Exp. 22	Sovuz 21/Evn 23	ICC III F	4/Exp. 23
					SVO-ZV	SVO-ZV	SVO-ZV	SVO-ZV	SVO-ZV	SVO-ZV	SVO-ZV	SVO-ZV
Sample Location			Potable Water		3VU-ZV	3 V O-Z V	3VU-ZV	300-20	3VU-ZV	3VU-ZV	3VU-ZV	3VU-ZV
Sample Location			Maximum	Maximum								
Sample Description		Test	Contaminant	Contaminant	Potable Water	Potable Water	Potable Water	Potable Water	Potable Water	Potable Water	Potable Water	Potable Water
Sample Date		Conducted	Level	Level	10/20/2009	11/10/2009	1/6/2010	2/3/2010	3/3/2010	3/31/2010	4/26/2010	5/18/2010
Analysis/Sample ID	Units	by	(MCL)	Source	20091214003	20091130010	20100222005	20100222007	20100326001	20100603002	20100527009	20100527013
1,1-Dichloropropanone	μg/L	U.S.			<2	<2	<2	<2	<2	<6	<2	<2
1,1-Dichloropropene	μg/L	U.S.			< 0.4	< 0.4	< 0.4	< 0.4	< 0.4	<1.2	< 0.4	< 0.4
cis-1,3-Dichloropropene	μg/L	U.S.			< 0.4	< 0.4	< 0.4	< 0.4	< 0.4	<1.2	< 0.4	< 0.4
trans-1,3-Dichloropropene	μg/L	U.S.			<2	<2	<2	<2	<2	<6	<2	<2
Diethyl ether	μg/L	U.S.			<2	<2	<2	<2	<2	<6	<2	<2
Ethylbenzene	μg/L	U.S.	700	EPA	< 0.4	< 0.4	< 0.4	< 0.4	< 0.4	<1.2	< 0.4	< 0.4
Ethyl methacrylate	ug/L	U.S.			<2	<2	<2	<2	<2	<6	<2	<2
Hexachlorobutadiene	μg/L	U.S.	1	EPA HA	<2	<2	<2	<2	<2	<6	<2	<2
Hexachloroethane	ug/L	U.S.	1	EPA HA	<2	<2	<2	<2	<2	<6	<2	<2
2-Hexanone	μg/L	U.S.	_		<2	<2	<2	<2	<2	<6	<2	<2
Iodomethane	μg/L	U.S.			<2	<2	<2	<2	<2	<6	<2	<2
Isopropylbenzene (Cumene)	ug/L	U.S.	4.000	EPA DWEL	<0.4	<0.4	<0.4	<0.4	<0.4	<1.2	<0.4	< 0.4
4-Isopropyltoluene (Cymene)	μg/L	U.S.	1,000	DITTO TYLL	<0.4	<0.4	<0.4	<0.4	<0.4	<1.2	<0.4	<0.4
Methacrylonitrile	μg/L	U.S.			<2	<2	<2	<2	<2	<6	<2	<2
Methyl acrylate	μg/L	U.S.			<2	<2	<2	<2	<2	<6	<2	<2
Methyl-t-butylether (MTBE)	μg/L	U.S.			<2	<2	<2	<2	<2	<6	<2	<2
Methylene chloride (Dichloromethane)	μg/L	U.S.	15,000/5	SWEG/EPA	<0.4	<0.4	<0.4	<0.4	<0.4	<1.2	<0.4	<0.4
Methyl methacrylate	μg/L	U.S.	10,000/0	SVIEG/EI/I	<2	<2	<2	<2	<2	<6	<2	<2
4-Methyl-2-pentanone	μg/L	U.S.			<0.4	<0.4	<0.4	<0.4	<0.4	<1.2	<0.4	<0.4
Naphthalene	μg/L	U.S.	100	ЕРА НА	<0.4	<0.4	<0.4	<0.4	<0.4	<1.2	<0.4	<0.4
Nitrobenzene	μg/L	U.S.	100	LIAIA	<2	<2	<2	<2	<2	<6	<2	<2
2-Nitropropane	μg/L	U.S.			<2	<2	<2	<2	<2	<6	<2	<2
Pentachloroethane	μg/L	U.S.			<2	<2	<2	<2	<2	<6	<2	<2
Propionitrile (Ethyl cyanide)	μg/L	U.S.			<10	<10	<10	<10	<10	<30	<10	<10
n-Propylbenzene	μg/L	U.S.			<0.4	<0.4	<0.4	<0.4	<0.4	<1.2	<0.4	<0.4
Styrene	μg/L	U.S.	100	EPA	<0.4	<0.4	<0.4	<0.4	<0.4	<1.2	<0.4	<0.4
1.1.1.2-Tetrachloroethane	μg/L	U.S.	70	EPA HA	<0.4	<0.4	<0.4	<0.4	<0.4	<1.2	<0.4	<0.4
1.1.2.2-Tetrachloroethane	μg/L	U.S.	0.3	EPA HA	<0.4	<0.4	<0.4	<0.4	<0.4	<1.2	<0.4	<0.4
Tetrachloroethene	μg/L	U.S.	5	EPA	<0.4	<0.4	<0.4	<0.4	<0.4	<1.2	<0.4	<0.4
Tetrahydrofuran	μg/L	U.S.	J	LIA	<2	<2	<2	<2	<2	<6	<2	<2
Toluene	μg/L	U.S.	1.000	EPA	<0.4	<0.4	<0.4	<0.4	<0.4	<1.2	<0.4	<0.4
1,2,3-Trichlorobenzene	μg/L	U.S.	1,000	LIA	<0.4	<0.4	<0.4	<0.4	<0.4	<1.2	<0.4	<0.4
1.2.4-Trichlorobenzene	μg/L μg/L	U.S.	70	EPA	<0.4	< 0.4	<0.4	<0.4	<0.4	<1.2	<0.4	<0.4
1.1.1-Trichloroethane	μg/L μg/L	U.S.	200	EPA	<0.4	<0.4	<0.4	<0.4	<0.4	<1.2	<0.4	<0.4
1.1.2-Trichloroethane	μg/L μg/L	U.S.	5	EPA	<0.4	<0.4	<0.4	<0.4	<0.4	<1.2	<0.4	<0.4
Trichloroethene	цg/L цg/L	U.S.	5	EPA	<0.4	<0.4	<0.4	<0.4	<0.4	<1.2	<0.4	<0.4
Trichlorofluoromethane	μg/L μg/L	U.S.	2,000	EPA HA	<0.4	<0.4	<0.4	<0.4	<0.4	<1.2	<0.4	<0.4
1,2,3-Trichloropropane	μg/L μg/L	U.S.	40	EPA HA	<0.4	<0.4	<0.4	<0.4	<0.4	<1.2	<0.4	<0.4
1,2,4-Trimethylbenzene	μg/L μg/L	U.S.	40	ЕГА ПА	<0.4	<0.4	<0.4	<0.4	<0.4	<1.2	<0.4	<0.4
1.3.5-Trimethylbenzene	μg/L μg/L	U.S.			<0.4	<0.4	<0.4	<0.4	<0.4	<1.2	<0.4	<0.4
Vinyl Acetate	μg/L μg/L	U.S.			<0.4 <2	<0.4 <2	<0.4 <2	<0.4 <2	<0.4 <2	<1.2 <6	<0.4 <2	<0.4 <2
Vinyl Chloride		U.S.	2	EDA	<2 <2		<2 <2					<2 <2
7	μg/L		-	EPA		<2		<2	<2	<6	<2	
m&p-Xylene	μg/L	U.S.	Total Xylenes 10,000	EPA	<0.4	<0.4	<0.4	<0.4	<0.4	<1.2	< 0.4	< 0.4
o-Xylene	μg/L	U.S.	Total Xylenes 10,000	EPA	<0.4	<0.4	< 0.4	<0.4	< 0.4	<1.2	<0.4	< 0.4
	<u> </u>				<u> </u>		<u> </u>	<u> </u>			<u> </u>	

NA=Not analyzed; MI=Matrix interference *MORD limit 1.5 mg/L (Russian method)

^{**}limit does not include contribution from formate #TDS allowable limit after mineralization SWEG - 1000 days (5-2006)

Mission	I		1		Cover 10/E-m 21	ICC III E2/E 21	100 20 4	/Evn 22	Cover 20/E 22	Cover 21/E-m 22	ICC III E	4/Evn 22
MINION						ISS ULF3/Exp. 21		/Exp. 22	Soyuz 20/Exp. 22			4/Exp. 23
Samula I agastian			Datable Water		SVO-ZV	SVO-ZV	SVO-ZV	SVO-ZV	SVO-ZV	SVO-ZV	SVO-ZV	SVO-ZV
Sample Location			Potable Water Maximum	Maximum								
Sample Description		Test	Contaminant	Contaminant	Potable Water	Potable Water	Potable Water	Potable Water	Potable Water	Potable Water	Potable Water	Potable Water
Sample Date		Conducted	Level	Level	10/20/2009	11/10/2009	1/6/2010	2/3/2010	3/3/2010	3/31/2010	4/26/2010	5/18/2010
Analysis/Sample ID	Units	by	(MCL)	Source	20091214003	20091130010	20100222005	20100222007	20100326001	20100603002	20100527009	20100527013
Volatile Organics - Non-Targets (Tentatively Ident					20002222000	20001100010				201000000	2010002:000	20100021010
Acetaldehyde	μg/L	U.S.	4	,	not found	not found	not found	not found	not found	not found	not found	not found
Butyraldehyde (Butanal)	μg/L	U.S.			not found	not found	not found	not found	not found	not found	not found	not found
1,8-Cineole	μg/L	U.S.			not found	not found	not found	not found	not found	not found	not found	not found
Cyclohexanone	μg/L	U.S.			not found	not found	not found	not found	not found	not found	not found	not found
Difluorodimethylsilane	μg/L	U.S.			not found	not found	not found	not found	not found	not found	not found	not found
Dimethoxymethane (Formal)	μg/L	U.S.			not found	not found	not found	not found	not found	not found	not found	not found
2,5-Dimethylfuran	μg/L	U.S.			not found	not found	not found	not found	not found	not found	not found	not found
2,6-Dimethyl-1,7-octadiene	μg/L	U.S.			not found	not found	not found	not found	not found	not found	not found	not found
Dimethyl sulfide (Thiobismethane)	μg/L	U.S.			not found	not found	not found	not found	not found	not found	not found	not found
p-Dioxane	μg/L	U.S.			not found	not found	not found	not found	not found	not found	not found	not found
1,3-Dioxolane (Ethylene glycol formal)	μg/L	U.S.			not found	not found	not found	not found	not found	not found	not found	not found
Ethyl acetate	μg/L	U.S.			not found	not found	not found	not found	not found	not found	not found	not found
Fluorotrimethylsilane	μg/L	U.S.			not found	not found	not found	not found	not found	not found	not found	not found
2-Heptanone	μg/L μg/L	U.S.			not found	not found	not found	not found	not found	not found	not found	not found
•												
3-Hexanone	μg/L	U.S.			not found	not found	not found	not found	not found	not found	not found	not found
3-Hydroxy-3-methylbutyric acid	μg/L	U.S.			not found	not found	not found	not found	not found	not found	not found	not found
Isobutyronitrile	μg/L	U.S.			not found	not found	not found	not found	not found	not found	not found	not found
Isooctanol	μg/L	U.S.			not found	not found	not found	not found	not found	not found	not found	not found
L-Menthol	μg/L	U.S.			not found	not found	not found	not found	not found	not found	not found	not found
Menthone	μg/L	U.S.			not found	not found	not found	not found	not found	not found	not found	not found
Methyl acetate	μg/L	U.S.			not found	not found	not found	not found	not found	not found	not found	not found
cis-1-Methyl-4-(1-methylethenyl)-cyclohexane	μg/L	U.S.			not found	not found	not found	not found	not found	not found	not found	not found
3-Methyl-2-pentanone	μg/L	U.S.			not found	not found	not found	not found	not found	not found	not found	not found
2-Methyl-1-propene	μg/L	U.S.			not found	not found	not found	not found	not found	not found	not found	not found
alpha-Methyl styrene	μg/L	U.S.			not found	not found	not found	not found	not found	not found	not found	not found
2-Nonanone	μg/L	U.S.			not found	not found	not found	not found	not found	not found	not found	not found
2-Octanone	μg/L	U.S.			not found	not found	not found	not found	not found	not found	not found	not found
2-Pentanone	μg/L	U.S.			not found	not found	not found	not found	not found	not found	not found	not found
1,1,2-Trichloro-1,2,2-trifluoroethane (Freon 113)	μg/L	U.S.			not found	not found	not found	not found	not found	not found	not found	not found
Trimethylsilanol	μg/L	U.S.			not found	not found	not found	not found	not found	not found	not found	not found
Entro stable Or double	-											
Extractable Organics		11.0			10						10	10
Acetophenone	μg/L	U.S.			<16	<8	<8	<8	<8	<32	<16	<16
Benzaldehyde	μg/L	U.S.			<8	<4	<4	<4	<4	<16	<8	<8
Benzoic acid	μg/L	U.S.			<24	<12	<12	<12	<12	<48	<24	<24
Benzothiazole Royard alcohol	μg/L	U.S.			<8	<4	<4	<4	<4	<16	<8	<8
Benzyl alcohol	μg/L	U.S.	7 000	EDA DITET	<8	<4	<4	<4	<4	<16	<8	<8
Benzyl butyl phthlate	μg/L	U.S.	7,000	EPA DWEL	<8	<4	<4	<4	<4	<16	<8	<8
2-Butoxyethanol	μg/L	U.S.			<16	<8	<8	<8	<8	<32	<16	<16
2-(2-Butoxyethoxy)ethanol	μg/L	U.S.			<16	<8	<8	<8	<8	<32	<16	<16
2-(2-Butoxyethoxy)ethyl acetate	μg/L	U.S.			<8	<4	<4 <8	<4	<4	<16	<8 <16	<8
n-Butylpalmitate Butylated hydroxyanisole (BHA)	μg/L	U.S.			<16	<8		<8	<8	<32		<16
J J J , , , , ,	μg/L	U.S.			<8	<4	<4	<4	<4	<16	<8	<8
N-Butylbenzenesulfonamide	μg/L	U.S. U.S.			<8 <24	<4 <12	<4 <12	<4 <12	<4	<16	<8 <24	<8 <24
3-tert-Butylphenol	μg/L								<12	<48		<24 <8
Caffeine	μg/L	U.S.			<8	<4	<4	<4	<4	<16	<8	
tris-2-Chloroethyl phosphate	μg/L	U.S.			<8	<4	<4	<4	<4	<16	<8	<8

NA=Not analyzed; MI=Matrix interference

SWEG - 1000 days (5-2006)

^{*}MORD limit 1.5 mg/L (Russian method)

^{**}limit does not include contribution from formate #TDS allowable limit after mineralization

Mission					Soyuz 19/Exp. 21	ISS ULF3/Exp. 21	ISS 20A	/Exp. 22	Soyuz 20/Exp. 22	Soyuz 21/Exp. 23	ISS ULF	4/Exp. 23
				•	SVO-ZV	SVO-ZV	SVO-ZV	SVO-ZV	SVO-ZV	SVO-ZV	SVO-ZV	SVO-ZV
Sample Location Sample Description		Test	Potable Water Maximum Contaminant	Maximum Contaminant	Potable Water	Potable Water	Potable Water	Potable Water	Potable Water	Potable Water	Potable Water	Potable Water
1												
Sample Date		Conducted	Level	Level	10/20/2009	11/10/2009	1/6/2010	2/3/2010	3/3/2010	3/31/2010	4/26/2010	5/18/2010
Analysis/Sample ID Cholesterol	Units	by	(MCL)	Source	20091214003	20091130010	20100222005	20100222007	20100326001	20100603002	20100527009	20100527013
	μg/L	U.S.			<64	<32	<32	<32	<32	<128	<64	<64
o-Cresol (2-Methylphenol)	μg/L	U.S.			<8	<4	<4	<4	<4	<16	<8	<8
Cyclododecane	μg/L	U.S.			<8	<4	<4	<4	<4	<16	<8	<8
Decamethylcyclopentasiloxane	μg/L	U.S.			<8	<4	<4	<4	<4	<16	<8	<8
Decanoic acid	μg/L	U.S.			<16	<8	<8	<8	<8	<32	<16	<16
2,6-Di-t-butyl-1,4-benzoquinone	μg/L	U.S.			<8	<4	<4	<4	<4	<16	<8	<8
2,4-Di-t-butylphenol	μg/L	U.S.			<8	<4	<4	<4	<4	<16	<8	<8
1,4 Diacetylbenzene	μg/L	U.S.			<8	<4	<4	<4	<4	<16	<8	<8
N,N-Dibutylformamide	μg/L	U.S.	40.000/4.000	CHEC/EDA DIJET	<8	<4	<4	<4	<4	<16	<8	<8
Dibutyl phthalate	μg/L	U.S.	-,,	SWEG/EPA DWEI	<8	<4	<4	<4	<4	<16	<8	<8
Dibutylamine	μg/L	U.S.	Dialkylamines 300	SWEG	<8	<4	<4	<4	<4	<16	<8	<8
N,N-Diethyl-m-toluamide	μg/L	U.S.	00.000	EDA DIJET	<8	<4	<4	<4	<4	<16	<8	<8
Diethylphthalate	μg/L	U.S.	30,000	EPA DWEL	<8	<4	<4	<4	<4	<16	<8	<8
Diethylene glycol monoethyl ether	μg/L	U.S.			<8	<4	<4	<4	<4	<16	<8	<8
N,N-Diethylformamide	μg/L	U.S.			<24	<12	<12	<12	<12	<48	<24	<24
Diiodomethane (Methyl iodide)	μg/L	U.S.			<8	<4	<4	<4	<4	<16	<8	<8
Diisopropyl adipate	μg/L	U.S.			<8	<4	<4	<4	<4	<16	<8	<8
Dimethyl phthalate	μg/L	U.S.			<8	<4	<4	<4	<4	<16	<8	<8
N,N-Dimethyl acetamide	μg/L	U.S.			<8	<4	<4	<4	<4	<16	<8	<8
N,N-Dimethylbenzylamine	μg/L	U.S.	Dialkylamines 300	SWEG	<8	<4	<4	<4	<4	<16	<8	<8
N,N-Dimethylformamide	μg/L	U.S.			<16	<8	<8	<8	<8	<32	<16	<16
Dipropylene glycol methyl ether	μg/L	U.S.			<8	<4	<4	<4	<4	<16	<8	<8
Dodecamethylcyclohexasiloxane	μg/L	U.S.			<8	<4	<4	<4	<4	<16	<8	<8
2-Ethoxyethanol	μg/L	U.S.			<8	<4	<4	<4	<4	<16	<8	<8
2-Ethyl-1-hexanol	μg/L	U.S.			<8	<4	<4	<4	<4	<16	<8	<8
2-Ethylhexanoic acid	μg/L	U.S.			<8	<4	<4	<4	<4	<16	<8	<8
bis-2-Ethylhexyl adipate	μg/L	U.S.	400	EPA	<8	<4	<4	<4	<4	<16	<8	<8
bis-2-Ethylhexyl phthalate (Dioctyl phthlate)	μg/L	U.S.	20,000/6	SWEG/EPA	<8	<4	<4	<4	<4	<16	<8	<8
4-Ethylmorpholine	μg/L	U.S.			<8	<4	<4	<4	<4	<16	<8	<8
1-Formylpiperidine	μg/L	U.S.			<8	<4	<4	<4	<4	<16	<8	<8
Heptanoic acid	μg/L	U.S.			<8	<4	<4	<4	<4	<16	<8	<8
2-Heptanone	μg/L	U.S.			<8	<4	<4	<4	<4	<16	<8	<8
gamma-Hexalactone	μg/L	U.S.			<8	<4	<4	<4	<4	<16	<8	<8
Hexanoic acid	μg/L	U.S.			<16	<8	<8	<8	<8	<32	<16	<16
2-Hexanol	μg/L	U.S.			<8	<4	<4	<4	<4	<16	<8	<8
2-Hydroxybenzothiazole	μg/L	U.S.			<8	<4	<4	<4	<4	<16	<8	<8
Ibuprofen Lodoform	μg/L	U.S.			<8	<4	<4	<4	<4	<16	<8	<8
Iodoform	μg/L	U.S.	100	EDA HA	<8	<4	<4	<4	<4	<16	<8	<8
Isophorone	μg/L	U.S.	100	EPA HA	<8	<4	<4	<4	<4	<16	<8	<8
4-Isopropylphenol	μg/L	U.S.			<8	<4	<4	<4	<4	<16	<8	<8
Lauramide	μg/L	U.S.			<8	<4	<4	<4	<4	<16	<8	<8
Lauric acid (Dodecanoic acid)	μg/L	U.S.			<240	<120	<120	<120	<120	<480	<240	<240
p-Menth-1-en-8-ol (alpha-Terpineol)	μg/L	U.S.	00.000	CHEC	<8	<4	<4	<4	<4	<16	<8	<8
2-Mercaptobenzothiazole	μg/L	U.S.	30,000	SWEG	<80	<40	<40	<40	<40	<160	<80	<80
2-Methyl-2,4-pentanediol	μg/L	U.S.			<8	<4	<4	<4	<4	<16	<8	<8
1-Methyl-2-pyrrolidinone	μg/L	U.S.			<8	<4	<4	<4	<4	<16	<8	<8
Methyl-4-hydroxybenzoate	μg/L	U.S.			<8	<4	<4	<4	<4	<16	<8	<8
Methyl sulfone	μg/L	U.S.			<8	<4	<4	<4	<4	<16	<8	<8

NA=Not analyzed; MI=Matrix interference

^{*}MORD limit 1.5 mg/L (Russian method)

^{**}limit does not include contribution from formate #TDS allowable limit after mineralization

Mission					Soyuz 19/Exp. 21	ISS ULF3/Exp. 21	ISS 20A	/Exp. 22	Soyuz 20/Exp. 22	Soyuz 21/Exp. 23	ISS ULF	4/Exp. 23
					SVO-ZV	SVO-ZV	SVO-ZV	SVO-ZV	SVO-ZV	SVO-ZV	SVO-ZV	SVO-ZV
Sample Location			Potable Water									
			Maximum	Maximum								
Sample Description		Test	Contaminant	Contaminant	Potable Water	Potable Water	Potable Water	Potable Water	Potable Water	Potable Water	Potable Water	Potable Water
Sample Date		Conducted	Level	Level	10/20/2009	11/10/2009	1/6/2010	2/3/2010	3/3/2010	3/31/2010	4/26/2010	5/18/2010
Analysis/Sample ID	Units	by	(MCL)	Source	20091214003	20091130010	20100222005	20100222007	20100326001	20100603002	20100527009	20100527013
2-Methyl butyric acid	μg/L	U.S.			<24	<12	<12	<12	<12	<48	<24	<24
2-Methylthiobenzothiazole	μg/L	U.S.			<8	<4	<4	<4	<4	<16	<8	<8
Monomethyl phthalate	μg/L	U.S.			<8	<4	<4	<4	<4	<16	<8	<8
Myristic acid	μg/L	U.S.			<48	<24	<24	<24	<24	<96	<48	<48
(+)-Neomenthol	μg/L	U.S.			<8	<4	<4	<4	<4	<16	<8	<8
Nicotine	μg/L	U.S.			<8	<4	<4	<4	<4	<16	<8	<8
Nonadecane	μg/L	U.S.			<8	<4	<4	<4	<4	<16	<8	<8
Nonanoic acid	μg/L	U.S.			<24	<12	<12	<12	<12	<48	<24	<24
1-Octadecanol	μg/L	U.S.			<24	<12	<12	<12	<12	<48	<24	<24
Octamethylcyclotetrasiloxane	μg/L	U.S.			<8	<4	<4	<4	<4	<16	<8	<8
Octanoic acid	μg/L	U.S.			<16	<8	<8	<8	<8	<32	<16	<16
4-tert-Octylphenol	μg/L	U.S.			<8	<4	<4	<4	<4	<16	<8	<8
Oleic acid	μg/L	U.S.			<80	<40	<40	<40	<40	<160	<80	<80
Oxindole	μg/L	U.S.			<8	<4	<4	<4	<4	<16	<8	<8
Palmitic acid	μg/L	U.S.			<240	<120	<120	<120	<120	<480	<240	<240
Palmitoleic acid	μg/L	U.S.			NA	NA	NA	NA	NA	NA	NA	NA
Pentacosane	μg/L	U.S.			<8	<4	<4	<4	<4	<16	<8	<8
sec-Phenethyl alcohol	μg/L	U.S.			<8	<4	<4	<4	<4	<16	<8	<8
Phenol	μg/L	U.S.	1,000/4,000	MORD/SWEG	<8	<4	<4	<4	<4	<16	<8	<8
2-Phenoxyethanol	μg/L	U.S.			<8	<4	<4	<4	<4	<16	<8	<8
N-Phenyl-2-naphthylamine	μg/L	U.S.	260,000	SWEG	<8	<4	<4	<4	<4	<16	<8	<8
2-Phenyl-2-propanol	μg/L	U.S.			<8	<4	<4	<4	<4	<16	<8	<8
2-Phenylacetic acid	μg/L	U.S.			<32	<16	<16	<16	<16	<64	<32	<32
Phenethyl alcohol	μg/L	U.S.			<8	<4	<4	<4	<4	<16	<8	<8
2-Phenylphenol	μg/L	U.S.			<8	<4	<4	<4	<4	<16	<8	<8
Salicyclic Acid	μg/L	U.S.			<64	<32	<32	<32	<32	<128	<64	<64
trans-Squalene	μg/L	U.S.			<16	<8	<8	<8	<8	<32	<16	<16
Stearic acid	μg/L	U.S.			NA	NA	NA	NA	NA	NA	NA	NA
1-Tetradecanol	μg/L	U.S.			<8	<4	<4	<4	<4	<16	<8	<8
Tetramethylsuccinonitrile	μg/L	U.S.			<8	<4	<4	<4	<4	<16	<8	<8
Tetramethyl thiourea	μg/L	U.S.			<8	<4	<4	<4	<4	<16	<8	<8
Tetramethylurea	μg/L	U.S.			<8	<4	<4	<4	<4	<16	<8	<8
Thymol	μg/L	U.S.			<8	<4	<4	<4	<4	<16	<8	<8
1,3,5-Triallyl-1,3,5-triazine-2,4,6(1H,3H,5H)-trione	μg/L	U.S.	m · II · · · · ·		<8	<4	<4	<4	<4	<16	<8	<8
Tributylamine	μg/L	U.S.	Trialkylamines 400	SWEG	<8	<4	<4	<4	<4	<16	<8	<8
Tributyl phosphate	μg/L	U.S.			<8	<4	<4	<4	<4	<16	<8	<8
Triethyl phosphate	μg/L	U.S.			<16	<8	<8	<8	<8	<32	<16	<16
2,2,4-Trimethyl-1,3-pentanediol diisobutyrate	μg/L	U.S.			<16	<8	<8	<8	<8	<32	<16	<16
Tripropylene glycol monomethyl ether	μg/L	U.S.			<8	<4	<4	<4	<4	<16	<8	<8
Undecanoic acid	μg/L	U.S.			<48	<24	<24	<24	<24	<96	<48	<48
2-Undecanone	μg/L	U.S.			<8	<4	<4	<4	<4	<16	<8	<8
Valeric acid (Pentanoic acid)	μg/L	U.S.			<48	<24	<24	<24	<24	<96	<48	<48
Vanillin	μg/L	U.S.			<16	<8	<8	<8	<8	<32	<16	<16
AL LI MAYGGREE								.			.	
Alcohols (DAI/GC/MS)					4	a						4
1-Butanol	μg/L	U.S.			<100	<100	<100	<100	<100	<100	<100	<100
2-Butanol	μg/L	U.S.			<100	<100	<100	<100	<100	<100	<100	<100
Ethanol	μg/L	U.S.			<100	<100	<100	<100	<100	<100	<100	<100

NA=Not analyzed; MI=Matrix interference

^{*}MORD limit 1.5 mg/L (Russian method)

^{**}limit does not include contribution from formate #TDS allowable limit after mineralization

Mission					Sovuz 19/Evn 21	ISS ULF3/Exp. 21	ISS 20A	/Exp. 22	Soyuz 20/Exp. 22	Sovuz 21/Evn 23	ISS III F	4/Exp. 23
				ļ	SVO-ZV	SVO-ZV	SVO-ZV	SVO-ZV	SVO-ZV	SVO-ZV	SVO-ZV	SVO-ZV
Sample Location			Potable Water		3V 0-ZV	3VU-ZV	3VU-ZV	3VU-ZV	3VU-ZV	3VU-ZV	3VU-ZV	3VU-ZV
Sample Location			Maximum	Maximum								
Sample Description		Test	Contaminant	Contaminant	Potable Water	Potable Water	Potable Water	Potable Water	Potable Water	Potable Water	Potable Water	Potable Water
Sample Date		Conducted	Level	Level	10/20/2009	11/10/2009	1/6/2010	2/3/2010	3/3/2010	3/31/2010	4/26/2010	5/18/2010
Analysis/Sample ID	Units	by	(MCL)	Source	20091214003	20091130010	20100222005	20100222007	20100326001	20100603002	20100527009	20100527013
Methanol	μg/L	U.S.	, ,		<100	<100	<100	<100	<100	<100	<100	<100
2-Methyl-1-butanol	μg/L	U.S.			<100	<100	<100	<100	<100	<100	<100	<100
2-Methyl-2-butanol	μg/L	U.S.			<100	<100	<100	<100	<100	<100	<100	<100
3-Methyl-1-butanol (Isopentanol)	μg/L	U.S.			<100	<100	<100	<100	<100	<100	<100	<100
2-Methyl-1-propanol	μg/L	U.S.			<100	<100	<100	<100	<100	<100	<100	<100
2-Methyl-2-propanol	μg/L	U.S.			<100	<100	<100	<100	<100	<100	<100	<100
1-Pentanol (Amyl alcohol)	μg/L	U.S.			<100	<100	<100	<100	<100	<100	<100	<100
2-Pentanol (sec-Amyl alcohol)	μg/L	U.S.			<100	<100	<100	<100	<100	<100	<100	<100
3-Pentanol	μg/L	U.S.			<100	<100	<100	<100	<100	<100	<100	<100
1-Propanol	μg/L	U.S.			<100	<100	<100	<100	<100	<100	<100	<100
2-Propanol (Isopropanol)	μg/L	U.S.			<100	<100	<100	<100	<100	<100	<100	<100
(wopropunos)	ия/ ш	0.0.			100	1100	100	100	100	100	100	1100
Glycols (DAI/GC/MS)							1					
1,2-Ethanediol (Ethylene glycol)	μg/L	U.S.	12000/14000	MORD/EPA HA	<1000	<1000	<1000	<1000	<1000	<1000	<1000	<1000
1,2-Propanediol (Propylene glycol)	μg/L	U.S.	12000/11000	MOND LITTIN	<500	<500	<500	<500	<500	<500	<500	<500
1,5 Frommentor (Fropyrene grycor)	µg/L	0.5.			1000	1000	1000	1000	1000	1000	1000	1000
Glycerol (LC/MS/MS)												
Glycerol (1,2,3-Propanetriol)	μg/L	U.S.			NA	NA	NA	NA	NA	NA	NA	NA
Cifyceror (1,5,0 1 ropunction)	µg/L	0.5.			IVA	III	IVA	IVA	IVA	INA	IVA	IVA
Silanols (LC/MS/MS) (R&D Method - NIST tracea	ble stan	dard not ava	ilable)									
Dimethylsilanediol (DMSD)	μg/L	U.S.			NA	NA	NA	NA	NA	NA	NA	NA
	rs-	0.01				****				1112	1112	1112
Carboxylates (CE)												
Acetate	μg/L	U.S.			<125	<125	<125	<125	<125	<125	<125	<125
Formate	μg/L	U.S.	2.500.000	SWEG	<125	<125	<125	<125	<125	<125	<125	<125
Glycolate	μg/L	U.S.	2,000,000	51124	<125	<125	<125	<125	<125	<125	<125	<125
Glyoxylate	μg/L	U.S.			<125	<125	<125	<125	<125	<125	<125	<125
Lactate	μg/L	U.S.			<1000	<1000	<1000	<1000	<1000	<1000	<1000	<1000
Oxalate	μg/L	U.S.			<125	<125	MI	<125	<125	<125	<125	<125
Propionate	μg/L	U.S.			<125	<125	<125	<125	<125	<125	<125	<125
	ra-											
Aldehydes												
Formaldehyde	ug/L	U.S.	12.000/1.000	SWEG/EPA HA	<2	20	<2	<2	<2	3	4	4
	rs-	0.0.	_2,000, 1,000	,		~~	12		120	, , ,	•	•
Amines (CE)							1					
Ethylamine	μg/L	U.S.	Monoalkylamines 2000	SWEG	<125	<125	<125	<125	<125	<125	<125	<125
Methylamine	μg/L μg/L	U.S.	Monoalkylamines 2000	SWEG	<125	<125	<125	<125	<125	<125	<125	<125
n-Propylamine	μg/L	U.S.	Monoalkylamines 2000	SWEG	<125	<125	<125	<125	<125	<125	<125	<125
Trimethylamine	μg/L	U.S.	Trialkylamines 400	SWEG	<125	<125	<125	<125	<125	<125	<125	<125
,/	ry-	0.0.	, , , , , , , , , , , , , , , , , , ,	51124	-1200	1280	1120	1220	1200	1200	1220	1200
Non-volatiles (LC/UV-VIS)							1					
Urea	μg/L	U.S.			<800	<800	<800	<800	<800	<800	<800	<800
Caprolactam	μg/L	U.S.	100.000	SWEG	<8	<4	<4	<4	<4	<16	<8	<8
	μ ₅ /	0.0.	100,000	DIILG	,,	``	``	``	` .	110	``	``
Organic Carbon Recovery	percent	U.S.			0.20	2.25	0.21	0.24	0.27	0.18	0.36	0.37
Unaccounted Organic Carbon	mg/L	U.S.			2.69	0.35	3.37	3.25	3.47	3.52	0.30	0.47
AMERICA OF BRIDE ONL DON	mg/L	0.5.	l l		2.00	0.00	3.31	3.23	0.71	0.02	0.10	0.71

SWEG - 1000 days (5-2006)

NA=Not analyzed; MI=Matrix interference

^{*}MORD limit 1.5 mg/L (Russian method)

^{**}limit does not include contribution from formate #TDS allowable limit after mineralization

Mission		1		1	Cover 1	9/Exp. 21	166 111 1	F3/Exp. 21	ISS 20A/Exp. 22	ISS 19A	/Evn 22	Soyuz 21/Exp. 23
Mission					WPA PWD Hot	WPA PWD	WPA PWD Hot	WPA PWD	WPA PWD	WPA PWD	WPA PWD Hot	WPA PWD Hot
					WITH WE HOL	Ambient	WITH WE HOL	Ambient	Ambient	Ambient	WITH WE HOL	WITTI WE HOU
Sample Location			Potable Water									
Comple Description			Maximum Contaminant	Maximum Contaminant	Potable Water	Potable Water	Datable Water	Potable Water	Potable Water	Potable Water	Potable Water	Potable Water
Sample Description Sample Date		Test Conducted	Contaminant Level	Level	10/20/2009	10/20/2009	Potable Water 11/10/2009	11/10/2009	1/6/2010	3/3/2010	3/3/2010	3/31/2010
Analysis/Sample ID	Units	by	(MCL)	Source	20091214001	20091214002	20091130012	20091130013	20100222001	20100422006	20100422007	20100603003
Analysis bampic 1D	Omts	Бу	(MCL)	Source	20031214001	20031214002	20031130012	20031130013	20100222001	20100422000	20100422007	20100003003
Physical Characteristics												
pH	pH units	U.S.	4.5-8.5	41000	7.38	7.09	7.06	6.93	6.91	6.74	6.55	7.15
Conductivity	μS/cm	U.S.			2	2	2	2	14	2	2	2
Turbidity	NTU	U.S.	1	41000	NA	< 0.1	<0.1	<0.1	<0.1	<0.1	0.2	0.2
Total Solids	mg/L	U.S.	100	41000	NA	NA	<5	<5	<5	<5	<5	NA
Iodine (LCV)												
. ,				41000 (tl I max/tl I at								
Total I	mg/L	U.S.	6/0.2	pt of consumption)	< 0.05	< 0.05	< 0.05	< 0.05	< 0.05	< 0.05	< 0.05	< 0.05
Iodine	mg/L	U.S.			< 0.05	<0.05	< 0.05	<0.05	<0.05	< 0.05	<0.05	<0.05
Iodide	mg/L	U.S.			< 0.05	< 0.05	< 0.05	< 0.05	< 0.05	< 0.05	< 0.05	< 0.05
Anions (IC/ISE)												
Bromide	mg/L	U.S.			< 0.5	< 0.5	< 0.5	< 0.5	< 0.5	< 0.5	< 0.5	< 0.5
Chloride	mg/L	U.S.			< 0.15	< 0.15	< 0.15	< 0.15	< 0.15	< 0.15	< 0.15	< 0.15
Fluoride	mg/L	U.S.			< 0.1	<0.1	<0.1	<0.1	<0.1	<0.1	<0.1	<0.1
Nitrate as Nitrogen (NO3-N)	mg/L	U.S.	10	41000	<0.11	<0.11	<0.11	<0.11	<0.11	<0.11	<0.11	<0.11
Nitrite as Nitrogen (NO2-N)	mg/L	U.S.			NA	NA	NA	NA	NA	NA	NA	NA
Phosphate as P (PO4-P) Sulfate	mg/L mg/L	U.S. U.S.	250	41000	<0.24 <0.75	<0.24 <0.75	<0.24 <0.75	<0.24 <0.75	<0.24 <0.75	<0.24 <0.75	<0.24 <0.75	<0.24 <0.75
Sunde	mg/L	U.S.	250	41000	<0.75	<0.75	<0.75	<0.75	<0.75	<0.75	<0.75	<0.75
Cations (IC)												
Ammonia as Nitrogen (NH3-N)	mg/L	U.S.	1	SWEG&41000	0.061	< 0.002	0.039	< 0.002	< 0.002	< 0.002	< 0.002	< 0.002
Lithium	mg/L	U.S.			< 0.002	< 0.002	< 0.002	< 0.002	< 0.002	< 0.002	< 0.002	< 0.002
Metals (ICP/MS)												
Calcium	mg/L	U.S.	30	41000	< 0.01	< 0.01	< 0.01	< 0.01	< 0.01	< 0.01	0.01	0.03
Magnesium	mg/L	U.S.	50	41000	<0.01	< 0.01	<0.01	< 0.01	< 0.01	<0.01	<0.01	<0.01
Potassium	mg/L	U.S.	340	41000	< 0.01	< 0.01	0.01	< 0.01	0.03	< 0.01	< 0.01	< 0.01
Sodium	mg/L	U.S.			< 0.01	< 0.01	0.02	< 0.01	0.02	< 0.01	< 0.01	< 0.01
Aluminum	μg/L	U.S.			<2	<2	<2	<2	<2	<2	<2	<2
Antimony	μg/L	U.S.	2,000	SWEG	<2	<2	<2	<2	<2	<2	<2	<2
Arsenic	μg/L	U.S.	10	41000	<1	<1	<1	<1	<1	<1	<1	<1
Barium	μg/L	U.S.	10,000	SWEG&41000	<1	<1	<1	<1	<1	<1	<1	<1
Beryllium	μg/L	U.S.		arrana :	<1	<1	<1	<1	<1	<1	<1	<1
Cadmium	μg/L	U.S.	22	SWEG&41000	<1	<1	<1	<1	<1	<1	<1	<1
Chromium	μg/L	U.S.	230	41000	<5	<5	<5	<5	<5	<5	<5	<5
Cobalt	μg/L	U.S.	1.000	41000	<1	<1	<1	<1	<1	<1	<1	<1
Copper Iron	μg/L μg/L	U.S. U.S.	1,000 300	41000 41000	<1 <5	<1 <5	<1 42	<1 <5	<1 <5	<1 <5		<1 <5
Lead	μg/L μg/L	U.S.	50	41000	<5 <1	<0 <1	42 <1	<3 <1	 <1 	<5 <1	<5 <1	<5 <1
Manganese	μg/L μg/L	U.S.	300	SWEG&41000	<1	<1	<1	<1	<1	<1	<1	<1
Mercury	μg/L μg/L	U.S.	2	41000	<0.5	<0.5	<0.5	<0.5	<0.5	<0.5	<0.5	<0.5
Molybdenum	μg/L	U.S.	~	11000	<1	<1	<1	<1	<1	<1	<1	<1
Nickel	μg/L	U.S.	300	SWEG&41000	17	24	25	33	10	34	26	40
Selenium	μg/L	U.S.	10	41000	<1	<1	<1	<1	<1	<1	1	<1
Silver	μg/L	U.S.	400	SWEG&41000	<2	<2	<2	<2	<2	<2	<2	<2

Mission					Sovuz 1	9/Exp. 21	ISS UL	F3/Exp. 21	ISS 20A/Exp. 22	ISS 19A	/Exp. 23	Soyuz 21/Exp. 23
					WPA PWD Hot	WPA PWD Ambient	WPA PWD Hot	WPA PWD Ambient	WPA PWD Ambient	WPA PWD Ambient	WPA PWD Hot	WPA PWD Hot
Sample Location Sample Description Sample Date		Test Conducted	Potable Water Maximum Contaminant Level	Maximum Contaminant Level	Potable Water 10/20/2009	Potable Water 10/20/2009	Potable Water 11/10/2009	Potable Water 11/10/2009	Potable Water 1/6/2010	Potable Water 3/3/2010	Potable Water 3/3/2010	Potable Water 3/31/2010
Analysis/Sample ID	Units	by	(MCL)	Source	20091214001	20091214002	20091130012	20091130013	20100222001	20100422006	20100422007	20100603003
Zinc	μg/L	U.S.	2,000	SWEG&41000	<1	<1	<1	<1	<1	<1	<1	<1
Silicon (ICP/MS)												
Silicon (ICP/MS)					NA	NA	NA	43	22	136	185	NA
Sincon (ICI7MS)		1			INA	IVA	IVA	10	22	130	100	IVA
Total Organic Carbon (Sievers)												
Total Inorganic Carbon	mg/L	U.S.			0.71	0.73	1.09	1.08	1.22	1.34	1.23	1.70
Total Organic Carbon	mg/L	U.S.	3	41000	0.30	0.22	0.22	0.14	0.16	0.17	0.16	0.26
Volatile Organics												
Acetone	μg/L	U.S.	15,000	SWEG	<2	<2	<2	<2	<2	39	11	<2
Acryloniltrile	μg/L	U.S.			<2	<2	<2	<2	<2	<2	<2	<2
Allyl chloride (3-Chloropropene)	μg/L	U.S.			<2	<2	<2	<2	<2	<2	<2	<2
Benzene	μg/L	U.S.	70/5	SWEG/EPA	< 0.4	< 0.4	< 0.4	< 0.4	< 0.4	< 0.4	< 0.4	< 0.4
Bromobenzene	μg/L	U.S.			< 0.4	< 0.4	< 0.4	< 0.4	< 0.4	< 0.4	< 0.4	< 0.4
Bromochloromethane	μg/L	U.S.	90	EPA HA	<4	<4	<4	<4	<4	<4	<4	<4
Bromodichloromethane	μg/L	U.S.	THM 80	EPA	< 0.4	< 0.4	< 0.4	< 0.4	< 0.4	< 0.4	< 0.4	< 0.4
Bromoform	μg/L	U.S.	THM 80	EPA	<2	<2	<2	<2	<2	<2	<2	<2
Bromomethane	μg/L	U.S.	10	EPA HA	<2	<2	<2	<2	<2	<2	<2	<2
2-Butanone (Methyl ethyl ketone)	μg/L	U.S.	54000/4000	SWEG/EPA	<2	<2	<2	<2	<2	<2	<2	<2
n-Butylbenzene	μg/L	U.S.			< 0.4	< 0.4	< 0.4	< 0.4	< 0.4	< 0.4	< 0.4	< 0.4
sec-Butylbenzene	μg/L	U.S.			< 0.4	< 0.4	< 0.4	< 0.4	< 0.4	< 0.4	< 0.4	< 0.4
tert-Butylbenzene	μg/L	U.S.			< 0.4	< 0.4	< 0.4	< 0.4	< 0.4	< 0.4	< 0.4	< 0.4
Carbon disulfide	μg/L	U.S.			<2	<2	<2	<2	<2	<2	<2	<2
Carbon tetrachloride	μg/L	U.S.	5	EPA	< 0.4	< 0.4	< 0.4	< 0.4	< 0.4	< 0.4	< 0.4	< 0.4
Chloroacetonitrile	μg/L	U.S.			<10	<10	<10	<10	<10	<10	<10	<10
Chlorobenzene	μg/L	U.S.	100	EPA	< 0.4	< 0.4	< 0.4	< 0.4	< 0.4	< 0.4	< 0.4	< 0.4
1-Chlorobutane (Butyl chloride)	μg/L	U.S.			< 0.4	< 0.4	< 0.4	< 0.4	< 0.4	< 0.4	< 0.4	< 0.4
Chloroethane	μg/L	U.S.			<2	<2	<2	<2	<2	<2	<2	<2
Chloroform	μg/L	U.S.	6,500/THM 80	SWEG/EPA	< 0.4	< 0.4	< 0.4	< 0.4	< 0.4	< 0.4	< 0.4	< 0.4
Chloromethane	μg/L	U.S.	30	EPA HA	<2	<2	<2	<2	<2	<2	<2	<2
2-Chlorotoluene	μg/L	U.S.	100	EPA HA	<0.4	< 0.4	< 0.4	<0.4	<0.4	<0.4	<0.4	<0.4
4-Chlorotoluene	μg/L	U.S.	100	EPA HA	<0.4	< 0.4	< 0.4	<0.4	<0.4	<0.4	<0.4	<0.4
Dibromochloromethane (DDGP)	μg/L	U.S.	THM 80	EPA	<0.4	<0.4	<0.4	<0.4	<0.4	<0.4	<0.4	<0.4
1,2-Dibromo-3-chloropropane (DBCP)	μg/L	U.S.	0.2	EPA	<2	<2	<2	<2	<2	<2	<2	<2
1,2-Dibromoethane (EDB) Dibromomethane	μg/L μg/L	U.S. U.S.	0.05	EPA	<0.4 <0.4	<0.4 <0.4	<0.4 <0.4	<0.4 <0.4	<0.4 <0.4	<0.4 <0.4	<0.4 <0.4	<0.4 <0.4
1.2-Dichlorobenzene	μg/L μg/L	U.S.	600	EPA	<0.4	<0.4	<0.4	<0.4	<0.4	<0.4	<0.4	<0.4
1.3-Dichlorobenzene	μg/L μg/L	U.S.	600	EPA HA	<0.4	<0.4	<0.4	<0.4	<0.4	<0.4	<0.4	<0.4
1.4-Dichlorobenzene	μg/L μg/L	U.S.	75	EPA HA	<0.4	<0.4	<0.4	<0.4	<0.4	<0.4	<0.4	<0.4
trans-1.4-Dichloro-2-butene	μg/L μg/L	U.S.	13	ыл	<0.4	<0.4	<0.4	<0.4	<0.4	<0.4	<0.4	<0.4
Dichlorodifluoromethane	μg/L μg/L	U.S.	1.000	EPA HA	<2	<2	<2	<2	<2	<2	<2	<2
1.1-Dichloroethane	μg/L	U.S.	2,500	22.11111	<0.4	<0.4	<0.4	<0.4	<0.4	<0.4	<0.4	<0.4
1,2-Dichloroethane	μg/L	U.S.	5	EPA	<0.4	<0.4	<0.4	<0.4	<0.4	<0.4	<0.4	<0.4
1,1-Dichloroethene	μg/L	U.S.	7	EPA	<0.4	<0.4	<0.4	<0.4	<0.4	<0.4	<0.4	<0.4
cis1,2-Dichloroethene	μg/L	U.S.	70	EPA	<0.4	<0.4	<0.4	<0.4	<0.4	<0.4	<0.4	<0.4
trans-1,2-Dichloroethene	μg/L	U.S.	100	EPA	<0.4	<0.4	< 0.4	<0.4	<0.4	<0.4	<0.4	<0.4
1,2-Dichloropropane	μg/L	U.S.	5	EPA	<0.4	<0.4	<0.4	<0.4	<0.4	<0.4	<0.4	<0.4

Mission					Soyuz 1	19/Exp. 21	ISS UL	F3/Exp. 21	ISS 20A/Exp. 22	ISS 19A	Exp. 23	Soyuz 21/Exp. 23
					WPA PWD Hot	WPA PWD	WPA PWD Hot	WPA PWD	WPA PWD	WPA PWD	WPA PWD Hot	WPA PWD Hot
						Ambient		Ambient	Ambient	Ambient		
Sample Location			Potable Water Maximum	Maximum								
Sample Description		Test	Contaminant	Contaminant	Potable Water	Potable Water	Potable Water	Potable Water				
Sample Date		Conducted	Level	Level	10/20/2009	10/20/2009	11/10/2009	11/10/2009	1/6/2010	3/3/2010	3/3/2010	3/31/2010
Analysis/Sample ID	Units	by	(MCL)	Source	20091214001	20091214002	20091130012	20091130013	20100222001	20100422006	20100422007	20100603003
1,3-Dichloropropane	μg/L	U.S.	,		<0.4	< 0.4	< 0.4	< 0.4	< 0.4	< 0.4	< 0.4	<0.4
2,2-Dichloropropane	μg/L	U.S.			< 0.4	< 0.4	< 0.4	< 0.4	< 0.4	< 0.4	< 0.4	< 0.4
1,1-Dichloropropanone	μg/L	U.S.			<2	<2	<2	<2	<2	<2	<2	<2
1,1-Dichloropropene	μg/L	U.S.			< 0.4	< 0.4	< 0.4	< 0.4	< 0.4	< 0.4	< 0.4	< 0.4
cis-1,3-Dichloropropene	μg/L	U.S.			< 0.4	< 0.4	< 0.4	< 0.4	< 0.4	< 0.4	< 0.4	< 0.4
trans-1,3-Dichloropropene	μg/L	U.S.			<2	<2	<2	<2	<2	<2	<2	<2
Diethyl ether	μg/L	U.S.			<2	<2	<2	<2	<2	<2	<2	<2
Ethylbenzene	μg/L	U.S.	700	EPA	< 0.4	< 0.4	< 0.4	< 0.4	< 0.4	< 0.4	< 0.4	<0.4
Ethyl methacrylate	μg/L	U.S.			<2	<2	<2	<2	<2	<2	<2	<2
Hexachlorobutadiene	μg/L	U.S.	1	EPA HA	<2	<2	<2	<2	<2	<2	<2	<2
Hexachloroethane	μg/L	U.S.	1	EPA HA	<2	<2	<2	<2	<2	<2	<2	<2
2-Hexanone Iodomethane	μg/L	U.S. U.S.			<2	<2	<2	<2	<2	<2	<2	<2
	μg/L		4.000	EDA DILET	<2 <0.4	<2	<2	<2 <0.4	<2 <0.4	<2 <0.4	<2 <0.4	<2 <0.4
Isopropylbenzene (Cumene) 4-Isopropyltoluene (Cymene)	μg/L	U.S. U.S.	4,000	EPA DWEL		<0.4	<0.4					
Methacrylonitrile	μg/L ug/L	U.S.			<0.4 <2	<0.4 <2	<0.4 <2	<0.4 <2	<0.4 <2	<0.4 <2	<0.4 <2	<0.4 <2
Methyl acrylate	μg/L μg/L	U.S.			<2	<2	<2	<2	<2 <2	<2	<2	<2 <2
Methyl-t-butylether (MTBE)	μg/L μg/L	U.S.			<2	<2	<2	<2	<2	<2	<2	<2
Methylene chloride (Dichloromethane)	μg/L μg/L	U.S.	15.000/5	SWEG/EPA	<0.4	<0.4	<0.4	<0.4	<0.4	<0.4	<0.4	<0.4
Methyl methacrylate	μg/L	U.S.	13,000/3	SWEG/EFA	<2	<2	<2	<2	<2	<2	<2	<2
4-Methyl-2-pentanone	μg/L	U.S.			<0.4	<0.4	<0.4	<0.4	<0.4	<0.4	<0.4	<0.4
Naphthalene	μg/L	U.S.	100	EPA HA	<0.4	<0.4	<0.4	<0.4	<0.4	<0.4	<0.4	<0.4
Nitrobenzene	μg/L	U.S.			<2	<2	<2	<2	<2	<2	<2	<2
2-Nitropropane	μg/L	U.S.			<2	<2	<2	<2	<2	<2	<2	<2
Pentachloroethane	μg/L	U.S.			<2	<2	<2	<2	<2	<2	<2	<2
Propionitrile (Ethyl cyanide)	μg/L	U.S.			<10	<10	<10	<10	<10	<10	<10	<10
n-Propylbenzene	μg/L	U.S.			< 0.4	< 0.4	< 0.4	< 0.4	< 0.4	< 0.4	< 0.4	< 0.4
Styrene	μg/L	U.S.	100	EPA	< 0.4	< 0.4	< 0.4	< 0.4	< 0.4	< 0.4	< 0.4	< 0.4
1,1,1,2-Tetrachloroethane	μg/L	U.S.	70	EPA HA	< 0.4	< 0.4	< 0.4	< 0.4	< 0.4	< 0.4	< 0.4	< 0.4
1,1,2,2-Tetrachloroethane	μg/L	U.S.	0.3	EPA HA	< 0.4	< 0.4	< 0.4	< 0.4	< 0.4	< 0.4	< 0.4	< 0.4
Tetrachloroethene	μg/L	U.S.	5	EPA	< 0.4	< 0.4	< 0.4	< 0.4	< 0.4	< 0.4	< 0.4	< 0.4
Tetrahydrofuran	μg/L	U.S.			<2	<2	<2	<2	<2	<2	<2	<2
Toluene	μg/L	U.S.	1,000	EPA	<0.4	<0.4	<0.4	<0.4	<0.4	<0.4	< 0.4	<0.4
1,2,3-Trichlorobenzene	μg/L	U.S.		TID 4	<0.4	<0.4	<0.4	<0.4	< 0.4	<0.4	<0.4	<0.4
1,2,4-Trichlorobenzene	μg/L	U.S.	70	EPA	<0.4	<0.4	<0.4	<0.4	<0.4	<0.4	<0.4	<0.4
1,1,1-Trichloroethane 1,1,2-Trichloroethane	μg/L	U.S. U.S.	200 5	EPA EPA	<0.4 <0.4	<0.4 <0.4	<0.4 <0.4	<0.4 <0.4	<0.4 <0.4	<0.4 <0.4	<0.4 <0.4	<0.4 <0.4
Trichloroethene	μg/L μg/L	U.S.	5	EPA EPA	<0.4	<0.4	<0.4	<0.4 <0.4	<0.4 <0.4	<0.4	<0.4	<0.4
Trichlorofluoromethane	μg/L μg/L	U.S.	2,000	EPA EPA HA	<0.4	<0.4 <0.4	<0.4 <0.4	<0.4 <0.4	<0.4 <0.4	<0.4	<0.4	<0.4 <0.4
1,2,3-Trichloropropane	μg/L ug/L	U.S.	2,000	EPA HA EPA HA	<0.4 <0.4	<0.4 <0.4	<0.4 <0.4	<0.4 <0.4	<0.4 <0.4	<0.4 <0.4	<0.4	<0.4 <0.4
1.2.4-Trimethylbenzene	μg/L ug/L	U.S.	40	ЕГА ПА	<0.4	<0.4 <0.4	<0.4 <0.4	<0.4 <0.4	<0.4 <0.4	<0.4	<0.4	<0.4 <0.4
1,3,5-Trimethylbenzene	μg/L μg/L	U.S.			<0.4	<0.4	<0.4	<0.4	<0.4	<0.4	<0.4	<0.4
Vinyl Acetate	μg/L μg/L	U.S.			<2	<2	<2	<2	<2	<2	<2	<0.4
Vinyl Acctate Vinyl Chloride	μg/L μg/L	U.S.	2	EPA	<2	<2	<2	<2	<2	<2	<2	<2
m&p-Xylene	μg/L μg/L	U.S.	Total Xylenes 10,000	EPA	<0.4	<0.4	<0.4	<0.4	<0.4	<0.4	<0.4	<0.4
o-Xylene	μg/L	U.S.	Total Xylenes 10,000	EPA	<0.4	<0.4	<0.4	<0.4	<0.4	<0.4	<0.4	<0.4
	F8-2		, , , , , , , , , , , , , , , , , , , ,		1	.012			.511	.012		.012

Mission					Corner 1	0/E-m 21	TCC III I	F3/Exp. 21	ICC 204/Erm 22	TCC 10.4	/Erm 22	Corner 21/Erm 22
WISSION					WPA PWD Hot	9/Exp. 21			ISS 20A/Exp. 22	ISS 19A		Soyuz 21/Exp. 23
					WPA PWD Hot	WPA PWD Ambient	WPA PWD Hot	WPA PWD Ambient	WPA PWD Ambient	WPA PWD Ambient	WPA PWD Hot	WPA PWD Hot
Sample Location			Potable Water			Ambient		Ambient	Ambient	Ambient		
Sample Librarion			Maximum	Maximum	Potable Water	Potable Water						
Sample Description		Test	Contaminant	Contaminant	Totable Water	Totable Water	Potable Water	Potable Water	Potable Water	Potable Water	Potable Water	Potable Water
Sample Date		Conducted	Level	Level	10/20/2009	10/20/2009	11/10/2009	11/10/2009	1/6/2010	3/3/2010	3/3/2010	3/31/2010
Analysis/Sample ID	Units	by	(MCL)	Source	20091214001	20091214002	20091130012	20091130013	20100222001	20100422006	20100422007	20100603003
Volatile Organics - Non-Targets (Tentatively Identi)% match quality)									
Acetaldehyde	μg/L	U.S.			not found	not found	not found	not found	not found	not found	not found	not found
Butyraldehyde (Butanal)	μg/L	U.S.			not found	not found	not found	not found	not found	not found	not found	not found
1,8-Cineole	μg/L	U.S.			not found	not found	not found	not found	not found	not found	not found	not found
Cyclohexanone	μg/L	U.S.			not found	not found	not found	not found	not found	not found	not found	not found
Difluorodimethylsilane	μg/L	U.S.			not found	not found	not found	not found	not found	not found	not found	not found
2,5-Dimethylfuran	μg/L	U.S.			not found	not found	not found	not found	not found	not found	not found	not found
2,6-Dimethyl-1,7-octadiene	μg/L	U.S.			not found	not found	not found	not found	not found	not found	not found	not found
Dimethyl sulfide (Thiobismethane)	μg/L				not found	not found	not found	not found	not found	not found	not found	not found
p-Dioxane	μg/L	U.S.			not found	not found	not found	not found	not found	not found	not found	not found
1,3-Dioxolane (Ethylene glycol formal)	μg/L	U.S.			not found	not found	not found	not found	not found	not found	not found	not found
Ethyl acetate	μg/L	U.S.			not found	not found	not found	not found	not found	not found	not found	not found
Fluorotrimethylsilane	μg/L	U.S.			not found	not found	not found	not found	not found	not found	not found	not found
2-Heptanone	μg/L	U.S.			not found	not found	not found	not found	not found	not found	not found	not found
4-Heptanone	μg/L	U.S.			not found	not found	not found	not found	not found	not found	not found	not found
3-Hexanone	μg/L	U.S.			not found	not found	not found	not found	not found	not found	not found	not found
3-Hydroxy-3-methylbutyric acid	μg/L	U.S.			not found	not found	not found	not found	not found	not found	not found	not found
Isobutyronitrile	μg/L	U.S.			not found	not found	not found	not found	not found	not found	not found	not found
Isooctanol	μg/L	U.S.			not found	not found	not found	not found	not found	not found	not found	not found
L-Menthol	μg/L	U.S.			not found	not found	not found	not found	not found	not found	not found	not found
Menthone	μg/L	U.S.			not found	not found	not found	not found	not found	not found	not found	not found
Methyl acetate	μg/L	U.S.			not found	not found	not found	not found	not found	not found	not found	not found
cis-1-Methyl-4-(1-methylethenyl)-cyclohexane	μg/L	U.S.			not found	not found	not found	not found	not found	not found	not found	not found
3-Methyl-2-pentanone	μg/L	U.S.			not found	not found	not found	not found	not found	not found	not found	not found
2-Methyl-1-propene	μg/L	U.S.			not found	not found	not found	not found	not found	not found	not found	not found
alpha-Methyl styrene	μg/L	U.S.			not found	not found	not found	not found	not found	not found	not found	not found
2-Nonanone	μg/L	U.S.			not found	not found	not found	not found	not found	not found	not found	not found
2-Octanone	μg/L	U.S.			not found	not found	not found	not found	not found	not found	not found	not found
2-Pentanone	μg/L	U.S.			not found	not found	not found	not found	not found	not found	not found	not found
1,1,2-Trichloro-1,2,2-trifluoroethane (Freon 113)	μg/L	U.S.			not found	not found	not found	not found	not found	not found	not found	not found
Trimethylsilanol	μg/L	U.S.			not found	not found	not found	not found	not found	not found	not found	not found
T												
Extractable Organics	ug/L	U.S.			<16	<16	<16	<8	<16	<16	<16	<16
Acetophenone Benzaldehyde	1.0	U.S.				_						<16 <8
Benzaidenyde Benzoic acid	μg/L	U.S. U.S.			<8 <24	<8 <24	<8 <24	<4 <12	<8 <12	<8 <24	<8 <24	<8 <24
Benzothiazole	μg/L μg/L	U.S.			<24 <8	<24 <8	<24 <8	<1z <4	<1z <8	<24 <8	<24 <8	<24 <8
Benzyl alcohol	μg/L μg/L	U.S.			<8 <8	<8 <8	<8	<4	<8	<8	<8	<8 <8
Benzyl butyl phthlate	μg/L μg/L	U.S.	7.000	EPA DWEL	<8 <8	<8 <8	<8	<4	<8	<8	<8	<8 <8
2-Butoxvethanol	μg/L μg/L	U.S.	7,000	ELADWEL	<16	<16	<16	<8	<16	<16	<16	<16
2-(2-Butoxyethoxy)ethanol	μg/L μg/L	U.S.			<16	<16	<16	<8	<16	<16	<16	<16
2-(2-Butoxyethoxy)ethyl acetate	μg/L μg/L	U.S.			<8	<8	<8	<4	<8	<8	<8	<8
n-Butylpalmitate	ug/L	U.S.			<16	<16	<16	<8	<16	<16	<16	<16
Butylated hydroxyanisole (BHA)	μg/L μg/L	U.S.			<8	<8	<8	<4	<8	<8	<8	<8
N-Butylbenzenesulfonamide	μg/L μg/L	U.S.			<8	<8	<8	<4	<8	<8	<8	<8
3-tert-Butylphenol	μg/L μg/L	U.S.			<24	<24	<24	<12	<24	<24	<24	<24
	r8~	U.S.			<8	<8	<8	<4	<8	<8		<8

Mission					Corner 1	0/E 21	ICC III I	F3/Exp. 21	ICC 204/E-m 22	TCC 10.4	/E 22	Corner 21/Erm 22
WISSION					WPA PWD Hot	9/Exp. 21 WPA PWD	WPA PWD Hot	WPA PWD	ISS 20A/Exp. 22 WPA PWD	ISS 19A WPA PWD	WPA PWD Hot	Soyuz 21/Exp. 23 WPA PWD Hot
					WPA PWD Hot	Ambient	WPA PWD Hot	Ambient	Ambient	WPA PWD Ambient	WPA PWD Hot	WPA PWD Hot
Sample Location			Potable Water			Ambient		Ambient	Ambient	Ambient		
Sample House			Maximum	Maximum	Potable Water	Potable Water						
Sample Description		Test	Contaminant	Contaminant	Totable Water	Totable Trater	Potable Water	Potable Water	Potable Water	Potable Water	Potable Water	Potable Water
Sample Date		Conducted	Level	Level	10/20/2009	10/20/2009	11/10/2009	11/10/2009	1/6/2010	3/3/2010	3/3/2010	3/31/2010
Analysis/Sample ID	Units	by	(MCL)	Source	20091214001	20091214002	20091130012	20091130013	20100222001	20100422006	20100422007	20100603003
tris-2-Chloroethyl phosphate	μg/L	U.S.			<8	<8	<8	<4	<8	<8	<8	<8
Cholesterol	μg/L	U.S.			<64	<64	<64	<32	<64	<64	<64	<64
o-Cresol (2-Methylphenol)	μg/L	U.S.			<8	<8	<8	<4	<8	<8	<8	<8
Cyclododecane	μg/L	U.S.			<8	<8	<8	<4	<8	<8	<8	<8
Decamethylcyclopentasiloxane	μg/L	U.S.			<8	<8	<8	<4	<8	<8	<8	<8
Decanoic acid	μg/L	U.S.			<16	<16	<16	<8	<16	<16	<16	<16
2,6-Di-t-butyl-1,4-benzoquinone	μg/L	U.S.			<8	<8	<8	<4	<8	<8	<8	<8
2,4-Di-t-butylphenol	μg/L	U.S.			<8	<8	<8	<4	<8	<8	<8	<8
1,4 Diacetylbenzene	μg/L	U.S.			<8	<8	<8	<4	<8	<8	<8	<8
N,N-Dibutylformamide	μg/L	U.S. U.S.	40.000/4.000	CMEC/EDA DIACE	<8	<8 <8	<8	<4	<8 <8	<8 <8	<8 <8	<8 <8
Dibutyl phthalate	μg/L		.,,	SWEG/EPA DWEL	<8		<8	<4				
Dibutylamine N.N-Diethyl-m-toluamide	μg/L μg/L	U.S. U.S.	Dialkylamines 300	SWEG	<8 <8	<8 <8	<8 <8	<4 <4	<8 <8	<8 <8	<8 <8	<8 <8
Diethylphthalate	μg/L μg/L	U.S.	30,000	EPA DWEL	<8 <8	<8	<8 <8	<4	<8	<8	<8	<8
Diethylene glycol monoethyl ether	μg/L μg/L	U.S.	30,000	EPA DWEL	<8	<8	<8	<4	<8	<8 <8	<8	<8
N,N-Diethylformamide	μg/L ug/L	U.S.			<8 <24	<8 <24	<8 <24	<12	<8 <24	<8 <24	<8 <24	<8 <24
Diiodomethane (Methyl iodide)	μg/L μg/L	U.S.			<24 <8	<24 <8	<24 <8	<1z <4	<24 <8	<24 <8	<24 <8	<24 <8
Diisopropyl adipate	μg/L μg/L	U.S.			<8	<8	<8	<4	<8	<8	<8	<8
Dimethyl phthalate	μg/L μg/L	U.S.			<8	<8	<8	<4	<8	<8	<8	<8
N,N-Dimethyl acetamide	μg/L μg/L	U.S.			<8	<8	<8	<4	<8	<8	<8	<8
N.N-Dimethylbenzylamine	μg/L		Dialkylamines 300	SWEG	<8	<8	<8	<4	<8	<8	<8	<8
N,N-Dimethylformamide	μg/L	U.S.	Drain/jamines 000	5 W.Eu	<16	<16	<16	<8	<16	<16	<16	<16
Dipropylene glycol methyl ether	μg/L	U.S.			<8	<8	<8	<4	<8	<8	<8	<8
Dodecamethylcyclohexasiloxane	μg/L	U.S.			<8	<8	<8	<4	<8	<8	<8	<8
2-Ethoxyethanol	μg/L	U.S.			<8	<8	<8	<4	<8	<8	<8	<8
2-Ethyl-1-hexanol	μg/L	U.S.			<8	<8	<8	<4	<8	<8	<8	<8
2-Ethylhexanoic acid	μg/L	U.S.			<8	<8	<8	<4	<8	<8	<8	<8
bis-2-Ethylhexyl adipate	μg/L	U.S.	400	EPA	<8	<8	<8	<4	<8	<8	<8	<8
bis-2-Ethylhexyl phthalate (Dioctyl phthlate)	μg/L	U.S.	20,000/6	SWEG/EPA	<8	<8	<8	<4	<8	<8	<8	<8
4-Ethylmorpholine	μg/L	U.S.			<8	<8	<8	<4	<8	<8	<8	<8
1-Formylpiperidine	μg/L	U.S.			<8	<8	<8	<4	<8	<8	<8	<8
Heptanoic acid	μg/L	U.S.			<8	<8	<8	<4	<8	<8	<8	<8
2-Heptanone	μg/L	U.S.			<8	<8	<8	<4	<8	<8	<8	<8
gamma-Hexalactone	μg/L	U.S.			<8	<8	<8	<4	<8	<8	<8	<8
Hexanoic acid	μg/L	U.S.			<16	<16	<16	<8	<16	<16	<16	<16
2-Hexanol	μg/L	U.S.			<8	<8	<8	<4	<8	<8	<8	<8
2-Hydroxybenzothiazole	μg/L	U.S.			<8	<8	<8	<4	<8	<8	<8	<8
Ibuprofen	μg/L	U.S.			<8	<8	<8	<4	<8	<8	<8	<8
Iodoform	μg/L	U.S.	100		<8	<8	<8	<4	<8	<8	<8	<8
Isophorone	μg/L	U.S.	100	EPA HA	<8	<8	<8	<4	<8	<8	<8	<8
4-Isopropylphenol Lauramide	μg/L	U.S.			<8	<8	<8	<4	<8	<8	<8	<8
	μg/L	U.S.			<8	<8	<8	<4	<8	<8	<8	<8
Lauric acid (Dodecanoic acid) p-Menth-1-en-8-ol (alpha-Terpineol)	μg/L	U.S. U.S.			<240	<240	<240	<120	<240	<240	<240	<240
p-Mentin-1-en-8-0i (aipna-1 erpineoi) 2-Mercaptobenzothiazole	μg/L μg/L	U.S.	30,000	SWEG	<8 <80	<8 <80	<8 <80	<4 <40	<8 <80	<8 <80	<8 <80	<8 <80
~ wici capionelizuiliazuie	րց/ Լ		30,000	SWEG				<40 <4	<80 <8	<80 <8	<80 <8	<80 <8
		HC										
2-Methyl-2,4-pentanediol 1-Methyl-2-pyrrolidinone	μg/L μg/L	U.S. U.S.			<8 <8	<8 <8	<8 <8	<4	<8	<8	<8 <8	<8

Mission		1			Sovuz 1	19/Exp. 21	111 221	F3/Exp. 21	ISS 20A/Exp. 22	ISS 19A	/Evp. 23	Sovuz 21/Exp. 23
					WPA PWD Hot	WPA PWD	WPA PWD Hot	WPA PWD	WPA PWD	WPA PWD	WPA PWD Hot	WPA PWD Hot
					WPA PWD HOL	Ambient	WPA PWD HOL	Ambient	Ambient	Ambient	WPA PWD Hot	WPA PWD HOL
Sample Location			Potable Water			Ambient		Ambient	Ambient	Ambient		
			Maximum	Maximum	Potable Water	Potable Water						
Sample Description		Test	Contaminant	Contaminant	1 otable Water	I otable Water	Potable Water	Potable Water	Potable Water	Potable Water	Potable Water	Potable Water
Sample Date		Conducted	Level	Level	10/20/2009	10/20/2009	11/10/2009	11/10/2009	1/6/2010	3/3/2010	3/3/2010	3/31/2010
Analysis/Sample ID	Units	by	(MCL)	Source	20091214001	20091214002	20091130012	20091130013	20100222001	20100422006	20100422007	20100603003
Methyl sulfone	μg/L	U.S.			111	102	66	34	34	66	87	41
2-Methyl butyric acid	μg/L	U.S.			<24	<24	<24	<12	<24	<24	<24	<24
2-Methylthiobenzothiazole	μg/L	U.S.			<8	<8	<8	<4	<8	<8	<8	<8
Monomethyl phthalate	μg/L	U.S.			<8	<8	<8	<4	<8	<8	<8	<8
Myristic acid	μg/L	U.S.			<48	<48	<48	<24	<48	<48	<48	<48
(+)-Neomenthol	μg/L	U.S.			<8	<8	<8	<4	<8	<8	<8	<8
Nicotine	μg/L	U.S.			<8	<8	<8	<4	<8	<8	<8	<8
Nonadecane	μg/L	U.S.			<8	<8	<8	<4	<8	<8	<8	<8
Nonanoic acid	μg/L	U.S.			<24	<24	<24	<12	<24	<24	<24	<24
1-Octadecanol	μg/L	U.S.			<24	<24	<24	<12	<24	<24	<24	<24
Octamethylcyclotetrasiloxane	μg/L	U.S.			<8	<8	<8	<4	<8	<8	<8	<8
Octanoic acid	μg/L	U.S.			<16	<16	<16	<8	<16	<16	<16	<16
4-tert-Octylphenol	μg/L	U.S.			<8	<8	<8	<4	<8	<8	<8	<8
Oleic acid	μg/L	U.S.			<80	<80	<80	<40	<80	<80	<80	<80
Oxindole	μg/L	U.S.			<8	<8	<8	<4	<8	<8	<8	<8
Palmitic acid	μg/L	U.S.			<240	<240	<240	<120	<240	<240	<240	<240
Palmitoleic acid	ug/L	U.S.			NA	NA	NA	NA	NA	NA	NA	NA
Pentacosane	μg/L	U.S.			<8	<8	<8	<4	<8	<8	<8	<8
sec-Phenethyl alcohol	μg/L	U.S.			<8	<8	<8	<4	<8	<8	<8	<8
Phenol	μg/L	U.S.	4,000	SWEG	<8	<8	<8	<4	<8	<8	<8	<8
2-Phenoxyethanol	μg/L	U.S.			<8	<8	<8	<4	<8	<8	<8	<8
N-Phenyl-2-naphthylamine	μg/L	U.S.	260,000	SWEG	<8	<8	<8	<4	<8	<8	<8	<8
2-Phenyl-2-propanol	μg/L	U.S.			<8	<8	<8	<4	<8	<8	<8	<8
2-Phenylacetic acid	μg/L	U.S.			<32	<32	<32	<16	<32	<32	<32	<32
Phenethyl alcohol	μg/L	U.S.			<8	<8	<8	<4	<8	<8	<8	<8
2-Phenylphenol	μg/L	U.S.			<8	<8	<8	<4	<8	<8	<8	<8
Salicyclic Acid	μg/L	U.S.			<64	<64	<64	<32	<64	<64	<64	<64
trans-Squalene	μg/L	U.S.			<16	<16	<16	<8	<16	<16	<16	<16
Stearic acid	μg/L	U.S.			NA	NA	NA	NA	NA	NA	NA	NA
1-Tetradecanol	μg/L	U.S.			<8	<8	<8	<4	<8	<8	<8	<8
Tetramethylsuccinonitrile	μg/L	U.S.			<8	<8	<8	<4	<8	<8	<8	<8
Tetramethyl thiourea	μg/L	U.S.			<8	<8	<8	<4	<8	<8	<8	<8
Tetramethylurea	μg/L	U.S.			<8	<8	<8	<4	<8	<8	<8	<8
Thymol	μg/L	U.S.			<8	<8	<8	<4	<8	<8	<8	<8
1,3,5-Triallyl-1,3,5-triazine-2,4,6(1H,3H,5H)-trione	μg/L	U.S.			<8	<8	<8	<4	<8	<8	<8	<8
Tributylamine	μg/L	U.S.	Trialkylamines 400	SWEG	<8	<8	<8	<4	<8	<8	<8	<8
Tributyl phosphate	μg/L	U.S.			<8	<8	<8	<4	<8	<8	<8	<8
Triethyl phosphate	μg/L	U.S.			<16	<16	<16	<8	<16	<16	<16	<16
2,2,4-Trimethyl-1,3-pentanediol diisobutyrate	μg/L	U.S.			<16	<16	<16	<8	<16	<16	<16	<16
Tripropylene glycol monomethyl ether	μg/L	U.S.			<8	<8	<8	<4	<8	<8	<8	<8
Undecanoic acid	μg/L	U.S.			<48	<48	<48	<24	<48	<48	<48	<48
2-Undecanone	μg/L	U.S.			<8	<8	<8	<4	<8	<8	<8	<8
Valeric acid (Pentanoic acid)	μg/L	U.S.			<48	<48	<48	<24	<48	<48	<48	<48
Vanillin	μg/L	U.S.			<16	<16	<16	<8	<16	<16	<16	<16

Mission				I	C 1	0/E 21	TCCTIT	F2/F 21	TCC 20 A /E 22	TCC 10.4	/E 22	C 21/E 22
1411991011					WPA PWD Hot	9/Exp. 21 WPA PWD	WPA PWD Hot	F3/Exp. 21 WPA PWD	ISS 20A/Exp. 22 WPA PWD	ISS 19A WPA PWD	WPA PWD Hot	Soyuz 21/Exp. 23 WPA PWD Hot
					WPA PWD Hot	Ambient	WPA PWD Hot	Ambient	Ambient	WPA PWD Ambient	WPA PWD Hot	WPA PWD Hot
Sample Location			Potable Water			Ambient		Ambient	Ambient	Ambient		
Dampio 200ation			Maximum	Maximum	Potable Water	Potable Water						
Sample Description		Test	Contaminant	Contaminant	Totable Trater	Totable Trater	Potable Water	Potable Water	Potable Water	Potable Water	Potable Water	Potable Water
Sample Date		Conducted	Level	Level	10/20/2009	10/20/2009	11/10/2009	11/10/2009	1/6/2010	3/3/2010	3/3/2010	3/31/2010
Analysis/Sample ID	Units	by	(MCL)	Source	20091214001	20091214002	20091130012	20091130013	20100222001	20100422006	20100422007	20100603003
Acid Extractables-EPA 625 List												
4-Chloro-3-methylphenol	μg/L	U.S.			<16	<16	<16	<8	<16	<16	<16	<16
2-Chlorophenol	μg/L	U.S.	40	EPA HA	<16	<16	<16	<8	<16	<16	<16	<16
2,4-Dichlorophenol	μg/L	U.S.	20	EPA HA	<16	<16	<16	<8	<16	<16	<16	<16
2,4-Dimethylphenol	μg/L	U.S.			<16	<16	<16	<8	<16	<16	<16	<16
2,4-Dinitrophenol	μg/L	U.S.			<16	<16	<16	<8	<16	<16	<16	<16
2-Methyl-4,6-dinitrophenol	μg/L	U.S.			<16	<16	<16	<8	<16	<16	<16	<16
2-Nitrophenol	μg/L	U.S.			<16	<16	<16	<8	<16	<16	<16	<16
4-Nitrophenol	μg/L	U.S.	60	EPA HA	<16	<16	<16	<8	<16	<16	<16	<16
Pentachlorophenol	μg/L	U.S.	1	EPA	<16	<16	<16	<8	<16	<16	<16	<16
Phenol	μg/L	U.S.	4,000/2,000	SWEG/EPA HA	<8	<8	<8	<4	<8	<8	<8	<8
2,4,5-Trichlorophenol	μg/L	U.S.			<16	<16	<16	<8	<16	<16	<16	<16
2,4,6-Trichlorophenol	μg/L	U.S.	10	EPA DWEL	<16	<16	<16	<8	<16	<16	<16	<16
4-Methylphenol	μg/L	U.S.			NA	NA	NA	NA	<4	<8	<8	<8
Base/Neutral Extractables - EPA 625 List												
Benzidine	μg/L	U.S.			<16	<16	<16	<8	<16	<16	<16	<16
3,3-Dichlorobenzidine	μg/L	U.S.			<16	<16	<16	<8	<16	<16	<16	<16
bis-(2-Ethylhexyl)phthalate	μg/L	U.S.	20,000/6	SWEG/EPA	<8	<8	<8	<4	<8	<8	<8	<8
Benzyl butyl phthalate	μg/L	U.S.	7,000	EPA DWEL	<8	<8	<8	<4	<8	<8	<8	<8
Dibutylphthalate	μg/L	U.S.	40,000/4,000	SWEG/EPA DWEL	<8	<8	<8	<4	<8	<8	<8	<8
Diethylphthalate	μg/L	U.S.	30,000	EPA DWEL	<8	<8	<8	<4	<8	<8	<8	<8
Dimethylphthalate	μg/L	U.S.			<8	<8	<8	<4	<8	<8	<8	<8
Di-n-octyl phthalate	μg/L	U.S.			<16	<16	<16	<8	<16	<16	<16	<16
N-Nitrosodimethylamine	μg/L	U.S.			<16	<16	<16	<8	<16	<16	<16	<16
N-Nitrosodiphenylamine	μg/L	U.S.			<16	<16	<16	<8	<16	<16	<16	<16
N-Nitrosodi-n-propylamine	μg/L	U.S.	100	EDA DILET	<16	<16	<16	<8	<16	<16	<16	<16
2,4-Dinitrotoluene	μg/L	U.S.	100	EPA DWEL	<16	<16	<16	<8	<16	<16	<16	<16
2,6-Dinitrotoluene	μg/L	U.S.	40	EPA DWEL	<16	<16	<16	<8	<16	<16	<16	<16
Isophorone Nitrobenzene	μg/L	U.S. U.S.	100	EPA HA	<8 <16	<8	<8 <16	<4	<8	<8	<8	<8
Acenaphthene	μg/L	U.S.	2000	EPA DWEL	<16 <16	<16	<16 <16	<8 <8	<16 <16	<16 <16	<16 <16	<16 <16
Acenaphthylene	μg/L	U.S.	2000	EPA DWEL	<16	<16 <16	<16	<8	NA	<16	<16	<16
Anthracene	μg/L μg/L	U.S.	10.000	EPA DWEL	<16 <16	<16 <16	<16 <16	<8 <8	NA <16	<16 <16	<16 <16	<16 <16
Benzo(a)anthracene	μg/L μg/L	U.S.	10,000	EPA DWEL	<16 <16	<16 <16	<16 <16	<8 <8	<16 <16	<16	<16 <16	<16 <16
Benzo(a)pyrene	μg/L μg/L	U.S.	0.2	EPA	<10	<10	<10	< 8	<10	<10	<10	<10
Benzo(b)fluoranthene	μg/L μg/L	U.S.	U.£	EFA	<10 <8	<10 <8	<10 <8	<3 <4	<10 <8	<10 <8	<10 <8	<10 <8
Benzo(ghi)perylene	μg/L μg/L	U.S.			<10	<10	<10	<5	<10	<10	<10	<10
Benzo(k)fluoroanthene	μg/L μg/L	U.S.			<10 <8	<8	<8	<4	<8	<8	<8	<8
Chrysene	μg/L μg/L	U.S.			<8 <20	<8 <20	<8 <20	<10	<20	<8 <20	<8 <20	<8
Dibenzo(a,h)anthracene	μg/L ug/L	U.S.			<20 <10	<10	<20 <10	<10 <5	<10	<10	<20 <10	<10
Fluoranthene	μg/L μg/L	U.S.			<8	<8	<8	<4	<8	<8	<8	<10 <8
Fluorene	μg/L μg/L	U.S.	1,000	EPA DWEL	<16	<16	<16	<8	<16	<16	<16	<16
Indeno(1,2,3-cd)pyrene	μg/L μg/L	U.S.	1,000	ELA DVVEL	<10	<10	<10	< 8	<10	<10	<10	<10
Naphthalene	μg/L μg/L	U.S.	100	EPA HA	<40	<40	<40	<20	<40	<40	<10 <40	<10 <40
Phenanthrene	μg/L μg/L	U.S.	100	EIAHA	<40 <8	<8	<40 <8	<4	<8	<8	<8	<8
Pyrene	μg/L μg/L	U.S.			<8	<8	<8	<4	<8	<8	<8	<8
- J	μg/L μg/L	U.S.		1	<16	<16	<16	<8	<16	<16	<16	<16

In an			1									
Mission					Soyuz 1	19/Exp. 21	ISS UL	F3/Exp. 21	ISS 20A/Exp. 22	ISS 19A	Exp. 23	Soyuz 21/Exp. 23
					WPA PWD Hot	WPA PWD Ambient	WPA PWD Hot	WPA PWD Ambient	WPA PWD Ambient	WPA PWD Ambient	WPA PWD Hot	WPA PWD Hot
Sample Location			Potable Water Maximum	Maximum	Potable Water	Potable Water						
Sample Description Sample Date		Test Conducted	Contaminant	Contaminant Level	10/20/2009	10/20/2009	Potable Water 11/10/2009	Potable Water 11/10/2009	Potable Water 1/6/2010	Potable Water 3/3/2010	Potable Water 3/3/2010	Potable Water 3/31/2010
			Level									
Analysis/Sample ID	Units	by	(MCL)	Source	20091214001	20091214002	20091130012	20091130013	20100222001	20100422006	20100422007	20100603003
bis(2-Chloroethoxy) methane	μg/L	U.S.			<16	<16	<16	<8	<16	<16	<16	<16
bis(2-Chloroisopropyl) ether	μg/L	U.S.	300	EPA HA	<16	<16	<16	<8	<16	<16	<16	<16
4-Bromophenyl phenyl ether	μg/L	U.S.			<16	<16	<16	<8	<16	<16	<16	<16
4-Chlorophenyl phenyl ether	μg/L	U.S.			<16	<16	<16	<8	<16	<16	<16	<16
2-Chloronaphthalene	μg/L	U.S.			<16	<16	<16	<8	<16	<16	<16	<16
1,2-Dichlorobenzene	μg/L	U.S.	600	EPA	<16	<16	<16	<8	<16	<16	<16	<16
1,3-Dichlorobenzene	μg/L	U.S.	600	EPA HA	<16	<16	<16	<8	<16	<16	<16	<16
1,4-Dichlorobenzene	μg/L	U.S.	75	EPA	<16	<16	<16	<8	<16	<16	<16	<16
Hexachlorobenzene	μg/L	U.S.	30	EPA DWEL	<16	<16	<16	<8	<16	<16	<16	<16
Hexachlorobutadiene	μg/L	U.S.	1	EPA HA	<16	<16	<16	<8	<16	<16	<16	<16
Hexachlorocyclopentadiene	μg/L	U.S.	50	EPA	<16	<16	<16	<8	<16	<16	<16	<16
Hexachloroethane	μg/L	U.S.	1	EPA HA	<16	<16	<16	<8	<16	<16	<16	<16
1,2,4-Trichlorobenzene	μg/L	U.S.	70	EPA	<16	<16	<16	<8	<16	<16	<16	<16
Alcohols (DAI/GC/MS)												
1-Butanol	μg/L	U.S.			<100	<100	<100	<100	<100	<100	<100	<100
2-Butanol	μg/L	U.S.			<100	<100	<100	<100	<100	<100	<100	<100
Ethanol	μg/L	U.S.			<100	<100	<100	<100	<100	<100	<100	<100
Methanol	μg/L	U.S.	40000	SWEG	<100	<100	<100	<100	<100	<100	<100	<100
2-Methyl-1-butanol	μg/L	U.S.			<100	<100	<100	<100	<100	<100	<100	<100
2-Methyl-2-butanol	μg/L	U.S.			<100	<100	<100	<100	<100	<100	<100	<100
3-Methyl-1-butanol (Isopentanol)	μg/L	U.S.			<100	<100	<100	<100	<100	<100	<100	<100
2-Methyl-1-propanol	μg/L	U.S.			<100	<100	<100	<100	<100	<100	<100	<100
2-Methyl-2-propanol	μg/L	U.S.			<100	<100	<100	<100	<100	<100	<100	<100
1-Pentanol (Amyl alcohol)	μg/L	U.S.			<100	<100	<100	<100	<100	<100	<100	<100
2-Pentanol (sec-Amyl alcohol)	μg/L	U.S.			<100	<100	<100	<100	<100	<100	<100	<100
3-Pentanol	μg/L	U.S.			<100	<100	<100	<100	<100	<100	<100	<100
1-Propanol	μg/L	U.S.			<100	<100	<100	<100	<100	<100	<100	<100
2-Propanol (Isopropanol)	μg/L	U.S.			<100	<100	<100	<100	<100	<100	<100	<100
Glycols (DAI/GC/MS)												
1,2-Ethanediol (Ethylene glycol)	μg/L	U.S.	12000/4000	MORD/SWEG	<1000	<1000	<1000	<1000	<1000	<1000	<1000	<1000
1,2-Propanediol (Propylene glycol)	μg/L	U.S.	1700000	SWEG	<500	< 500	< 500	< 500	<500	< 500	<500	< 500
Glycerol (LC/MS/MS)												
Glycerol (1,2,3-Propanetriol)	μg/L	U.S.			NA	NA	NA	NA	NA	NA	NA	NA
Silanols (GC/MS & LC/MS/MS) (R&D Method -			not available)									
Dimethylsilanediol (DMSD)	μg/L	U.S.			NA	NA	NA	<400	<400	<400	<400	NA
Carboxylates (CE)												
Acetate	μg/L	U.S.			<125	<125	<125	<125	<125	<125	<125	<125
Formate	μg/L	U.S.	2,500,000	SWEG	<125	<125	<125	<125	<125	<125	<125	<125
Glycolate	μg/L	U.S.			<125	<125	<125	<125	<125	<125	<125	<125
Glyoxylate	μg/L	U.S.			<125	<125	<125	<125	<125	<125	<125	<125
Lactate	μg/L	U.S.			<1000	<1000	<1000	<1000	<1000	<1000	<1000	<1000
Oxalate	μg/L	U.S.			<125	<125	<125	<125	<125	<125	<125	<125
Propionate	μg/L	U.S.			<125	<125	<125	<125	<125	<125	<125	<125

Mission					Soyuz 1	9/Exp. 21	ISS UL	F3/Exp. 21	ISS 20A/Exp. 22	ISS 19A/	Exp. 23	Soyuz 21/Exp. 23
					WPA PWD Hot	WPA PWD	WPA PWD Hot	WPA PWD	WPA PWD	WPA PWD	WPA PWD Hot	WPA PWD Hot
						Ambient		Ambient	Ambient	Ambient		
Sample Location Sample Description		Test	Potable Water Maximum Contaminant	Maximum Contaminant	Potable Water	Potable Water	Potable Water	Potable Water				
Sample Date		Conducted	Level	Level	10/20/2009	10/20/2009	11/10/2009	11/10/2009	1/6/2010	3/3/2010	3/3/2010	3/31/2010
Analysis/Sample ID	Units	by	(MCL)	Source	20091214001	20091214002	20091130012	20091130013	20100222001	20100422006	20100422007	20100603003
Aldehydes												
Formaldehyde	μg/L	U.S.	12,000/1,000	SWEG/EPA HA	<2	2	<2	<2	4	5	2	4
Amines (CE)												
Ethylamine	μg/L	U.S.	Monoalkylamines 2000	SWEG	<125	<125	<125	<125	<125	<125	<125	<125
Methylamine	μg/L	U.S.	Monoalkylamines 2000	SWEG	<125	<125	<125	<125	<125	<125	<125	<125
n-Propylamine	μg/L	U.S.	Monoalkylamines 2000	SWEG	<125	<125	<125	<125	<125	<125	<125	<125
Trimethylamine	μg/L	U.S.	Trialkylamines 400	SWEG	<125	<125	<125	<125	<125	<125	<125	<125
						·				·		
Non-volatiles (LC/UV-VIS)												
Urea	μg/L	U.S.			<800	<800	<800	<800	<800	<800	<800	<800
Caprolactam	μg/L	U.S.	100,000	SWEG	<8	<8	<8	<4	<8	<8	<8	<8
Organic Carbon Recovery	percent	U.S.			9.60	12.05	7.59	6.33	6.46	24.73	18.64	4.73
Unaccounted Organic Carbon	mg/L	U.S.			0.27	0.19	0.21	0.13	0.15	0.13	0.13	0.24

Mission		ICC III E	4/E 22	ı	C 22/E 2/	1	C 2	2/E 25	TCC 10.4	/F 22	C 22/E 24
1411221011		WPA PWD	4/Exp. 23 WPA PWD	WPA PWD	Soyuz 22/Exp. 24 WPA PWD Hot	WPA PWD		3/Exp. 25 WPA PWD Hot	PWD Aux Port	/Exp. 23 WPA RIP	Soyuz 22/Exp. 24 WPA RIP
		Ambient	Ambient	Ambient	WYA PWD Hot	MPA PWD Ambient	WPA PWD Hot	WPA PWD Hot	P WD Aux Port	WPA KIP	WPA KIP
Sample Location		Ambient	Ambient	Ambient		Ambient			•		
_									Processed	Processed	
Sample Description		Potable Water	Potable Water	Potable Water	Potable Water	Potable Water	Potable Water	Potable Water	Water	Water	Processed Water
Sample Date		4/26/2010	5/18/2010	7/14/2010	8/25/2010	9/15/2010	10/19/2010	11/23/2010	2/25/2010	2/25/2010	7/29/2010
Analysis/Sample ID	Units	20100527006	20100527007	20100926001	20100926002	20100926003	20101129001	20101129002	20100422009	20100422008	20100926004
Physical Characteristics									•		
pH	pH units	7.35	7.25	7.95	7.88	7.56	6.86	6.95	6.05	6.18	7.31
Conductivity	μS/cm	2	1	1	1	1	4	2	3	3	3
Turbidity	NTU	<0.1	<0.1	<0.1	<0.1	<0.1	NA.	<0.1	0.1	<0.1	NA.
Total Solids	mg/L	<5	<5	NA	NA	NA	NA	NA	<5	<5	NA
Iodine (LCV)											
Total I	mg/L	<0.05	< 0.05	< 0.05	< 0.05	< 0.05	< 0.05	<0.05	2.24	2.72	2.63
Iodine	mg/L	< 0.05	< 0.05	< 0.05	< 0.05	< 0.05	< 0.05	< 0.05	1.11	2.03	2.05
Iodide	mg/L	< 0.05	< 0.05	< 0.05	< 0.05	< 0.05	< 0.05	< 0.05	1.13	0.69	0.58
Anions (IC/ISE)											
Bromide	mg/L	< 0.5	< 0.5	< 0.5	< 0.5	< 0.5	< 0.5	< 0.5	< 0.5	< 0.5	< 0.5
Chloride	mg/L	< 0.15	< 0.15	< 0.15	< 0.15	< 0.15	< 0.15	< 0.15	< 0.15	< 0.15	< 0.15
Fluoride	mg/L	<0.1	<0.1	<0.1	<0.1	<0.1	< 0.1	<0.1	<0.1	<0.1	<0.1
Nitrate as Nitrogen (NO3-N)	mg/L	<0.11	<0.11	<0.11	<0.11	<0.11	<0.11	<0.11	<0.11	<0.11	<0.11
Nitrite as Nitrogen (NO2-N)	mg/L	NA	NA	NA	NA	NA	NA	NA	NA	NA	NA
Phosphate as P (PO4-P)	mg/L	<0.24	<0.24	<0.02	< 0.02	<0.02	<0.01	<0.01	<0.24	<0.24	<0.02
Sulfate	mg/L	<0.75	< 0.75	<0.75	< 0.75	< 0.75	< 0.75	<0.75	< 0.75	< 0.75	< 0.75
Cations (IC)											
Ammonia as Nitrogen (NH3-N)	mg/L	< 0.002	< 0.002	< 0.002	< 0.002	< 0.002	< 0.002	< 0.002	< 0.002	< 0.002	< 0.002
Lithium	mg/L	< 0.002	<0.002	<0.002	<0.002	<0.002	<0.002	<0.002	<0.002	< 0.002	<0.002
Metals (ICP/MS)									\ \		
Calcium	mg/L	< 0.01	< 0.01	< 0.02	0.03	0.03	0.03	< 0.01	< 0.01	< 0.01	0.03
Magnesium	mg/L	< 0.01	< 0.01	< 0.02	< 0.02	< 0.02	< 0.01	< 0.01	< 0.01	< 0.01	< 0.02
Potassium	mg/L	< 0.01	< 0.01	< 0.02	< 0.02	< 0.02	0.03	< 0.01	0.01	< 0.01	< 0.02
Sodium	mg/L	< 0.01	0.02	< 0.02	< 0.02	< 0.02	0.02	< 0.01	< 0.01	< 0.01	< 0.02
Aluminum	μg/L	<2	<2	<4	<4	<4	<2	<2	<2	<2	<4
Antimony	μg/L	<2	<2	<4	<4	<4	<2	<2	<2	<2	<4
Arsenic	μg/L	<1	<1	<2	<2	<2	<1	<1	<1	<1	<2
Barium	μg/L	<1	<1	<2	<2	15	<1	<1	<1	<1	<2
Beryllium	μg/L	<1	<1	<2	<2	<2	<1	<1	<1	<1	<2
Cadmium	μg/L	<1	<1	<2	<2	<2	<1	<1	<1	<1	<2
Chromium	μg/L	<5	<5	<10	<10	<10	<5	<5	<5	<5	<10
Cobalt	μg/L	<1	<1	<2	<2	<2	<1	<1	<1	<1	<2
Copper	μg/L	<1	<1 <5	<2	<2	<2	<1	<1	<1 <5	<1	<2
Iron	μg/L	<5		<10	<10 <2	<10	<5	<5		<5	<10
Lead	μg/L ug/I	<1 <1	<1 <1	<2 <2	<2 <2	<2 <2	<1 <1	<1 <1	<1 <1	<1 <1	<2 <2
Manganese Mercury	μg/L μg/L	<0.5	<0.5	<z <1</z 	<z <1</z 	<z <1</z 	<0.5	<0.5	<0.5	<0.5	<z <1</z
Molybdenum	μg/L μg/L	<1	<0.5	<2	<2	<2	<1	<1	<1	<1	<2
Nickel	μg/L μg/L	27	29	29	17	4	17	23	270	42	53
Selenium	μg/L	<1	<1	<2	<2	<2	<1	<1	<1	<1	<2
Silver	μg/L	<2	<2	<4	<4	<4	<2	<2	<2	<2	<4

Mission		ISS III E	4/Exp. 23		Soyuz 22/Exp. 24	1	Sovna 2	3/Exp. 25	ISS 10A	/Exp. 23	Soyuz 22/Exp. 24
		WPA PWD	WPA PWD	WPA PWD	WPA PWD Hot	WPA PWD		WPA PWD Hot	PWD Aux Port	WPA RIP	WPA RIP
		Ambient	Ambient	Ambient	WPA PWD Hot	Ambient	WPA PWD Hot	WPA PWD HOL	PWD Aux Pon	WPA KIP	WPA KIP
Sample Location		rimbient	rimbient	rimbient		rimbient			N		
-									Processed	Processed	
Sample Description		Potable Water	Potable Water	Potable Water	Potable Water	Potable Water	Potable Water	Potable Water	Water	Water	Processed Water
Sample Date		4/26/2010	5/18/2010	7/14/2010	8/25/2010	9/15/2010	10/19/2010	11/23/2010	2/25/2010	2/25/2010	7/29/2010
Analysis/Sample ID	Units	20100527006	20100527007	20100926001	20100926002	20100926003	20101129001	20101129002	20100422009	20100422008	20100926004
Zinc	μg/L	<1	<1	<2	<2	2	<1	<1	2	<1	<2
Silicon (ICP/MS)									` 		
Silicon (ICP/MS)		NA	NA	1530	2130	1210	2360	223	NA	NA	1770
,									1		
Total Organic Carbon (Sievers)	İ										
Total Inorganic Carbon	mg/L	0.87	0.84	1.17	0.93	1.05	0.98	0.73	1.21	1.31	0.94
Total Organic Carbon	mg/L	0.16	0.22	1.51	2.19	1.11	2.51	0.15	0.24	0.18	2.18
Voletile Outenier									1		
Volatile Organics Acetone	μg/L	11	<2	<2	<2	<2	40	<2	25	<2	<6
Acryloniltrile	μg/L	<2	<2	<2	<2	<2	<6	<2	<2	<2	<6
Allyl chloride (3-Chloropropene)	μg/L	<2	<2	<2	<2	<2	<6	<2	<2	<2	<6
Benzene	μg/L	<0.4	<0.4	<0.4	<0.4	<0.4	<1.2	<0.4	< 0.4	< 0.4	<1.2
Bromobenzene	μg/L	<0.4	<0.4	<0.4	<0.4	<0.4	<1.2	<0.4	< 0.4	<0.4	<1.2
Bromochloromethane	μg/L	<4	<4	<4	<4	<4	<12	<4	<4	<4	<12
Bromodichloromethane	μg/L	< 0.4	< 0.4	< 0.4	< 0.4	< 0.4	<1.2	< 0.4	< 0.4	< 0.4	<1.2
Bromoform	μg/L	<2	<2	<2	<2	<2	<6	<2	<2	<2	<6
Bromomethane	μg/L	<2	<2	<2	<2	<2	<6	<2	<2	<2	<6
2-Butanone (Methyl ethyl ketone)	μg/L	<2	<2	<2	4	4	24	<2	<2	<2	<6
n-Butylbenzene	μg/L	< 0.4	< 0.4	< 0.4	< 0.4	< 0.4	<1.2	< 0.4	< 0.4	< 0.4	<1.2
sec-Butylbenzene	μg/L	< 0.4	< 0.4	< 0.4	< 0.4	< 0.4	<1.2	< 0.4	< 0.4	< 0.4	<1.2
tert-Butylbenzene	μg/L	< 0.4	< 0.4	< 0.4	< 0.4	< 0.4	<1.2	< 0.4	< 0.4	< 0.4	<1.2
Carbon disulfide	μg/L	<2	<2	<2	<2	<2	<6	<2	<2	<2	<6
Carbon tetrachloride	μg/L	< 0.4	< 0.4	< 0.4	< 0.4	< 0.4	<1.2	<0.4	< 0.4	< 0.4	<1.2
Chloroacetonitrile	μg/L	<10	<10	<10	<10	<10	<30	<10	<10	<10	<30
Chlorobenzene	μg/L	< 0.4	< 0.4	< 0.4	< 0.4	< 0.4	<1.2	< 0.4	< 0.4	< 0.4	<1.2
1-Chlorobutane (Butyl chloride)	μg/L	< 0.4	< 0.4	< 0.4	< 0.4	< 0.4	<1.2	< 0.4	< 0.4	< 0.4	<1.2
Chloroethane	μg/L	<2	<2	<2	<2	<2	<6	<2	<2	<2	<6
Chloroform	μg/L	< 0.4	< 0.4	< 0.4	< 0.4	< 0.4	<1.2	< 0.4	< 0.4	< 0.4	<1.2
Chloromethane	μg/L	<2	<2	<2	<2	<2	<6	<2	<2	<2	<6
2-Chlorotoluene	μg/L	<0.4	<0.4	<0.4	<0.4	<0.4	<1.2	<0.4	<0.4	<0.4	<1.2
4-Chlorotoluene	μg/L	<0.4	<0.4	<0.4	<0.4	<0.4	<1.2	<0.4	<0.4	<0.4	<1.2
Dibromochloromethane (PD CD)	μg/L	<0.4	<0.4	<0.4	<0.4	<0.4	<1.2	<0.4	<0.4	<0.4	<1.2
1,2-Dibromo-3-chloropropane (DBCP)	μg/L	<2	<2	<2	<2	<2	<6	<2	<2	<2	<6
1,2-Dibromoethane (EDB)	μg/L	<0.4	<0.4	<0.4	<0.4	<0.4	<1.2	<0.4	<0.4	<0.4	<1.2
Dibromomethane	μg/L	<0.4	<0.4	<0.4	<0.4	<0.4	<1.2	<0.4	<0.4	<0.4	<1.2
1,2-Dichlorobenzene	μg/L	<0.4	<0.4	<0.4	<0.4	<0.4	<1.2	<0.4	<0.4	<0.4	<1.2
1,3-Dichlorobenzene	μg/L	<0.4	<0.4	<0.4	<0.4	<0.4	<1.2	<0.4	<0.4	<0.4	<1.2
1,4-Dichlorobenzene	μg/L	<0.4	<0.4	<0.4	<0.4	<0.4	<1.2	<0.4	<0.4	<0.4	<1.2
trans-1,4-Dichloro-2-butene Dichlorodifluoromethane	μg/L	<0.4	<0.4	<0.4	<0.4	<0.4	<1.2	<0.4	<0.4	<0.4	<1.2
1.1-Dichloroethane	μg/L	<2	<2	<2	<2	<2	<6	<2	<2	<2	<6
1,2-Dichloroethane	μg/L	<0.4	<0.4	<0.4	<0.4	<0.4	<1.2	<0.4	< 0.4	<0.4	<1.2
1,1-Dichloroethane	μg/L	<0.4	<0.4	<0.4	<0.4	<0.4	<1.2	<0.4	<0.4	<0.4	<1.2
cis1,2-Dichloroethene	μg/L	<0.4	<0.4	<0.4 <0.4	<0.4	<0.4 <0.4	<1.2 <1.2	<0.4	<0.4	<0.4	<1.2 <1.2
trans-1,2-Dichloroethene	μg/L	<0.4	<0.4		<0.4	<0.4 <0.4		<0.4	<0.4	<0.4	
1,2-Dichloropropane	μg/L μg/L	<0.4 <0.4	<0.4 <0.4	<0.4 <0.4	<0.4 <0.4	<0.4 <0.4	<1.2 <1.2	<0.4 <0.4	<0.4 <0.4	<0.4 <0.4	<1.2 <1.2

Mission		ISS ULF	4/Exp. 23		Soyuz 22/Exp. 24	ı	Soyuz 2	3/Exp. 25	ISS 19A	/Exp. 23	Soyuz 22/Exp. 24
		WPA PWD	WPA PWD	WPA PWD	WPA PWD Hot	WPA PWD		WPA PWD Hot	PWD Aux Port	WPA RIP	WPA RIP
		Ambient	Ambient	Ambient		Ambient					
Sample Location									,		
Sample Description		Potable Water	Datable Water	Detable Western	Datable Water	Datable Water	Datable Mister	Datable Water	Processed Water	Processed	D
Sample Description Sample Date		4/26/2010	Potable Water 5/18/2010	Potable Water 7/14/2010	Potable Water 8/25/2010	Potable Water 9/15/2010	Potable Water 10/19/2010	Potable Water 11/23/2010	2/25/2010	Water 2/25/2010	Processed Water 7/29/2010
Analysis/Sample ID	Units	20100527006	20100527007	20100926001	20100926002	20100926003	20101129001	20101129002	20100422009	20100422008	20100926004
1,3-Dichloropropane		<0.4	<0.4	<0.4	<0.4	<0.4	<1.2	<0.4	< 0.4	< 0.4	<1.2
2,2-Dichloropropane	μg/L μg/L	<0.4	<0.4	<0.4	<0.4	<0.4	<1.2	<0.4	<0.4	<0.4	<1.2
1,1-Dichloropropanoe		<0.4 <2	<0.4	<0.4	<0.4 <2	<0.4	<6	<0.4 <2	<0.4	<0.4	<6
1,1-Dichloropropene	μg/L	<0.4	<0.4	<0.4	<0.4		<1.2	<0.4	<0.4	<0.4	<1.2
cis-1,3-Dichloropropene	μg/L	<0.4	<0.4	<0.4	<0.4	<0.4 <0.4	<1.2	<0.4	<0.4	<0.4	<1.2
trans-1,3-Dichloropropene	μg/L μg/L	<0.4 <2	<0.4	<0.4	<0.4 <2	<0.4	<6	<0.4 <2	<0.4	<0.4	<6
		<2 <2	<2 <2	<2 <2	<2 <2	<2 <2		<2 <2	<2 <2	<2 <2	<6
Diethyl ether Ethylbenzene	μg/L	<2 <0.4	<2. <0.4	<2 <0.4	<2 <0.4	<2. <0.4	<6 <1.2	<2 <0.4	<2. <0.4	<2 <0.4	<0 <1.2
Ethyl methacrylate	μg/L		<0.4 <2		<0.4 <2	<0.4 <2	<1.2 <6	<0.4 <2	<0.4 <2	<0.4	
	μg/L	<2		<2	+						<6
Hexachlorobutadiene Hexachloroethane	μg/L	<2	<2	<2 <2	<2 <2	<2 <2	<6	<2 <2	<2 <2	<2 <2	<6 <6
	μg/L	<2	<2 <2	<2 <2	<2 <2	<2 <2	<6	<2 <2	<2 <2	<2 <2	-
2-Hexanone	μg/L	<2			+		<6		•		<6
Iodomethane Isopropylbenzene (Cumene)	μg/L	<2 <0.4	<2 <0.4	<2 <0.4	<2 <0.4	<2 <0.4	<6 <1.2	<2 <0.4	5 <0.4	<2 <0.4	<6 <1.2
	μg/L										
4-Isopropyltoluene (Cymene)	μg/L	<0.4	<0.4	<0.4	<0.4	<0.4	<1.2	<0.4	<0.4	<0.4	<1.2
Methacrylonitrile	μg/L	<2	<2 <2	<2 <2	<2	<2 <2	<6	<2 <2	<2 <2	<2 <2	<6
Methyl acrylate	μg/L	<2			<2		<6				<6
Methyl-t-butylether (MTBE)	μg/L	<2	<2	<2	<2	<2	<6	<2	<2	<2	<6
Methylene chloride (Dichloromethane)	μg/L	<0.4	<0.4	<0.4	<0.4	<0.4	<1.2	<0.4	<0.4	<0.4	<1.2
Methyl methacrylate	μg/L	<2	<2	<2	<2	<2	<6	<2	<2	<2	<6
4-Methyl-2-pentanone Naphthalene	μg/L	<0.4	<0.4	<0.4	<0.4	<0.4	<1.2	<0.4	<0.4	<0.4	<1.2
Nitrobenzene	μg/L	<0.4 <2	<0.4 <2	<0.4 <2	<0.4 <2	<0.4 <2	<1.2 <6	<0.4 <2	<0.4 <2	<0.4 <2	<1.2 <6
2-Nitropropane	μg/L	<2 <2	<2 <2	<2 <2	<2 <2	<2 <2	<6	<2 <2	<2 <2	<2 <2	<6
Pentachloroethane	μg/L	<2 <2	<2	<2	<2	<2	<6	<2 <2	<2	<2	<6
Propionitrile (Ethyl cyanide)	μg/L	<10	<10	<2 <10	<10	<10	<30	<10	<2 <10	<10	<30
n-Propylbenzene	μg/L μg/L	<0.4	<0.4	<0.4	<0.4	<0.4	<1.2	<0.4	<0.4	<0.4	<1.2
Styrene	μg/L	<0.4	<0.4	<0.4	<0.4	<0.4	<1.2	<0.4	<0.4	<0.4	<1.2
1,1,1,2-Tetrachloroethane	μg/L μg/L	<0.4	<0.4	<0.4	<0.4	<0.4	<1.2	<0.4	<0.4	<0.4	<1.2
1.1.2.2-Tetrachloroethane	μg/L	<0.4	<0.4	<0.4	<0.4	<0.4	<1.2	<0.4	<0.4	<0.4	<1.2
Tetrachloroethene	μg/L μg/L	<0.4	<0.4	<0.4	<0.4	<0.4	<1.2	<0.4	<0.4	<0.4	<1.2
Tetrahydrofuran	μg/L	<2	<2	<2	<2	<2	<6	<2	<2	<2	<6
Toluene	μg/L	<0.4	<0.4	<0.4	<0.4	<0.4	<1.2	<0.4	<0.4	<0.4	<1.2
1,2,3-Trichlorobenzene	μg/L	<0.4	<0.4	<0.4	< 0.4	<0.4	<1.2	<0.4	<0.4	<0.4	<1.2
1,2,4-Trichlorobenzene	μg/L	<0.4	<0.4	<0.4	<0.4	<0.4	<1.2	<0.4	<0.4	<0.4	<1.2
1.1.1-Trichloroethane	μg/L	<0.4	<0.4	<0.4	<0.4	<0.4	<1.2	<0.4	<0.4	<0.4	<1.2
1,1,2-Trichloroethane	μg/L	<0.4	<0.4	<0.4	<0.4	<0.4	<1.2	<0.4	<0.4	<0.4	<1.2
Trichloroethene	μg/L	<0.4	<0.4	<0.4	<0.4	<0.4	<1.2	<0.4	<0.4	<0.4	<1.2
Trichlorofluoromethane	μg/L	<0.4	<0.4	<0.4	<0.4	<0.4	<1.2	<0.4	<0.4	<0.4	<1.2
1,2,3-Trichloropropane	μg/L	<0.4	<0.4	<0.4	<0.4	<0.4	<1.2	< 0.4	<0.4	<0.4	<1.2
1,2,4-Trimethylbenzene	μg/L μg/L	<0.4	<0.4	<0.4	<0.4	<0.4	<1.2	<0.4	<0.4	<0.4	<1.2
1,3,5-Trimethylbenzene	μg/L μg/L	<0.4	<0.4	< 0.4	<0.4	<0.4	<1.2	<0.4	<0.4	<0.4	<1.2
Vinvl Acetate	μg/L	<2	<2	<2	<2	<2	<6	<0.4	<2	<2	<6
Vinyl Chloride	μg/L	<2	<2	<2	<2	<2	<6	<2	<2	<2	<6
m&p-Xylene	μg/L μg/L	<0.4	<0.4	<0.4	<0.4	<0.4	<1.2	<0.4	<0.4	<0.4	<1.2
o-Xylene	μg/L	<0.4	<0.4	<0.4	<0.4	<0.4	<1.2	<0.4	<0.4	1.5	<1.2
o rejiene	μg/L	\U.4	\U. 1	\U.4	\U.4	\U. 4	1.6	\U. 1	V.12	1.3	\1.ω

Mission		ICC III E	4/E 22	r	C 22/E 24	1	G 2	2/E 25	TCC 10.4	/E 22	C 22/E 24
1/11/2010(1			4/Exp. 23		Soyuz 22/Exp. 24			3/Exp. 25	ISS 19A		Soyuz 22/Exp. 24
		WPA PWD Ambient	WPA PWD Ambient	WPA PWD Ambient	WPA PWD Hot	WPA PWD Ambient	WPA PWD Hot	WPA PWD Hot	PWD Aux Port	WPA RIP	WPA RIP
Sample Location		Ambient	Ambient	Ambient		Ambient			4		
Sample Location					1				Processed	Processed	
Sample Description		Potable Water	Potable Water	Potable Water	Potable Water	Potable Water	Potable Water	Potable Water	Water	Water	Processed Water
Sample Date		4/26/2010	5/18/2010	7/14/2010	8/25/2010	9/15/2010	10/19/2010	11/23/2010	2/25/2010	2/25/2010	7/29/2010
Analysis/Sample ID	Units	20100527006	20100527007	20100926001	20100926002	20100926003	20101129001	20101129002	20100422009	20100422008	20100926004
Volatile Organics - Non-Targets (Tentatively Iden	tified Comp								1		
Acetaldehyde	μg/L	not found	not found	not found	not found	not found	not found	not found	not found	not found	not found
Butyraldehyde (Butanal)	μg/L	not found	not found	not found	not found	not found	not found	not found	not found	not found	not found
1,8-Cineole	μg/L	not found	not found	not found	not found	not found	not found	not found	not found	not found	not found
Cyclohexanone	μg/L	not found	not found	not found	not found	not found	not found	not found	not found	not found	not found
Difluorodimethylsilane	μg/L	not found	not found	not found	not found	not found	not found	not found	not found	not found	not found
2,5-Dimethylfuran	μg/L	not found	not found	not found	not found	not found	not found	not found	not found	not found	not found
2,6-Dimethyl-1,7-octadiene	μg/L	not found	not found	not found	not found	not found	not found	not found	not found	not found	not found
Dimethyl sulfide (Thiobismethane)	μg/L	not found	not found	not found	not found	not found	not found	not found	not found	not found	not found
p-Dioxane	μg/L	not found	not found	not found	not found	not found	not found	not found	not found	not found	not found
1,3-Dioxolane (Ethylene glycol formal)	μg/L	not found	not found	not found	not found	not found	not found	not found	not found	not found	not found
Ethyl acetate	μg/L	not found	not found	not found	not found	not found	not found	not found	not found	not found	not found
Fluorotrimethylsilane	μg/L	not found	not found	not found	not found	not found	not found	not found	not found	not found	not found
2-Heptanone	μg/L	not found	not found	not found	not found	not found	not found	not found	not found	not found	not found
4-Heptanone	μg/L	not found	not found	not found	not found	not found	not found	not found	not found	not found	not found
3-Hexanone	μg/L	not found	not found	not found	not found	not found	not found	not found	not found	not found	not found
3-Hvdroxv-3-methylbutyric acid	μg/L	not found	not found	not found	not found	not found	not found	not found	not found	not found	not found
Isobutyronitrile	μg/L μg/L	not found	not found	not found	not found	not found	not found	not found	not found	not found	not found
Isooctanol	μg/L μg/L	not found	not found	not found	not found	not found	not found	not found	not found	not found	not found
L-Menthol	μg/L μg/L	not found	not found	not found	not found	not found	not found	not found	not found	not found	not found
Menthone	μg/L μg/L	not found	not found	not found	not found	not found	not found	not found	not found	not found	not found
Methyl acetate	μg/L	not found	not found	not found	not found	not found	not found	not found	not found	not found	not found
cis-1-Methyl-4-(1-methylethenyl)-cyclohexane	μg/L	not found	not found	not found	not found	not found	not found	not found	not found	not found	not found
3-Methyl-2-pentanone	μg/L	not found	not found	not found	not found	not found	not found	not found	not found	not found	not found
2-Methyl-1-propene	μg/L	not found	not found	not found	not found	not found	not found	not found	not found	not found	not found
alpha-Methyl styrene	μg/L	not found	not found	not found	not found	not found	not found	not found	not found	not found	not found
2-Nonanone	μg/L	not found	not found	not found	not found	not found	not found	not found	not found	not found	not found
2-Octanone	μg/L	not found	not found	not found	not found	not found	not found	not found	not found	not found	not found
2-Pentanone	μg/L	not found	not found	not found	not found	not found	not found	not found	not found	not found	not found
1,1,2-Trichloro-1,2,2-trifluoroethane (Freon 113)	μg/L	not found	not found	not found	not found	not found	not found	not found	not found	not found	not found
Trimethylsilanol	μg/L	not found	not found	not found	not found	not found	not found	not found	not found	not found	not found
									1		
Extractable Organics											
Acetophenone	μg/L	<16	<8	<8	<8	<16	<16	<16	<16	<16	<40
Benzaldehyde	μg/L	<8	<4	<4	<4	<8	<8	<8	<8	<8	<20
Benzoic acid	μg/L	<24	<12	<12	<12	<24	<48	<48	<24	<24	<60
Benzothiazole	μg/L	<8	<4	<4	<4	<8	<8	<8	<8	<8	<20
Benzyl alcohol	μg/L	<8	<4	<4	<4	<8	<8	<8	<8	<8	<20
Benzyl butyl phthlate	μg/L	<8	<4	<4	<4	<8	<8	<8	<8	<8	<20
2-Butoxyethanol	μg/L	<16	<8	<8	<8	<16	<16	<16	<16	<16	<40
2-(2-Butoxyethoxy)ethanol	μg/L	<16	<8	<8	<8	<16	<16	<16	<16	<16	<40
2-(2-Butoxyethoxy)ethyl acetate	μg/L	<8	<4	<4	<4	<8	<8	<8	<8	<8	<20
n-Butylpalmitate	μg/L	<16	<8	<8	<8	<16	<16	<16	<16	<16	<40
Butylated hydroxyanisole (BHA)	μg/L	<8	<4	<4	<4	<8	<8	<8	<8	<8	<20
N-Butylbenzenesulfonamide	μg/L	<8	<4	<4	<4	<8	<8	<8	<8	<8	<20
3-tert-Butylphenol	μg/L	<24	<12	<12	<12	<24	<24	<24	<24	<24	<60
Caffeine	μg/L	<8	<4	<4	<4	<8	<8	<8	<8	<8	<20

Mission		ISS III F	4/Exp. 23		Soyuz 22/Exp. 24	1	Sovuz 2	3/Exp. 25	ISS 19A	/Fvn 23	Soyuz 22/Exp. 24
		WPA PWD	WPA PWD	WPA PWD	WPA PWD Hot	WPA PWD		WPA PWD Hot	PWD Aux Port	WPA RIP	WPA RIP
		Ambient	Ambient	Ambient	WITT WD Hot	Ambient	WITH WE HOL	WITTI WE HOL	T WE HEAT ON	***********	***************************************
Sample Location											
									Processed	Processed	
Sample Description		Potable Water	Potable Water	Potable Water	Potable Water	Potable Water	Potable Water	Potable Water	Water	Water	Processed Water
Sample Date	TT	4/26/2010	5/18/2010	7/14/2010	8/25/2010	9/15/2010	10/19/2010	11/23/2010	2/25/2010	2/25/2010	7/29/2010
Analysis/Sample ID	Units	20100527006	20100527007	20100926001	20100926002	20100926003	20101129001	20101129002	20100422009	20100422008	20100926004
tris-2-Chloroethyl phosphate	μg/L	<8	<4	<4	<4	<8	<8	<8	<8	<8	<20
Cholesterol o-Cresol (2-Methylphenol)	μg/L	<64 <8	<32	<32 <4	<32	<64 <8	<64 <8	<64	<64 <8	<64 <8	<160 <20
Cyclododecane	μg/L μg/L	<8 <8	<4	<4	<4 <4	<8 <8	<8 <8	<8 <8	<8 <8	<8 <8	<20 <20
Decamethylcyclopentasiloxane	μg/L μg/L	<8 <8	<4 <4	<4	<4 <4	<8 <8	<8 <8	<8 <8	<8 <8	<8 <8	<20 <20
Decamenyicyciopentasnoxane Decanoic acid	μg/L μg/L	<8 <16	<4 <8	<4 <8	<4 <8	<8 <16	<8 <24	<8 <24	<8 <16	<8 <16	<40
2,6-Di-t-butyl-1,4-benzoquinone	μg/L μg/L	<10 <8	<8 <4	< 6	<8 <4	<16 <8	<24 <8	<24 <8	<10 <8	<10 <8	<20
2.4-Di-t-butylphenol	μg/L μg/L	<8	<4	<4	<4	<8	<8	<8	<8	<8	<20
1.4 Diacetylbenzene	μg/L μg/L	<8	<4	<4	<4	<8	<8	<8	<8	<8	<20
N,N-Dibutylformamide	μg/L μg/L	<8	<4	<4	<4 <4	<8	<8	<8	<8 <8	<8 <8	<20
Dibutyl phthalate	μg/L μg/L	<8	<4	<4	<4	<8	<8	<8	<8	<8 <8	<20
Dibutylamine	μg/L μg/L	<8	<4	<4	<4 <4	<8	<8	<8	<8 <8	<8 <8	<20
N,N-Diethyl-m-toluamide	μg/L μg/L	<8	<4	<4	<4	<8	<8	<8	<8	<8	<20
Diethylphthalate	μg/L μg/L	<8	<4	<4	<4	<8	<8	<8	<8	<8	<20
Diethylene glycol monoethyl ether	μg/L μg/L	<8	<4	<4	<4	<8	<8	<8	<8	<8	<20
N,N-Diethylformamide	μg/L μg/L	<24	<12	<12	<12	<24	<24	<24	<24	<24	<60
Diiodomethane (Methyl iodide)	μg/L μg/L	<8	<4	<4	<4	<8	<8	<24 <8	<8	<8	<20
Diisopropyl adipate	μg/L μg/L	<8	<4	<4	<4	<8	<8	<8	<8	<8	<20
Dimethyl phthalate	μg/L μg/L	<8	<4	<4	<4	<8	<8	<8	<8	<8	<20
N,N-Dimethyl acetamide	μg/L	<8	<4	<4	<4	<8	<8	<8	<8	<8	<20
N.N-Dimethyl dectainide N.N-Dimethylbenzylamine	μg/L μg/L	<8	<4	<4	<4	<8	<8	<8	<8	<8	<20
N,N-Dimethylformamide	μg/L	<16	<8	<8	<8	<16	<16	<16	<16	<16	<40
Dipropylene glycol methyl ether	μg/L	<8	<4	<4	<4	<8	<8	<8	<8	<8	<20
Dodecamethylcyclohexasiloxane	μg/L	<8	<4	<4	<4	<8	<8	<8	<8	<8	<20
2-Ethoxyethanol	μg/L	<8	<4	<4	<4	<8	<16	<16	<8	<8	<20
2-Ethyl-1-hexanol	μg/L	<8	<4	<4	<4	<8	<8	<8	<8	<8	<20
2-Ethylhexanoic acid	μg/L	<8	<4	<4	<4	<8	<16	<16	<8	<8	<20
bis-2-Ethylhexyl adipate	μg/L	<8	<4	<4	<4	<8	<8	<8	<8	<8	<20
bis-2-Ethylhexyl phthalate (Dioctyl phthlate)	μg/L	<8	<4	<4	<4	<8	<8	<8	<8	<8	<20
4-Ethylmorpholine	μg/L	<8	<4	<4	<4	<8	<8	<8	<8	<8	<20
1-Formylpiperidine	μg/L	<8	<4	<4	<4	<8	<8	<8	<8	<8	<20
Heptanoic acid	μg/L	<8	<4	<4	<4	<8	<24	<24	<8	<8	<20
2-Heptanone	μg/L	<8	<4	<4	<4	<8	<8	<8	<8	<8	<20
gamma-Hexalactone	μg/L	<8	<4	<4	<4	<8	<8	<8	<8	<8	<20
Hexanoic acid	μg/L	<16	<8	<8	<8	<16	<24	<24	<16	<16	<40
2-Hexanol	μg/L	<8	<4	<4	<4	<8	<8	<8	<8	<8	<20
2-Hydroxybenzothiazole	μg/L	<8	<4	<4	<4	<8	<8	<8	<8	<8	<20
Ibuprofen	μg/L	<8	<4	<4	<4	<8	<24	<24	<8	<8	<20
Iodoform	μg/L	<8	<4	<4	<4	<8	<8	<8	<8	<8	<20
Isophorone	μg/L	<8	<4	<4	<4	<8	<8	<8	<8	<8	<20
4-Isopropylphenol	μg/L	<8	<4	<4	<4	<8	<8	<8	<8	<8	<20
Lauramide	μg/L	<8	<4	<4	<4	<8	<8	<8	<8	<8	<20
Lauric acid (Dodecanoic acid)	μg/L	<240	<120	<120	<120	<240	<240	<240	<240	<240	<600
p-Menth-1-en-8-ol (alpha-Terpineol)	μg/L	<8	<4	<4	<4	<8	<8	<8	<8	<8	<20
2-Mercaptobenzothiazole	μg/L	<80	<40	<40	<40	<80	<80	<80	<80	<80	<200
2-Methyl-2,4-pentanediol	μg/L	<8	<4	<4	<4	<8	<8	<8	<8	<8	<20
1-Methyl-2-pyrrolidinone	μg/L	<8	<4	<4	<4	<8	<8	<8	<8	<8	<20
Methyl-4-hydroxybenzoate	μg/L	<8	<4	<4	<4	<8	<8	<8	<8	<8	<20

Mission		ISS III E	4/Exp. 23		Sovuz 22/Exp. 24	1	Sovuz 2	3/Exp. 25	ISS 19A	/Evn 23	Sovuz 22/Exp. 24
		WPA PWD	WPA PWD	WPA PWD	WPA PWD Hot	WPA PWD		WPA PWD Hot	PWD Aux Port	WPA RIP	WPA RIP
		Ambient	Ambient	Ambient	WFA F WD HOL	Ambient	WFAFWD Hot	WFAFWD HOL	T WD Aux Foit	WFARIF	WIAKII
Sample Location		Ambient	Ambient	Ambient		Ambient					
									Processed	Processed	
Sample Description		Potable Water	Potable Water	Potable Water	Potable Water	Potable Water	Potable Water	Potable Water	Water	Water	Processed Water
Sample Date		4/26/2010	5/18/2010	7/14/2010	8/25/2010	9/15/2010	10/19/2010	11/23/2010	2/25/2010	2/25/2010	7/29/2010
Analysis/Sample ID	Units	20100527006	20100527007	20100926001	20100926002	20100926003	20101129001	20101129002	20100422009	20100422008	20100926004
Methyl sulfone	μg/L	40	34	36	<4	<8	48	50	88	98	150
2-Methyl butyric acid	μg/L	<24	<12	<12	<12	<24	<24	<24	<24	<24	<60
2-Methylthiobenzothiazole	μg/L	<8	<4	<4	<4	<8	<8	<8	<8	<8	<20
Monomethyl phthalate	μg/L	<8	<4	<4	<4	<8	<8	<8	<8	<8	<20
Myristic acid	μg/L	<48	<24	<24	<24	<48	<64	<64	<48	<48	<120
(+)-Neomenthol	μg/L	<8	<4	<4	<4	<8	<8	<8	<8	<8	<20
Nicotine	μg/L	<8	<4	<4	<4	<8	<8	<8	<8	<8	<20
Nonadecane	μg/L	<8	<4	<4	<4	<8	<8	<8	<8	<8	<20
Nonanoic acid	μg/L	<24	<12	<12	<12	<24	<24	<24	<24	<24	<60
1-Octadecanol	μg/L	<24	<12	<12	<12	<24	<24	<24	<24	<24	<60
Octamethylcyclotetrasiloxane	μg/L	<8	<4	<4	<4	<8	<8	<8	<8	<8	<20
Octanoic acid	μg/L	<16	<8	<8	<8	<16	<48	<48	<16	<16	<40
4-tert-Octylphenol	μg/L	<8	<4	<4	<4	<8	<8	<8	<8	<8	<20
Oleic acid	μg/L	<80	<40	<40	<40	<80	<80	<80	<80	<80	<200
Oxindole	μg/L	<8	<4	<4	<4	<8	<8	<8	<8	<8	<20
Palmitic acid	ug/L	<240	<120	<120	<120	<240	<240	<240	<240	<240	<600
Palmitoleic acid	ug/L	NA	NA	NA	NA	NA	NA	NA	NA	NA	NA
Pentacosane	μg/L	<8	<4	<4	<4	<8	<8	<8	<8	<8	<20
sec-Phenethyl alcohol	μg/L	<8	<4	<4	<4	<8	<8	<8	<8	<8	<20
Phenol	μg/L	<8	<4	<4	<4	<8	<8	<8	<8	<8	<20
2-Phenoxyethanol	μg/L	<8	<4	<4	<4	<8	<8	<8	<8	<8	<20
N-Phenyl-2-naphthylamine	μg/L	<8	<4	<4	<4	<8	<8	<8	<8	<8	<20
2-Phenyl-2-propanol	μg/L	<8	<4	<4	<4	<8	<8	<8	<8	<8	<20
2-Phenylacetic acid	μg/L	<32	<16	<16	<16	<32	<32	<32	<32	<32	<80
Phenethyl alcohol	μg/L	<8	<4	<4	<4	<8	<8	<8	<8	<8	<20
2-Phenylphenol	μg/L	<8	<4	<4	<4	<8	<8	<8	<8	<8	<20
Salicyclic Acid	μg/L	<64	<32	<32	<32	<64	<64	<64	<64	<64	<160
trans-Squalene	μg/L	<16	<8	<8	<8	<16	<16	<16	<16	<16	<40
Stearic acid	μg/L	NA	NA	NA	NA	NA	NA	NA	NA	NA	NA
1-Tetradecanol	μg/L	<8	<4	<4	<4	<8	<8	<8	<8	<8	<20
Tetramethylsuccinonitrile	μg/L	<8	<4	<4	<4	<8	<8	<8	<8	<8	<20
Tetramethyl thiourea	μg/L	<8	<4	<4	<4	<8	<8	<8	<8	<8	<20
Tetramethylurea	μg/L	<8	<4	<4	<4	<8	<8	<8	<8	<8	<20
Thymol	μg/L	<8	<4	<4	<4	<8	<8	<8	<8	<8	<20
1,3,5-Triallyl-1,3,5-triazine-2,4,6(1H,3H,5H)-trione	μg/L	<8	<4	<4	<4	<8	<8	<8	<8	<8	<20
Tributylamine	μg/L	<8	<4	<4	<4	<8	<8	<8	<8	<8	<20
Tributyl phosphate	μg/L	<8	<4	<4	<4	<8	<8	<8	<8	<8	<20
Triethyl phosphate	μg/L	<16	<8	<8	<8	<16	<16	<16	<16	<16	<40
2,2,4-Trimethyl-1,3-pentanediol diisobutyrate	ug/L	<16	<8	<8	<8	<16	<16	<16	<16	<16	<40
Tripropylene glycol monomethyl ether	ug/L	<8	<4	<4	<4	<8	<8	<8	<8	<8	<20
Undecanoic acid	μg/L	<48	<24	<24	<24	<48	<48	<48	<48	<48	<120
2-Undecanone	μg/L μg/L	<8	<4	<4	<4	<8	<8	<8	<8	<8	<20
Valeric acid (Pentanoic acid)	μg/L	<48	<24	<24	<24	<48	<48	<48	<48	<48	<120
Vanillin	μg/L	<16	<8	<8	<8	<16	<16	<16	<16	<16	<40
	rs/L	-10		, i	-0	-10	-10	-10	-10	-10	-10

Mission	1	ICC III E	4/Exp. 23		Sovuz 22/Exp. 24	1	Corner 2	2/E-m 25	TCC 10.4	/Exp. 23	Sovuz 22/Exp. 24
and and a second		WPA PWD	WPA PWD	WPA PWD	WPA PWD Hot	WPA PWD		3/Exp. 25 WPA PWD Hot	PWD Aux Port	WPA RIP	WPA RIP
		WPA PWD Ambient	WPA PWD Ambient	WPA PWD Ambient	WYA YWD HOU	WPA PWD Ambient	WPA PWD Hot	WPA PWD HOT	r WD Aux Port	WYA KIP	WYA KIP
Sample Location		Ambient	Ambient	Ambient		Ambient			1		
									Processed	Processed	
Sample Description		Potable Water	Potable Water	Potable Water	Potable Water	Potable Water	Potable Water	Potable Water	Water	Water	Processed Water
Sample Date		4/26/2010	5/18/2010	7/14/2010	8/25/2010	9/15/2010	10/19/2010	11/23/2010	2/25/2010	2/25/2010	7/29/2010
Analysis/Sample ID	Units	20100527006	20100527007	20100926001	20100926002	20100926003	20101129001	20101129002	20100422009	20100422008	20100926004
Acid Extractables-EPA 625 List									1		
4-Chloro-3-methylphenol	μg/L	<16	<8	<8	<8	<16	<16	<16	<16	<16	<40
2-Chlorophenol	μg/L	<16	<8	<8	<8	<16	<16	<16	<16	<16	<40
2,4-Dichlorophenol	μg/L	<16	<8	<8	<8	<16	<16	<16	<16	<16	<40
2,4-Dimethylphenol	μg/L	<16	<8	<8	<8	<16	<16	<16	<16	<16	<40
2,4-Dinitrophenol	μg/L	<16	<8	<8	<8	<16	<16	<16	<16	<16	<40
2-Methyl-4,6-dinitrophenol	μg/L	<16	<8	<8	<8	<16	<16	<16	<16	<16	<40
2-Nitrophenol	μg/L	<16	<8	<8	<8	<16	<16	<16	<16	<16	<40
4-Nitrophenol	μg/L	<16	<8	<8	<8	<16	<16	<16	<16	<16	<40
Pentachlorophenol	μg/L	<16	<8	<8	<8	<16	<16	<16	<16	<16	<40
Phenol	μg/L	<8	<4	<4	<4	<8	<8	<8	<8	<8	<20
2,4,5-Trichlorophenol	μg/L	<16	<8	<8	<8	<16	<16	<16	<16	<16	<40
2,4,6-Trichlorophenol	μg/L	<16	<8	<8	<8	<16	<16	<16	<16	<16	<40
4-Methylphenol	μg/L	<8	<4	<4	<4	<8	<8	<8	<8	<8	<20
D Al									1		
Base/Neutral Extractables - EPA 625 List							40	4.0	40		
Benzidine	μg/L	<16	<8	<8	<8	<16	<16	<16	<16	<16	<40
3,3-Dichlorobenzidine	μg/L	<16	<8	<8	<8	<16	<16	<16	<16	<16	<40
bis-(2-Ethylhexyl)phthalate	μg/L	<8	<4	<4.0	<4.0	<8	<8	<8	<8	<8	<20
Benzyl butyl phthalate Dibutylphthalate	μg/L	<8	<4	<4	<4	<8	<8	<8	<8	<8	<20
Diethylphthalate	μg/L	<8 <8	<4	<4 <4	<4	<8 <8	<8	<8	<8	<8 <8	<20
Dimethylphthalate	μg/L μg/L	<8	<4 <4	<4	<4 <4	<8	<8 <8	<8 <8	<8 <8	<8	<20 <20
Di-n-octyl phthalate	μg/L μg/L	<8 <16	<4 <8	<4 <8	<4 <8	<8 <16	<8 <16	<8 <16	<8 <16	<8 <16	<20 <40
N-Nitrosodimethylamine	μg/L μg/L	<16	<8	<8	<8	<16	<16	<16	<16	<16	<40
N-Nitrosodiphenylamine	μg/L μg/L	<16	<8	<8	<8	<16	<16	<16	<16	<16	<40
N-Nitrosodi-n-propylamine	μg/L μg/L	<16	<8	<8	<8	<16	<16	<16	<16	<16	<40
2,4-Dinitrotoluene	μg/L μg/L	<16	<8	<8	<8	<16	<16	<16	<16	<16	<40
2,6-Dinitrotoluene	μg/L μg/L	<16	<8	<8	<8	<16	<16	<16	<16	<16	<40
Isophorone	μg/L μg/L	<8	<4	<4	<4	<8	<8	<8	<8	<8	<20
Nitrobenzene	μg/L μg/L	<16	<8	<8	<8	<16	<16	<16	<16	<16	<40
Acenaphthene	μg/L	<16	<8	NA NA	NA NA	NA.	<16	<16	<16	<16	NA NA
Acenaphthylene	μg/L	<16	<8	<8	<8	<16	<16	<16	<16	<16	<40
Anthracene	μg/L	<16	<8	<8	<8	<16	<16	<16	<16	<16	<40
Benzo(a)anthracene	μg/L	<16	<8	<8	<8	<16	<16	<16	<16	<16	<40
Benzo(a)pyrene	μg/L	<10	<5	<5	<5	<10	<10	<10	<10	<10	<25
Benzo(b)fluoranthene	μg/L	<8	<4	<4	<4	<8	<8	<8	<8	<8	<20
Benzo(ghi)perylene	μg/L	<10	<5	<5	<5	<10	<10	<10	<10	<10	<25
Benzo(k)fluoroanthene	μg/L	<8	<4	<4	<4	<8	<8	<8	<8	<8	<20
Chrysene	μg/L	<20	<10	<10	<10	<20	<20	<20	<20	<20	<50
Dibenzo(a,h)anthracene	μg/L	<10	<5	<5	<5	<10	<10	<10	<1	<10	<25
Fluoranthene	μg/L	<8	<4	<4	<4	<8	<8	<8	<8	<8	<20
Fluorene	μg/L	<16	<8	<8	<8	<16	<16	<16	<16	<16	<40
Indeno(1,2,3-cd)pyrene	μg/L	<10	<5	<5	<5	<10	<10	<10	<10	<10	<25
Naphthalene	μg/L	<40	<20	<20	<20	<40	<40	<40	<40	<40	<100
Phenanthrene	μg/L	<8	<4	<4	<4	<8	<8	<8	<8	<8	<20
Pyrene	μg/L	<8	<4	<4	<4	<8	<8	<8	<8	<8	<20
bis(2-Chloroethyl) ether	μg/L	<16	<8	<8	<8	<16	<16	<16	<16	<16	<40

Mission	1	ICC III E	4/Exp. 23		Sorma 22/Evm 2/	1	Sovar 2	3/Exp. 25	ISS 19A	/Evn. 22	Soyuz 22/Exp. 24
1 TALLES AND A STATE OF THE STA					Soyuz 22/Exp. 24				PWD Aux Port		WPA RIP
		WPA PWD Ambient	WPA PWD Ambient	WPA PWD Ambient	WPA PWD Hot	WPA PWD Ambient	WPA PWD Hot	WPA PWD Hot	PWD Aux Port	WPA RIP	WPA RIP
Sample Location		Ambient	Ambient	Ambient		Ambient					
									Processed	Processed	
Sample Description		Potable Water	Potable Water	Potable Water	Potable Water	Potable Water	Potable Water	Potable Water	Water	Water	Processed Water
Sample Date		4/26/2010	5/18/2010	7/14/2010	8/25/2010	9/15/2010	10/19/2010	11/23/2010	2/25/2010	2/25/2010	7/29/2010
Analysis/Sample ID	Units	20100527006	20100527007	20100926001	20100926002	20100926003	20101129001	20101129002	20100422009	20100422008	20100926004
bis(2-Chloroethoxy) methane	μg/L	<16	<8	<8	<8	<16	<16	<16	<16	<16	<40
bis(2-Chloroisopropyl) ether	μg/L	<16	<8	<8	<8	<16	<16	<16	<16	<16	<40
4-Bromophenyl phenyl ether	μg/L	<16	<8	<8	<8	<16	<16	<16	<16	<16	<40
4-Chlorophenyl phenyl ether	μg/L	<16	<8	<8	<8	<16	<16	<16	<16	<16	<40
2-Chloronaphthalene	μg/L	<16	<8	<8	<8	<16	<16	<16	<16	<16	<40
1,2-Dichlorobenzene	μg/L	<16	<8	<8	<8	<16	<16	<16	<16	<16	<40
1,3-Dichlorobenzene	μg/L	<16	<8	<8	<8	<16	<16	<16	<16	<16	<40
1,4-Dichlorobenzene	μg/L	<16	<8	<8	<8	<16	<16	<16	<16	<16	<40
Hexachlorobenzene	μg/L	<16	<8	<8	<8	<16	<16	<16	<16	<16	<40
Hexachlorobutadiene	μg/L	<16	<8	<8	<8	<16	<16	<16	<16	<16	<40
Hexachlorocyclopentadiene	μg/L	<16	<8	<8	<8	<16	<16	<16	<16	<16	<40
Hexachloroethane	μg/L	<16	<8	<8	<8	<16	<16	<16	<16	<16	<40
1,2,4-Trichlorobenzene	μg/L	<16	<8	<8	<8	<16	<16	<16	<16	<16	<40
Alcohols (DAI/GC/MS)											
1-Butanol	μg/L	<100	<100	<100	<100	<100	<100	<100	<100	<100	<100
2-Butanol	μg/L	<100	<100	<100	<100	<100	<100	<100	<100	<100	<100
Ethanol	μg/L	<100	<100	<100	<100	<100	<100	<100	<100	<100	<100
Methanol	μg/L	<100	<100	<100	<100	<100	<100	<100	<100	<100	<100
2-Methyl-1-butanol	μg/L	<100	<100	<100	<100	<100	<100	<100	<100	<100	<100
2-Methyl-2-butanol	μg/L	<100	<100	<100	<100	<100	<100	<100	<100	<100	<100
3-Methyl-1-butanol (Isopentanol)	μg/L	<100	<100	<100	<100	<100	<100	<100	<100	<100	<100
2-Methyl-1-propanol	μg/L	<100	<100	<100	<100	<100	<100	<100	<100	<100	<100
2-Methyl-2-propanol	μg/L	<100	<100	<100	<100	<100	<100	<100	<100	<100	<100
1-Pentanol (Amyl alcohol)	μg/L	<100	<100	<100	<100	<100	<100	<100	<100	<100	<100
2-Pentanol (sec-Amyl alcohol)	μg/L	<100	<100	<100	<100	<100	<100	<100	<100	<100	<100
3-Pentanol	μg/L	<100	<100	<100	<100	<100	<100	<100	<100	<100	<100
1-Propanol	μg/L	<100	<100	<100	<100	<100	<100	<100	<100	<100	<100
2-Propanol (Isopropanol)	μg/L	<100	<100	<100	<100	<100	<100	<100	<100	<100	<100
Glycols (DAI/GC/MS)											
1,2-Ethanediol (Ethylene glycol)	μg/L	<1000	<1000	<1000	<1000	<1000	<1000	<1000	<1000	<1000	<1000
1,2-Propanediol (Propylene glycol)	μg/L	< 500	< 500	< 500	< 500	< 500	< 500	< 500	<500	< 500	<500
Glycerol (LC/MS/MS)											
Glycerol (1,2,3-Propanetriol)	μg/L	NA	NA	<300	<300	<300	<300	<300	NA	NA	<300
Silanols (GC/MS & LC/MS/MS) (R&D Metho	nd - NIST traces										
Dimethylsilanediol (DMSD)	μg/L	NA	NA	6100	8240	4730	8490	<400	NA	NA	7340
Carboxylates (CE)					 						
Acetate	μg/L	<125	<125	<125	<125	<125	<125	<125	<125	<125	<125
Formate	μg/L μg/L	<125	<125	<125	<125	<125	<125	<125	<125	<125	<125
Glycolate	μg/L μg/L	<125	<125	<125	<125	<125	<125	<125 <125	<125	<125	<125
Glyoxylate	μg/L μg/L	<125	<125	<125	<125	<125	<125	<125	<125	<125	<125
Lactate	μg/L μg/L	<1000	<1000	<1000	<1000	<1000	<1000	<1000	<1000	<1000	<1000
Oxalate	μg/L μg/L	<125	<125	<125	<125	<125	<125	<125	<125	<125	<125
Propionate	μg/L μg/L	<125	<125	<125	<125	<125	<125	<125	<125	<125	<125

Mission		ISS ULF	4/Exp. 23		Soyuz 22/Exp. 24	1	Soyuz 2	3/Exp. 25	ISS 19A	/Exp. 23	Soyuz 22/Exp. 24
		WPA PWD	WPA PWD	WPA PWD	WPA PWD Hot	WPA PWD	WPA PWD Hot	WPA PWD Hot	PWD Aux Port	WPA RIP	WPA RIP
		Ambient	Ambient	Ambient		Ambient					
Sample Location									Processed	Processed	
Sample Description		Potable Water	Potable Water	Potable Water	Potable Water	Potable Water	Potable Water	Potable Water	Water	Water	Processed Water
Sample Date		4/26/2010	5/18/2010	7/14/2010	8/25/2010	9/15/2010	10/19/2010	11/23/2010	2/25/2010	2/25/2010	7/29/2010
Analysis/Sample ID	Units	20100527006	20100527007	20100926001	20100926002	20100926003	20101129001	20101129002	20100422009	20100422008	20100926004
Aldehydes											
Formaldehyde	μg/L	6	3	8	5	6	<5	<5	3	6	8
Amines (CE)											
Ethylamine	μg/L	<125	<125	<125	<125	<125	<125	<125	<125	<125	<125
Methylamine	μg/L	<125	<125	<125	<125	<125	<125	<125	<125	<125	<125
n-Propylamine	μg/L	<125	<125	<125	<125	<125	<125	<125	<125	<125	<125
Trimethylamine	μg/L	<125	<125	<125	<125	<125	<125	<125	<125	<125	<125
Non-volatiles (LC/UV-VIS)											
Urea	μg/L	<800	<800	<800	<800	<800	<800	<800	<800	<800	<800
Caprolactam	μg/L	<8	<4	<4	<4	<8	<16	<16	<8	<8	<20
Organic Carbon Recovery	percent	12.07	4.43	106.01	98.19	111.42	90.19	8.68	16.63	16.44	89.58
Unaccounted Organic Carbon	mg/L	0.14	0.21	0.00	0.04	0.00	0.25	0.13	0.20	0.15	0.23